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Tracheal Tubes

1. General recommendations

Attention: Read these instructions carefully before using your tracheal tube for the first time.

Before initial use, every subsequent use and before sending the product back for repair, maintenance or service, the tracheal tube (hereinafter simply named "tube") must first be cleaned and disinfected in accordance with these instructions.

2. Intended use

Every tube is intended for only one patient. Please bear in mind that if several tubes are used, the components of the individual tubes may not be exchanged with one another. The treating physician is responsible for selecting the suitable type and size of tube and correct tube insertion, as well as for providing adequate patient training and information for handling and using the tube.

This product allows access to, and control of, the respiratory tract. The windowed tracheotomy tubes are indicated as well if an effective and gentle weaning of the patient from the tracheotomy tube is desired. If the sealing cap is used with the windowed tracheotomy tubes, the air is passed through the window and along the tube to the upper respiratory tract and through the vocal chords.

Attention: The use of windowless speaking tubes with a valve (so-called speaking valve tubes with a round or flat valve) can lead to death by suffocation. Therefore, only speaking valve tubes with window(s) (i.e., a window both in the outer and inner tube, or with a mesh filter in the outer tube and a window in the inner tube) may be used.

Attention: The subsequent processing of surface-treated and silver-plated tubes, including tube perforation, can damage the protective coating or silver-plating or create unfinished edges. No warranty is given for the improper execution of such works.

3. Testing before re-use

Every time before re-using the tube, it should be checked for fractures, cracks, deformations, damage and functionality. Special attention should be given to the tube end leading into the respiratory tract, as it should always be smooth and undamaged. Otherwise, injuries to the respiratory tract cannot be ruled out. Worn, corroded, deformed and porous tubes, or those with a changed surface or otherwise damaged, must be taken out of service or sent for repair.

4. Tube structure

HEBU tracheal tubes are manufactured in various models and sizes.

Apart from various lengths and diameters, special models with a window in both the outer and inner tube, or with a mesh filter in the outer tube and a window in the inner tube, or with a mesh filter in both the outer and inner tube are offered. Our medical product specialists are ready to answer your questions regarding details about tracheal tube application, shape and design. Manufacturing from a seamless tube facilitates cleaning.

5. Regular maintenance of the tube

Every six months, the tube should be sent to technical service for maintenance and overhaul.

6. Care and cleaning

So the product can be cleaned without time pressure, we always recommend using two tubes interchangeably. As long as one tube is being worn, the other one is being cleaned. A tip for initial use:

Simulate all exchange and cleaning steps in a "dry exercise" to ensure that you forgot no materials. You will gain self-confidence and become familiarized with the routine.

To minimize the risk of infection, every tube must be cleaned at least twice a day, and in case of stronger secretion, even more often. Tubes are manufactured from different materials according to their intended use and owing to their various desired properties.

Attention: Never, under any circumstances, clean the tubes with abrasive household cleaners, denture cleaners or high-percentage alcohol.

6.1 Tubes made of silver and German silver

Apart from an aesthetically pleasing aspect, a silver tube has the advantage of being particularly thin-walled, long-lived and easy to clean.

For cleaning silver tubes, warm water, a soft brush, and, if need be, some mild soap solution are sufficient.

Since silver typically turns brownish or black, these discolorations are no reason for reclamation. If necessary, the tube can be cleaned in a silver immersion bath, but the manufacturer's instructions should be followed. After such a cleaning, the tube must be thoroughly rinsed several times under warm, running water.

Tubes made of silver should never be immersed in a physiological salt solution, since prolonged contact may lead to the formation of silver chloride and possible ugly discolorations.

6.2 Tubes made of stainless steel

The types of stainless, high-grade steel used in the manufacturing of tubes form specific passive layers that have a protective function due to their alloy. These steels are only partly resistant to the attack of chloride ions and aggressive media and fluids!

In addition to the manufacturer's efforts to select the right materials and to carry out their careful processing, the user must take good and continuous care of his/her tubes, which may include their disinfection or sterilization.

We recommend the following processes for the reprocessing of our reusable tubes:

6.3 Manual cleaning

If possible, the tube should be disinfected and cleaned immediately upon use. Dirt should not be allowed to dry on the objects because it would make cleaning even more difficult.

Please take into account the following:

- The solutions for manual cleaning must be used according to the manufacturer's instructions.
- Use a suitable, soft brush for cleaning the tube, so every spot is reached.
- Remove blood or other incrustations or residues with a soft brush and a neutral or mildly alkaline cleaner. If hard crusts have formed, place them in warm water to soften them up before cleaning the tube with the brush.
- Never use brushes or sponges made of metal for manual cleaning.
- Take apart the tube as much as possible before cleaning.
- Residues from the cleaning stage must be fully removed in the subsequent final rinse.
- Afterwards, the external and internal tubes are dried. Before the internal tube is introduced into the external tube, we recommend rubbing off the outer surface of the internal tube with paraffin oil in order to facilitate the gliding of both tubes.

6.4 Machine cleaning

The cleaning and disinfection in a cleaning machine should always be given priority to manual cleaning, since machine processes can be standardized.

Always observe the operation and loading instructions of the respective manufacturer. Use only cleaners recommended by the manufacturer of your cleaning machine for the respective application purposes.

- Always load the tubes into the machine in disassembled condition in such a way that water can flow out.
- Make sure that the tubes are correctly positioned so rinsing reaches all parts.
- When removing the tubes from the cleaning machine, check them for visible dirt. If need be, repeat the cleaning cycle or carry out a manual pre-cleaning in the first place.
- Residues from the cleaning cycle must be completely removed by the subsequent final rinsing. Please use only fully demineralized or distilled water for this rinsing step.

6.5 Ultrasound treatment

To clean in the ultrasound bath, the tubes must be stored in an open position on wire baskets specially made for cleaning. Since warm water without any extra ingredients does not clean very well, a suitable cleaner must be added to the water. The manufacturer's recommendations with regard to concentration and temperature of the cleaning solutions in the ultrasound tank must be observed. Excessive dirt will lead to insufficient cleaning, so exchange the cleaning solution regularly according to the manufacturer's specifications. In addition, the bath and ultrasound exposure times according to the manufacturer's specifications for the cleaner used must be observed.

Basically, tubes cleaned by ultrasound must be subjected to a final rinsing. After the ultrasound treatment, it is necessary to inspect the tubes for loose components and to secure them again if need be. Rinsing should be carried out, if possible, with fully demineralized or distilled water to prevent water spots.

6.6 Thermal disinfection

For thermal disinfection, the tube can be disinfected in boiling water for at least 15 minutes. To prevent air bubbles, it is important to place the tube in the bath with the curved section (neck) pointing downwards. After disinfection has been completed, drain off the water. Afterwards, touch the tube only at the "shield" or neck-plate or with sterile, single-use gloves.

6.7 Chemical disinfection

1. The chemical disinfection solutions must be used according to the solvent manufacturer's specifications.
 2. Prior to use, the chemicals must be diluted with pure water to prepare the working solution. The addition of cleaners is not allowed. When using chemicals, make sure to observe the manufacturers' specifications (exposure time and concentration).
 3. Disinfecting solutions should be prepared daily for immediate use. In case of multiple use, the following problems may arise:
 - a. Increase in concentration owing to evaporation (danger of corrosion) or
 - b. Excessive dirt (corrosion and reduced effectiveness).
 4. After disinfection, always rinse again with enough clean, running water. To prevent water spots, use fully demineralized water.6.8 Drying of the tube
- Make sure to dry the tubes very well right after completing the cleaning or rinsing cycles. Afterwards, touch the tube only at the neck-plate or with sterile, single-use gloves.

7. Testing after cleaning

After every cleaning, the tube should be macroscopically clean, i.e., there should be no visible dirt and no traces of disinfectants or cleaners.

Test all moving parts very carefully.

If you notice that the tube malfunctions or is damaged, withdraw it immediately from use.

8. Processing limitations, disposal

Frequent reprocessing has no significant effect on the tube. Wear and damage during use are usually the main factors determining the product's useful life. Once the product has reached the end of its useful life, please take the tube to the proper place of disposal or to a recycling center.

9. How to insert the tube

The following recommendations are only of a general nature. Please consult the treating physician or the hospital's specialized staff for the exact handling of your tube.

Before attempting to use the product on your own, please rehearse its handling under professional supervision.

It is best to insert the tube standing straight in front of a mirror and grasping the tube's neck-plate with both hands. Then, place the tube's tip on the tracheostomy. The tube's neck-plate should point downwards. Insert the tube carefully into the tracheostomy, following the course of its curve. The use of an obturator (available as an optional accessory for your tube) and of stoma oil may facilitate the insertion of the tube.

Make sure that the tube has a tracheostomy protective cover.

To use the speaking insert, turn the tube upside down to avoid possible valve damage. Please make sure that by using additional aids and creating suitable environmental conditions the air that you breathe in is within the indicated temperature and humidity ranges and that no foreign bodies can be breathed in.

10. Contraindications

- ⇒ Unconscious and/or comatose patients
- ⇒ Severe respiratory tract blockage that can prevent adequate expiration
- ⇒ Severe and increased mucous secretion
- ⇒ Significantly reduced pulmonary elasticity that can cause air trapping
- ⇒ Difficult aspiration

In these cases, please consult your treating physician.

11. Exclusion of liability

Safety notice: Responsibility for the proper cleaning, disinfection and, if need be, sterilization of the tubes lies with the product user. Be sure to observe all local regulations!

HEBU medical, as manufacturer of the products, shall not be liable for direct or consequential damage caused by improper use or handling or by incorrect cleaning, sterilization or maintenance. If companies or individuals repair the tubes without the authorization of HEBU medical, the warranty shall cease to apply.

Failure to observe the recommendations given above, an improper handling or the incorrect use of the products delivered by us shall lead to the exclusion of any warranty claims. HEBU medical shall not be held liable for damage resulting from any such failure or misuse.

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HEBU medical GmbH