**HF 9501**

# HBS 100i

**Elektrochirurgie-Gerät / Electrosurgical Unit / Appareil  
d'électrochirurgie / Aparato de electrocirugía /  
Apparecchie di elettrochirurgia**

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Seite 3

## **Operating Manual**

Page 31

## **Mode d'emploi**

Page 59

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Página 87

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Página 115

**DEUTSCH**

**ENGLISH**

**FRANÇAIS**

**ESPAÑOL**

**ITALIANO**

## HF 9501

# Inhalt

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## HF 9501

# Contents

<b>1</b>	<b>Symbol descriptions</b> .....	<b>32</b>
<b>2</b>	<b>Introduction</b> .....	<b>33</b>
2.1	Monopolar operation .....	34
2.2	Bipolar operation .....	34
<b>3</b>	<b>Scope of delivery</b> .....	<b>34</b>
<b>4</b>	<b>Scope</b> .....	<b>35</b>
4.1	Inteded Use .....	35
4.2	Contraindications.....	35
<b>5</b>	<b>Structure and basic components of the system</b> .....	<b>36</b>
5.1	Main menu of HBS 100i .....	38
5.2	Overview of symbols for operating modes and settings .....	41
<b>6</b>	<b>Specifications</b> .....	<b>42</b>
<b>7</b>	<b>Accessories</b> .....	<b>43</b>
7.1	Accessories for the monopolar operating mode .....	43
	Neutral electrode.....	43
7.2	Accessories for the bipolar operating mode.....	44
7.3	Foot switch .....	44
<b>8</b>	<b>Commissioning the unit</b> .....	<b>45</b>
8.1	Before starting work .....	45
8.2	Switching on .....	45
8.3	Foot switch settings.....	46
8.4	Handles with buttons .....	46
8.5	Neutral electrode settings and attachment.....	46
	Application information for single-use neutral electrodes .....	47
	Application information for reusable neutral electrodes .....	47
	Areas for attachment of the neutral electrode.....	48
8.6	Setting the operating mode and power level.....	49
8.7	Saving the program and access.....	49
8.8	Switching off .....	49
8.9	Operating errors .....	49
<b>9</b>	<b>Guideline value for power settings</b> .....	<b>51</b>
<b>10</b>	<b>Protective measures and warnings</b> .....	<b>52</b>
<b>11</b>	<b>Review of the technical status</b> .....	<b>54</b>
11.1	Mechanical faults.....	54
11.2	Recurring safety inspections .....	54
11.3	Servicing and manufacturer address .....	55
<b>12</b>	<b>System maintenance work</b> .....	<b>55</b>
12.1	Cleaning .....	55
12.2	Fuses.....	55
12.3	Transport .....	55
12.4	Environmental directives .....	56
<b>13</b>	<b>Warranty card</b> .....	<b>57</b>

# 1 Symbol descriptions

Symbol	Definition
	CE-labelling
	Attention
	Manufacturer
	Lot-description
	Reference code
	Medical device / FDA Prescription device
	Medical device
	Non sterile
	Keep away from sunlight
	Dry storage required
	Temperature
	(Electronic) instruction for use
	No intervention on patients with cardiac pacemakers.
	WEEE-labelling
	Green dot

## 2 Introduction

An electrosurgical unit is a generator used to transform electronic energy into high-frequency alternating current (HF current).

If this current flows through human tissue, it creates heat which is used for cutting and for coagulation of the tissue.

The HBS 100i is a versatile and compact electrosurgical unit which complies to a high safety standard and unites economy with operative flexibility.

It offers the following application modes:

For monopolar surgical use

- Pure cut
- Blend cut
- Soft/Contact coagulation (no sparks/desiccate mode)
- Forced coagulation

For bipolar use:

- Bipolar cutting
- Bipolar Coagulation

Benefits of the HBS 100i:

- Simple, intuitive and safe operation
- Input through membrane keys
- One monopolar multifunction terminal for 3-pin standard or Martin plug
- One bipolar multifunction terminal for 3-pin standard or Martin plug
- Tone signal with scope for control of the volume and an alarm tone in case of faults without facility for switching off
- Adaptation (automatically) to different tissue impedances
- Neutral electrode monitoring
- Bipolar coagulation with 100W power
- Nine storage locations for individual program settings
- Power control in each operation mode, 1-25W in 1W steps for micro applications
- HEBU generators have an automatic output dosage by impedance with voltage regulation and arcing regulation and an automatic power adjustment with dynamic control and stop in case of anomalies.
- Worldwide application possible through connection scope from 100 to 260 volts, 50 to 60 Hz

## 2.1 Monopolar operation

In the monopolar mode, HF current is applied to the tissue using an active electrode. The cutting or coagulation effect is caused by a high concentration of current, i.e. a high current density over the small surface of the active electrode. This generates an increase in temperature, causing the water from the tissue immediately surrounding the active electrode to vaporize. Depending on the intensity of the HF current, this causes bleeding to stop or generates a cut in the tissue.

The HF current flows from the active electrode to the neutral electrode and is distributed over its large surface. This reduces the density of the current in this area and generally prevents unwanted thermal effects at the position of the neutral electrode. The current flows back to the electrosurgical unit through the neutral electrode.

## 2.2 Bipolar operation

In the bipolar mode, the current path is limited only to tissue located between the two electrodes of a bipolar instrument.

Use of the neutral electrode is not necessary. Bipolar application does not entail any risk of the current flowing through the patient's body. Consequently there is no danger of burns due to a poorly or incorrectly positioned neutral electrode.

Bipolar coagulation is therefore safer than the monopolar method and is particularly recommended for patients with pacemakers or for operations on organs with a small diameter.

## 3 Scope of delivery

HF 9501	HBS 100i
HF 9979-XX	Power cable
	Operating instructions in five languages
	Device manual in five languages
	Service manual in five languages

## 4 Scope

The use of electrosurgical equipment is predominantly in the following areas:

- General surgery
- Traumatology
- Cardiac and Thoracic Surgery
- Vascular Surgery
- Urology
- ENT
- Gynecology

### 4.1 Inteded Use

An electrosurgical unit is a generator that converts electronic energy into high frequency alternating current (RF) power. If this current flow through biological tissue, heat is generated there, which is used to cut and coagulate the tissue.

The device has monopolar and bipolar application modes.

In monopolar mode, RF current is applied to the tissue with an active electrode. The cutting or coagulation effect is caused by a high concentration of current, i. a high current density over the small area of the active electrode. As a result, an increase in temperature is generated and the water evaporates from the tissue immediately surrounding the active electrode. Depending on the intensity of the HF current, this leads to a standstill of bleeding or to a cut in the tissue. The HF current flows from the active electrode to the neutral electrode and is distributed on its large surface area. Thus, there reduces the density of the current and there are usually no unintended thermal effects at the position of the neutral electrode. Through the neutral electrode, the current flows back to the electrosurgery device.


In bipolar mode, the current path is limited only to tissue located between the two electrodes of a bipolar instrument. The use of the neutral electrode is not required.

Optionally, the HF units can be modularly expanded with an argon unit. The argon device is an argon unit designed for cutting and coagulation of biological tissue with high-frequency currents under argon protection.


### 4.2 Contraindications

HF surgical devices should not be used with simulators such as bone stimulators. It also discourages the use of electrosurgical devices when pacemakers, such as cardiac pacemakers, are in use.

## 5 Structure and basic components of the system

Front view of the HBS 100i

Rear view of the HBS 100i

The housing of the HBS 100i is made of metal and plastic. As no openings are available for ventilation, the unit can be cleaned easily (cf. 12.1)

**Located on the front of the unit are:**


- **(1)** Terminal for a neutral electrode
- Connections for
  - **(2)** one bipolar instrument: 2-pin standard plug or MARTIN plug
  - **(3)** one monopolar instrument: 3-pin standard or MARTIN plug
- **(4)** Main switch (on/off)

**Located on the back of the housing are:**



- **(5)** Terminal for a power supply cable (100-260 V, 50-60 Hz)
- **(6)** Earthing pins
- **(7)** Volume controller
- **(8)** Rating plate with manufacturer's information regarding serial number, mains voltage, output and fuses
- **(9)** Connection for foot switch of monopolar or bipolar mode (mono / bi)

## 5.1 Main menu of HBS 100i

### Monopolar


- (10) A yellow LED indicates active RF power during the operation in a mode of CUT
- (11) A blue LED indicates active RF power during the operation in a mode of COAG
- (12) Performance indication CUT
- (13) Performance indication COAG

- (14) Displays the selected mode of CUT
- (15) Button for pure cut
- (16) Button for blend cut
- (17) Button for reduce reduce power (CUT)
- (18) Button for increase power (CUT)
- (19) Displays the selected mode of COAG
- (20) Button for soft coagulation
- (21) Button for forced coagulation
- (22) Button for reduce power (Coag)
- (23) Button for increase power (Coag)


## Bipolar

- **(24)** Yellow LED indicates active RF current while operating in a BICUT mode
- **(25)** Performance indication CUT
- **(26)** Button for increasing power (BICUT)
- **(27)** Button for reducing power (BICUT)
- **(28)** Blue LED indicates active RF current while operating in a BICOAG mode
- **(29)** Performance indication BICOAG
- **(30)** Button for increasing power (BICOAG)
- **(31)** Button for reducing power (BICOAG)
- **(32)** Button BICOAG-Autostart


## Neutral electrode

- **(33)** Indication by red LED: Neutral electrode not plugged in, or not recognized
- **(34)** Indication by green LED: split neutral electrode detected
- **(35)** Indication by green LED: not split neutral electrode detected


## Foot switch

- (36) Button: setting of the foot switch
- (37) Indication by blue LED: footswitch set to monopolar coagulation
- (38) Indication by blue LED: footswitch set to bipolar

## Program

- (39) Displays the program number
- (40) Button to select the program

## 5.2 Overview of symbols for operating modes and settings

Button	Meaning
	Increase output
	Reduce output
	Pure Cut
	Blend Cut
	Soft/Contact coagulation (no sparks/desiccate mode)
	Forced Coagulation
	Bipolar Cut
	Bipolar Coagulation
	BICOAG Autostart
	Foot switch setting
	Program selection

Status	Meaning
	Undivided neutral electrode detected
	Divided neutral electrode detected
<b>ALARM</b>	Neutral electrode not detected or not connected
MONOPOLAR	The foot switch is set to <b>MONOPOLAR</b>
BIPOLAR	The foot switch is set to <b>BIPOLAR</b>
BIAUTOSTART	The foot switch is set to <b>BIAUTOSTART</b>

## 6 Specifications

<b>Mains connection</b>	
Supply voltage	100-260 Volt, 50 - 60 Hz.
Nominal frequency	500 kHz
<b>Output power in monopolar operation</b>	
Pure cut	100 W at 500 $\Omega$
Blend cut	100 W at 500 $\Omega$
Soft/Contact coagulation (no sparks/desiccate mode)	100 W at 50 $\Omega$
Forced coagulation	100 W at 500 $\Omega$
<b>Output power in bipolar operation</b>	
Bipolar cut	100 W at 300 $\Omega$
Bipolar Coagulation	100 W at 100 $\Omega$
BICOAG-AUTOSTART	Yes
<b>Safety measures</b>	
Type	CF
Protection class	I
Type Equipotential connection	Yes
Safety	In accordance with DIN EN ISO 60601-1
Electromagnetic compatibility	In accordance with DIN EN ISO 60601-1-2
Safety of high-frequency surgical units	In accordance with DIN EN ISO 60601-2-2
<b>Miscellaneous</b>	  
Operating temperature	+10° C to +40°C
Humidity during operation	30-75%
Air pressure during operating	700-1060hPa
Storage temperature	-10° C to +60°C
Storage humidity	10-85%
Air pressure in storage	500-1060hPa
Dimensions (L/W/H)	310 x 310 x 145 mm
Weight	6,9 kg

## 7 Accessories

 **Apply for a copy of our extensive electro-surgery accessories catalogue today!**

### 7.1 Accessories for the monopolar operating mode


Handles for monopolar electrodes can be connected to the unit using a 2-pin standard or MARTIN plug or using the 4 mm plug (banana plug) with only one contact.

3-pin standard plug

MARTIN monopolar plug

Monopolar 4mm 1-pin plug (banana plug)

### Neutral electrode

Valleylab Connecting plug for the neutral electrode

The HBS 100i reliably detects neutral electrodes with undivided and divided contact surface. A comparison of the electrode detected by the machine and the used electrode provides a direct indication of correct adhesion of the neutral electrode. Both single-use and reusable neutral electrodes can be used.

 **HEBUMedical recommends the use of the single-use safety neutral electrode with double divided contact surface (HF9564H) for maximum safety!**

## 7.2 Accessories for the bipolar operating mode

Bipolar instruments are connected to the electrosurgical unit using a 2-pin standard or Martin plug.

MARTIN bipolar plug

Connecting plug for bipolar instruments

## 7.3 Foot switch

Foot switches with 6-pin connector can be used:

6-pin connector

**HF 9533-02:**  
Foot switch,  
cable length 4.5m,  
with 6-pin connector

**HF 9535-01:**  
Double foot switch,  
cable length 3m,  
with 6-pin connector

## 8 Commissioning the unit

### 8.1 Before starting work

Before starting surgery, the HBS 100i electro-surgical unit is connected via the supply cable to the mains socket with PE conductor. The current cable may only be attached or removed when the device is switched off. The connecting socket for the power cable is located at the back of the housing.

Connect foot switch, the terminal is located at the back of the housing.

Handles and electrodes are connected correspondingly at the front of the housing to "MONOPOLAR" or "BIPOLAR".

**Before first using the electro-surgical unit, we recommend familiarizing employees with the operating characteristics of the different operating modes and power settings by carrying out practice exercises using fresh meat (beef, pork).**

### 8.2 Switching on

When using an electro-surgical unit which produces HF current, always be aware of two fundamental principles:

- **Current flows through everything located between the active and the neutral electrode (monopolar)**

The system is switched on by actuating the main switch (ON/OFF). However, high-frequency current only flows once a button is activated at the handle or foot switch.

Consequently, accessories can also be connected to the unit when it is switched on.

However, extreme care should always be taken when doing so!

**Always ensure that the system is not inadvertently activated during installation by pressing the foot switch or the button in the handle.**

### 8.3 Foot switch settings

After connecting the foot switch at the back, the required mode must be set at the front. To set the mode, press the button of the foot switch, until the LED is lit next to the desired setting.

Double and single footswitches can be used:

- With the double foot switch you can activate the cutting or coagulation power (if set) by pressing the pedal
- With the single footswitch you can only activate the set mode. It must be selected by pressing the foot switch symbol.

### 8.4 Handles with buttons

Activating the handle with buttons takes place independently of the foot switch settings. These cannot be activated using the foot switch.

### 8.5 Neutral electrode settings and attachment

In the monopolar operating mode, a neutral electrode is required.

After attaching the neutral electrode to the patient, the HF unit recognizes the neutral electrode automatically (cf. illustration on the left and chapter 5.2 Overview of symbols for operating modes and settings). The HBS 100 recognizes neutral electrodes with divided and undivided contact surfaces.

**If the displayed neutral electrode type does not agree with the used neutral electrode type, this is indicative of poor adhesion of the neutral electrode on the patient, a defective neutral electrode or a defective cable!**

**Please note the following information when attaching the neutral electrode**

## Application information for single-use neutral electrodes

- When using single-use products, always pay attention to the use-by date. Do not use electrodes which are over this date!
- Neutral electrodes designed for single use must not be used a second time!
- The electrode must not be attached to surfaces with a high amount of hair growth. If necessary, shave the relevant area before attaching.
- The neutral electrode must be attached in such a way that its entire surface adheres to the patient's body.
- After every change of the patient's position, check that the neutral electrode is adhering correctly and also check the cable routing.
- As soon as the electrode has been attached, it must not be pulled off and attached again. Should the electrode have to be repositioned, use a new neutral electrode.
- The neutral electrode must not come into contact with fluids and must not be coiled round.
- Do not apply any additional electrically conductive gel to the neutral electrode.
- When performing surgery on small children, special electrodes suitable for children or babies must be used (see manufacturer's instructions)
- Never under any circumstances use damaged products!
- When removing the neutral electrode, never pull on the cable!
- When pulling off single-use electrodes, ensure that the patient's skin is not damaged. Avoid any abrupt movements.


## Application information for reusable neutral electrodes

- The neutral electrode must be wiped with disinfectant before repeated use. To disinfect reusable rubber neutral electrodes, we recommend the products Incidin perfekt, Minutil and Incidur F from Ecolab.
- The neutral electrode must be attached in such a way that its entire surface adheres to the patient's body. A rubber band provided for this purpose can be used to ensure better adhesion.
- When attaching the neutral electrode, ensure that the long side is pointing towards the operation field.
- After each change in the patient's position, check the correct adhesion of the electrode and the connecting cable!
- The electrode must not be attached to surfaces with a high amount of hair growth. If necessary, shave the relevant area before attaching.
- The electrode must not be wet or coiled round.
- Fluids between the patient's skin and the neutral electrode must be avoided.
- Do not apply any additional electrically conductive gel to the neutral electrode.
- When performing surgery on small children, special electrodes suitable for children or babies must be used (see manufacturer's instructions)
- When removing the neutral electrode, never pull on the cable!
- Never under any circumstances carry out repairs on the neutral electrode!


Please note that rubber electrodes lose their conductive capability if active substances are eroded due to cleaning of the material. This type of electrode considerably increases the risk of burns. Consequently ensure that not only the unit but also the reusable neutral electrode is included in regular maintenance inspections.

### Areas for attachment of the neutral electrode


- Attach the neutral electrode to clean, dry skin.
- The electrode must be attached close to the operation area, but with a minimum direct distance from the operating field of 20 cm.
- Attach over tissue with a good supply of blood without dips or curves in the skin, for instance on the upper arm or thigh. (cf. picture below)
- Do not attach the neutral electrode to projecting surfaces.
- Do not attach above bones, scars, cuts or scratches.
- The electrode must not be attached to surfaces with a high amount of hair growth. If necessary, shave the relevant area before attaching.
- Do not attach above particularly adipose areas such as the abdomen or buttocks.
- The neutral electrode must not be attached above implants.
- Electrodes of monitoring systems must be placed at a distance of at least 20 cm from the operation area and from the position of the attached neutral electrode
- Correct attachment of the neutral electrode in adult patients

- Correct attachment of the neutral electrode in children

- Attachment in patients with pacemakers  
**Note! Consult the treating cardiologist!**

**RIGHT**

**WRONG!**

## 8.6 Setting the operating mode and power level

HBS 100i has four different monopolar modes and one for bipolar applications (for details, look chapter 5.2). The selection of all modes takes place according to the same principle: the button with the picture of the operating mode must be pressed. The LED next to the button indicates the selected mode. Use the arrow buttons in order to set the required power setting.


## 8.7 Saving the program and access

Any change of settings of modes or power is automatically saved in the program in which you are. To save individual settings first you have to select a different program. Therefore, the button with the program symbol must be pressed. Even after turning off and on again the device will keep these values.

## 8.8 Switching off

After terminating the operation, switch off the system at the main switch and pull out the plug. After switching off the system, all electrodes and forceps must be disconnected from the cable and the electrode cable from the unit.

## 8.9 Operating errors

The following error messages can occur:

Error	No.	Cause of error	Remedying errors
Err 1 Total Failure	1	"Power timeout error"	Check the power terminal and power supply unit, exchange if applicable. Check the controller and the power supply unit, exchange if applicable.
	10	"Power Error"	The output measured by the monitoring module is 20% greater than the set output. Check the controller, the interface and the power supply unit, exchange if applicable.
Err 2 Modules Failure	2	"NE electrode error"	Check the terminal, and the plate on the patient, where applicable exchange NE and/or connecting cable.
	6	Mono - Cut hand /FS1 Error"	Check the CUT terminals, check the settings, exchange accessories if necessary. If error persists, inform specialist personnel.
	7	"Mono – Coag (Bi) hand /FS1 Error"	Check the COAG terminals, check the settings, exchange accessories if necessary. If error persists, inform specialist personnel.
Err 03 PS Failure	4	"Power Supply 1 Error"	The necessary voltage from the block for regulating the power voltage for the output generator in idle running is missing. Check the power terminal and power supply unit, exchange if applicable.
	5	"Power Supply 2 Error"	The necessary voltage from the block is missing for correction of the output factor. Check the power terminal and power supply unit, exchange if applicable
Err 5 CPU Failure	11	"EEROR in CPU ROM"	Falsche Check-Summe des ROM – Speichers vom Steuerprozessor. Kontroller überprüfen, ggf. auswechseln.
	12	"EEROR in CPU RAM"	Incorrect checksum of the RAM memory from control processor. Incorrect checksum of the ROM memory from control processor. Check controller, exchange if applicable.
Err 6 Keys Failure	0-17	„Button stick“	Check your settings and the buttons.

If one or more numerical codes are displayed, or if the device is indicating any other recognizable defects, please refer to [service@hebumedical.de](mailto:service@hebumedical.de)

■ **HF current does not flow when the foot switch is actuated. The system only reacts with a tone signal and the display changes colour.**

Please check:

- ➔ That all terminals have been correctly connected
- ➔ The settings of the foot switch
- ➔ Whether the neutral electrode has been recognized

- **The neutral electrode is not detected or not correctly**  
Check all terminals. If the neutral electrode is still not recognized, exchange it where appropriate.
- **It is not possible to switch the machine on or off.**  
Check whether the power cable has been correctly connected and that the socket is intact. If the error persists, inform the specialist personnel.
- **The system works perfectly, but the settings cannot be saved.**  
Inform the specialist personnel.

## 9 Guideline value for power settings

When performing the power settings, the user should generally take into account the cutting speed, geometry of the electrode and properties of the tissue.

		Recommended settings
Monopolar cutting	Smooth cuts with thin needle electrode	5 - 40 W
	Cutting with knife electrode, lancet or wire loop	50 - 100 W
Monopolar coagulation	Coagulation of large surfaces	40 - 100 W
	Coagulation of small surfaces	20 - 50 W
Bipolar cut	Cut	40 – 60 W
Bipolar application	Coagulation	30-60 W

The recommended settings serve as guideline values only and must be adjusted to the relevant situation. In the event of unfavourable conditions, these values may deviate!

## 10 Protective measures and warnings

Please observe also the instructions provided on working with the neutral electrode in chapter 8.5, the operating instructions with the accessories and the applicable safety measures!

**Before any procedure, check the cables, accessories (handles, electrodes, forceps) and foot switches for any signs of visible irregularity. Do not use cables or instruments with broken or defective insulation. These can result in burns if they come into contact with the skin!**

The unit should only be operated at a mains socket with protective earth contact. When performing electrosurgical procedures, the risk of burns must be kept as small as possible. For this, the following safety measures must be observed:

All plug-in connections may only be carried out with compatible accessories.

Use the greatest care when working with the neutral conductor!

Protect the patient from contact with earthed metallic elements (OP table, infusion stand etc.) during use.

Avoid contact between supply cables and patients, other cables or instruments! Ensure that there are no loops in the cables of handles and electrodes or fixture to the OP table.

Avoid small-area skin-to-skin contact in the patient (e.g. a finger making contact with the thigh).

Avoid contact between HF instruments and non-insulated instruments.

Prevent the penetration of fluids between the patient's body and the neutral electrode.

The output power should not be set higher than is absolutely necessary for the procedure.

Do not touch the patient during activation of the coagulation or cutting current. Should this be necessary, a piece of fabric can be used for insulation purposes.

In patients connected to monitoring devices (ECG), the monitoring electrodes must be attached at the greatest possible distance to the application areas of the electrosurgical electrodes.

The use of highly flammable narcotic substances and oxidizing gases (N<sub>2</sub>O or O<sub>2</sub>) should be avoided. Should this be necessary, take particular precautions.

Exclusively non-flammable disinfectants must be used. If alcohol-based solutions are required, wait before performing surgery until these have completely evaporated. Bear in mind that disinfectants can flow down and collect in recesses of the body. A spark from the active electrode can cause combustion of the fluid in normal use.

Sparks from the active electrode can set bandages (textiles), metabolic gases and combustible fluids alight.

In the case of procedures involving high-frequency current applied to body parts with a minimal cross-section, the bipolar technique should be used in order to prevent the risk of unintentional coagulation.

In the case of procedures on body parts with a close link to the remainder of the body, the bipolar technique should be used for safety reasons.

The frequently used indirect coagulation mode should only be used with well insulated forceps. We recommend forceps with specially insulated grip areas (ask your HEBUmedical contact).

Surgical gloves do not afford the surgeon sufficient protection from burns!

During surgery, pacemakers are at risk of destruction. We recommend consulting the treating cardiologist and/or an authorized representative of the pacemaker manufacturer. As a rule, in such cases bipolar application technology is advisable. Electrosurgical procedures are prohibited on outpatients with pacemakers.

A marked drop in power at a normal setting can be an indication of incorrect application of the neutral electrode.

Soiled electrodes can bring about a drop in the quality of the unit. This applies particularly to gentle and bipolar coagulation. To prevent the tissue adhering to the active electrode (monopolar or bipolar), it is advisable to wait for a moment after coagulation of the tissue and to interrupt the current flow to allow the active tip to cool before separating the electrode from the tissue.

High-frequency creep currents can cause burns in areas distant from the application site through conductive connections.

To remove the cable from the unit, always hold the plug! Never pull on the cable, as this can result in damage to the cable insulators, burns on the skin or also fire.

Instruments which are temporarily not required during the procedure must be placed away from the patient on the instrument table or a similar surface.

Reusable electrosurgical instruments which are connected to the unit are generally not delivered in a sterile condition!

If the products are used on patients with transmissible spongiform encephalopathy or HIV infection, we decline any responsibility for their reuse.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## ■ Electromagnetic safety

In the development of electrosurgical units, HEBUmedical attaches particular importance to the stringent electromagnetic emission regulations. Consequently, the solutions we supply have a correspondingly minimal emission level for the required output. Measurements conform the high level of electromagnetic safety offered by HEBUmedical electrosurgical units.

Under typical working conditions based on an 8-hour day, a field occurs at a distance of 5 to 15 cm from the cables. At a distance of 20-40 cm, the value drops well below the prescribed limit.

Electromagnetic fields occur primarily around the cable.  
The electrosurgical unit itself does not constitute a significant emitting element.

## 11 Review of the technical status

Before any procedure, check the cables, accessories (handles, electrodes, forceps) and foot switches for any signs of visible irregularity. Do not use cables or instruments with broken or defective insulation. These can result in burns if they come into contact with the skin!

### 11.1 Mechanical faults

If damage occurs to switches, connecting sockets, houses or the touchscreen, or if the unit is dropped on the floor, please contact your authorized dealer for advice on further use.  
The dealer can arrange for a detailed technical check to be performed.

The manufacturer accepts no liability for calibration work or repairs to the unit carried out by the user.

### 11.2 Recurring safety inspections

A safety inspection of the electrosurgical system is required once a year. The manufacturer only permits the use of units for which an up-to-date inspection by suitably trained personnel can be verified.

All information relating to the inspections (date, inspecting staff member, result etc.) must be recorded in the unit's log book.

## 11.3 Servicing and manufacturer address

Should you require the instructions for use in paper form, please use the contact details below. The instructions for use in paper form will be made available to you within seven calendar days of receipt of the request.

Alternatively, you can print out the electronic instructions for use yourself.

HEBUmedical GmbH  
Badstraße 8  
78532 Tuttlingen / Germany

Tel. +49 7461 94 71 - 0  
Fax +49 7461 94 71 - 22  
eMail: [service@HEBUmedical.de](mailto:service@HEBUmedical.de)  
Web: [www.HEBUmedical.de](http://www.HEBUmedical.de)

## 12 System maintenance work

### 12.1 Cleaning

The electrosurgical unit can be wiped using a cloth dampened with standard disinfectant, and so kept clean using the simplest possible method. Avoid allowing cleaning fluid and disinfectant to penetrate into the unit.

### 12.2 Fuses

The fuse sockets are located on the back of the unit. The unit works with two fuses of 6.3 A each.

Exchanging the fuses:

- Disconnect the unit from the mains
- Release the fuse from the retainer using a flat blade screwdriver
- Exchange the fuse

After an exchange of fuses, the unit must be checked by suitably qualified personnel.

### 12.3 Transport

Prevent mechanical damage and moisture.

If the unit has been stored for a long period in a cold environment, wait before removing the protective packaging until the unit has reached room temperature.

During transportation, the standard safety measures additionally apply.

Only ever transport the unit in its original packaging.

## 12.4 Environmental directives

From the implementation date of the European directive 2002/96/EU as national legislation, the following regulations apply:

- Electrical and electronic units must not be disposed of with normal household waste.
- The consumer is required by law to return electrical and electronic devices at the end of their service life to designated public collection points or to the place of purchase.

Disposal details are set out by the applicable local legislation in the relevant country. The symbol on the product, the instructions for use or the packaging indicates the relevant applicable legislation. By sending end-of-life devices for reuse, material recycling or other forms of utilization you are making a major contribution to environmental protection.

Please also send the unit packaging for environmentally responsible recycling after the end of the product's life.

## 13 Warranty card

# Warranty card

We grant on all electrical medical products

We grant a warranty in accordance with the statutory and country-specific regulations on all electrical medical products (verification required in the form of an invoice or delivery note). The warranty period is two years. In the event of material or production errors, we provide free repair of the unit. Any attempt to carry out independent repairs will render the warranty cover null and void.

Damage caused by incorrect handling, overloading or normal wear and tear is not covered under the terms of the warranty.

In the event of a warranty claim and resulting repair, please enclose a copy of the invoice or delivery note (or batch number) with the instrument.

**Manufacturer and Aftersales Service Address:**

HEBUmedical GmbH  
Badstraße 8  
78532 Tuttlingen  
Tel. +49 7461 9471 - -0  
Fax +49 7461 9471 -22  
[service@HEBUmedical.de](mailto:service@HEBUmedical.de)

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