HF 9518-00

HEBU HFIT[®]-T

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HF 9518-00

Content

1	Symbol descriptions.....	37
2	Introduction.....	38
2.1	Functional principle.....	38
3	Scope of delivery	39
4	Scope.....	39
4.1	Intended use.....	39
4.2	Contraindications	39
5	Structure and basic components of the system.....	40
5.1	Main menu.....	41
6	Specifications	42
7	Commissioning.....	43
7.1	Installation	43
7.2	Switching on / off	43
7.3	Insertion of the syringe.....	43
7.4	Creating the connection between HF generator and perfusion pump via RS 232.....	44
7.5	Connection of the active electrode, neutral electrode and foot switch	44
7.6	Adjusting the treatment parameters	45
7.7	Creating and accessing the program	47
8	Application information.....	49
8.1	Attaching the neutral electrode	49
8.2	Bipolar thermotherapy.....	52
8.3	Start of treatment	53
8.3.1	Activation of the flow rate.....	53
8.3.2	Activation of the HF power	53
8.4	End of the treatment	54
8.4.1	Switching off the flow rate and the HF power	54
8.5	Synchronization mode	54
8.6	Automatic mode.....	55
8.7	Real-time monitoring.....	56

8.8 Setting recommendations56

9 Accessories.....57

10 Operating errors58

11 Performance diagram.....60

12 Protective measures and warnings.....61

13 Review of the technical status63

13.1 Mechanical faults63

13.2 Recurring safety inspections63

13.3 Servicing and manufacturer address.....63

14 System maintenance work64

14.1 Cleaning64

14.2 Fuses.....64

14.3 Transport64

14.4 Environmental directives64

15 Warranty card65

1 Symbol descriptions


Symbol	Definition
	CE-labelling
	Attention
	Manufacturer
	Lot-description
	Reference code
	Medical device / FDA Prescription device
	Medical device
	Non sterile
	Keep away from sunlight
	Dry storage required
	Temperature
	(Electronic) instruction for use
	No intervention on patients with cardiac pacemakers.
	WEEE-labelling
	Green dot

2 Introduction

Thank you for purchasing a product of HEBUmedical GmbH. The electrosurgical system has been designed and manufactured according to the latest development in science and technology. Please read these instructions carefully and keep them for future reference. Please note the safety and maintenance instructions to achieve good long-term results with this device.

2.1 Functional principle

While observing through imaging systems (ultrasound, CT), the high-frequency electric current will be applied via a needle electrode to the tissue thereby the tissue will be heated and denatured. Due to the use of saline solution heat accumulation and exsiccation of the tissue are avoided furthermore the electrical conduction and the temperature around the tip of the needle electrode are stabilized.

It is distinguished between two methods of HF induced thermotherapy:

The monopolar thermotherapy: In this case the HF current flows between the needle electrode and neutral electrode which is attached to the body of the patient. In principle, the current can flow through the entire body.

The bipolar thermotherapy: The induces HF current flows between two needle electrodes, which must be placed at a distance of about 6mm from each other.

Both methods can be realized with the HFiT®-T system. Additionally there is the possibility to monitor the temperature at the needle tip by a thermal probe.

3 Scope of delivery

In the delivery condition, the system consists of the following components:

HF 9518	HF Generator HEBU HFIT®-T
HF 9518-01	Perfusion pump
HF 9537-02	Double foot switch
HF 9529	Trolley
HF 9979-XX	Mains cable
HF 9418-02	RS 232 data cable
	Power connection cable
	Operating Manual
	Logbook
	Service Manual

Please check immediately the system for any transportation damages and completeness of all components after receiving the equipment. In case of a return, please use the original packaging.

4 Scope

The electrosurgical system HFIT®-T is used for percutaneous minimally-invasive treatment of inoperable liver cancer (HCC) and other pathogenic structures

- Oncology
- Urology
- General Surgery


4.1 Intended use

HFIT®-T needle electrodes are used as applicators in the radiofrequency thermal ablation of inoperable tumors of the liver and other pathogenic bio-tissues. High-frequency current is introduced into the tissue via the HFIT®-T needle electrode and perfused with physiological saline solution. HFIT®-T needle electrodes are to be used exclusively with the HEBU HFIT®-T electrosurgical system. It is a disposable product.

4.2 Contraindications


The HF monopolar electrodes must not be used in conjunction with the central nervous or circulatory system. Patients with pacemakers and in vivo cardiac defibrillators may have restrictions on use. HF monopolar electrodes should not be used in the immediate vicinity of bone, tooth root and metal.

5 Structure and basic components of the system

Front view of: HEBU HFiT®-T


- **(1)** Terminal for needle electrode
- **(2)** Connection for temperature sensor
- **(3)** Terminal for neutral electrode
- **(4)** Main switch (on/off)
- **(5)** LCD touchscreen display

Rear view HEBU HFiT®-T

- **(6)** Earthing pins
- **(7)** USB-connection for service
- **(8)** Connection RS 232
- **(9)** Volume controller
- **(10)** Connection for foot switch
- **(11)** Current terminal for perfusion pump
- **(12)** Main terminal / fuse

5.1 Main menu

- (13) Display patient identification
- (14) Display power / Menu power setting
- (15) Display flow rate / Menu flow rate
- (16) Display time / Menu time setting
- (17) Menu program selection
- (18) Menu program name
- (19) Lock: settings lock / unlock
- (20) Symbol synchronization mode active
- (21) Symbol automatic mode active
- (22) HF-Energy activate / deactivate
- (23) Perfusion activate / deactivate
- (24) Connection with perfusion pump establish / interrupt (RS 232)
- (25) Display temperature / Menu temperature setting
- (26) Display „introduced Energy “ in watt seconds / Resetting the display
- (27) Recognition and status of the neutral electrode
- (28) Display tissue impedance

The arrangement of terminals and controls of the perfusion pump please refer from the corresponding operating manual.

6 Specifications

HF Generator HEBU HFIT	
Supply voltage	100 – 260 V ac, 50/60 Hz
Power consumption	174 VA (100 V ac) / 165 VA (230 V ac)
Maximum HF-power output	100 W (150 Ohm) 345 V _{ss} / 122 V _{eff}
RF nominal frequency	500 kHz
Dimensions (L/B/H)	405 x 395 x 145 mm
Weight	8,1 kg
Safety	
Type	CF
Protection class	I
Equipotential connection	Yes
Safety	In accordance with DIN EN ISO 60601-1
Electromagnetic compatibility	In accordance with DIN EN ISO 60601-1-2
Safety of high-frequency surgical units	In accordance with DIN EN ISO 60601-2-2
Perfusion pump	
Power consumption	power consumption 23 VA
Delivery rate	0,1 – 200 ml/h for 50 ml syringe 0,1 – 120 ml/h for 20 ml syringe
Dimensions (L/B/H)	330 x 155 x 120 mm
Weight	2,4 kg
Safety	
Type	CF
Protection class	II
Equipotential connection	No
Safety	In accordance with DIN EN ISO 60601-1
Miscellaneous   	
Operating temperature	+10°C to +40°C
Humidity during operation	30-75%
Air pressure during operating	700-1060hPa
Storage temperature	-10°C to +60°C
Storage humidity	10-85% nicht kondensierend
Air pressure in storage	500-1060hPa

7 Commissioning

7.1 Installation

- The HF device should be placed on an inclined surface with a maximum incline of 10°. The surface should be at least equivalent to the base of the HF device.

It is recommended to place the HF device on the trolley HF 9529. This ensures the necessary stability and offers more mobility.

- The air intakes on the bottom of the HF device must remain free to ensure a good air circulation.

7.2 Switching on / off


Make sure that the main switch (4) is turned off. Position "0"

- Connect the power cable to the main power terminal (12) of the HF device and connect it to a mains socket with protective earth contact.
- Connect the power extension cable to the appropriate terminal of the perfusion pump and to the power supply terminal (13) of the HF device.
- Connect the RS232 to the perfusion pump in accordance with the operating manual and connect it to the RS232 terminal (8) of the HF device.
- Press the main switch (4) of the HF device.
- Start the perfusion pump in accordance with the operating manual.

7.3 Insertion of the syringe

- Connect the syringe (in coordination with the syringe selection) as usual to a Luer-Lock perfusions line.
- Now put the syringe into the syringe shaft. The syringe wing must lie correctly in the syringe hollow.
- Turn the syringe down holder into a closed position and secure the syringe plunger with the drive lock of the syringe pusher.

7.4 Creating the connection between HF generator and perfusion pump via RS 232

- Press the button “Perfusor“(24) in the main menu.

As soon as the HF device via RS232 is connected to the perfusion pump, on the screen of the perfusion pump flashes alternately “Stop” and the set flow rate.

The shown flow rate of the perfusion pump must match with the set flow rate of the HF device.

7.5 Connection of the active electrode, neutral electrode and foot switch

- Connect the plug of the needle electrode to the terminal (1) on the front of the HF device.
- Connect the 6-pin plug of the foot switch to the connection (10) on the back of the HF device. Tighten the sleeve by hand.
- Connect the plug of the neutral electrode to the terminal (3). Depending on the type and the status of the electrode the following statuses are indicated (27).

Statuses of the neutral electrode	
Status	Meaning
	Undivided neutral electrode detected
	Divided neutral electrode detected
 alarm	Neutral electrode not detected or not connected

7.6 Adjusting the treatment parameters


 For safety reasons, the activation of all control panels happens with a time delay.
For activation, the buttons must be pressed for about 0,5 s.

Power:


- Actuate the power display (14) in the main menu to reach the power settings menu.

- Select the required power level with the arrow keys (max. 100Watt)

- Actuate the “Exit” key to save the selection and return to the main menu.


 The operation of the keys „Sync P-Z-F“ and „Auto P-F“ is described in Chapters 7.5 and 7.6.

Flow rate:


- Actuate the flow rate display (15) in the main menu to reach the flow rate menu.

- Select the required flow rate with the arrow keys.

- Actuate the “Exit” key to save the selection and return to the main menu.


Treatment time:


- Actuate the time display (16) in the main menu to reach the time settings menu.

- Select the required treatment time with the arrow keys.

- Actuate the “Exit” key to save the selection and return to the main menu.


Treatment temperature:


- Actuate the temperature display (25) in the main menu to reach the temperature settings menu.

- Select the required treatment temperature with the arrow keys.

- Actuate the “Exit” key to save the selection and return to the main menu.

If the treatment is without a temperature sensor, the temperature must be set to "OFF".

7.7 Creating and accessing the program

The electrosurgical system offers the possibility to create and save nine different programs with individual treatment parameters.

Accessing the program:

- The active program with program number and program name is always displayed in the main menu.

Example: PROGRAM1 DR MULLER

To choose another program, actuate the program key in the main menu and it appears a program list.


- To select the required program actuate the key “Program name”.
- The check mark under the lock icon indicates that the particular program is secured.
- The column of the active program is highlighted in gray.


Nr	Program Name	Lock Icon
1	DR. MULLER	✓
2	Name 2	
3	Name 3	
4	Name 4	
5	Name 5	
6	Name 6	
7	Name 7	
8	Name 8	

Creating the program:


- To create a new program, actuate the key “Program” in the main menu.


- Select a program number by directly actuating the program name, e.g. no. 4.

- Unlock the selected program in the main menu by actuating the lock symbol (19).


- Then press the program name key (18).

- Now enter your desired name for the program. Actuate the “Enter” key to save the program name and return to the main menu.

- Actuate the “Esc” key to return to the main menu without saving.

- Now configure the parameters as described in chapter 4.5 and lock the program by pressing the lock symbol (19).


Changing the program:

If you want to change an existing program, the program must be unlocked. Subsequently, the corresponding change can be made. After changing, the program must be locked again.

8 Application information

8.1 Attaching the neutral electrode

In the monopolar operating mode, a neutral electrode is required.

After attaching the neutral electrode to the patient, the HF unit recognizes the neutral electrode automatically. The HEBU HFIT®-T recognizes neutral electrodes with divided and undivided contact surfaces.

As far as possible, neutral electrodes with divided contact surfaces are preferable.

If the displayed neutral electrode type does not agree with the used neutral electrode type, this is indicative of poor adhesion of the neutral electrode on the patient, a defective neutral electrode or a defective cable!

Please note the following information when attaching the neutral electrode

Application information for single-use neutral electrodes

- When using single-use products, always pay attention to the use-by date. Do not use electrodes which are over this date!
- Neutral electrodes designed for single use must not be used a second time!
- The electrode must not be attached to surfaces with a high amount of hair growth. If necessary, shave the relevant area before attaching.
- The neutral electrode must be attached in such a way that its entire surface adheres to the patient's body.
- After every change of the patient's position, check that the neutral electrode is adhering correctly and also check the cable routing.
- As soon as the electrode has been attached, it must not be pulled off and attached again. Should the electrode have to be repositioned, use a new neutral electrode.
- The neutral electrode must not come into contact with fluids and must not be coiled round.
- Do not apply any additional electrically conductive gel to the neutral electrode.
- When performing surgery on small children, special electrodes suitable for children or babies must be used (see manufacturer's instructions)
- Never under any circumstances use damaged products!
- When removing the neutral electrode, never pull on the cable!
- When pulling off single-use electrodes, ensure that the patient's skin is not damaged. Avoid any abrupt movements.


Application information for reusable neutral electrodes

- The neutral electrode must be wiped with disinfectant before repeated use. To disinfect reusable rubber neutral electrodes, we recommend the products Incidin perfekt, Minutil and Incidur F from Ecolab.
- The neutral electrode must be attached in such a way that its entire surface adheres to the patient's body. A rubber band provided for this purpose can be used to ensure better adhesion.
- When attaching the neutral electrode, ensure that the long side is pointing towards the operation field.
- After each change in the patient's position, check the correct adhesion of the electrode and the connecting cable!
- The electrode must not be attached to surfaces with a high amount of hair growth. If necessary, shave the relevant area before attaching.
- The electrode must not be wet or coiled round.
- Fluids between the patient's skin and the neutral electrode must be avoided.
- Do not apply any additional electrically conductive gel to the neutral electrode.
- When performing surgery on small children, special electrodes suitable for children or babies must be used (see manufacturer's instructions)
- When removing the neutral electrode, never pull on the cable!
- Never under any circumstances carry out repairs on the neutral electrode!


Please note that rubber electrodes lose their conductive capability if active substances are eroded due to cleaning of the material. This type of electrode considerably increases the risk of burns. Consequently ensure that not only the unit but also the reusable neutral electrode is included in regular maintenance inspections.

Areas for attachment of the neutral electrode


- Attach the neutral electrode to clean, dry skin.
- The electrode must be attached close to the operation area, but with a minimum direct distance from the operating field of 20 cm.
- Attach over tissue with a good supply of blood without dips or curves in the skin, for instance on the upper arm or thigh. (cf. picture below)
- Do not attach the neutral electrode to projecting surfaces.
- Do not attach above bones, scars, cuts or scratches.
- The electrode must not be attached to surfaces with a high amount of hair growth. If necessary, shave the relevant area before attaching.
- Do not attach above particularly adipose areas such as the abdomen or buttocks.
- The neutral electrode must not be attached above implants.
- Electrodes of monitoring systems must be placed at a distance of at least 20 cm from the operation area and from the position of the attached neutral electrode
- Correct attachment of the neutral electrode in adult patients

- Correct attachment of the neutral electrode in children

- Attachment in patients with pacemakers
Note! Consult the treating cardiologist!

RIGHT

WRONG

8.2 Bipolar thermotherapy

During the bipolar thermotherapy, the following point should be observed:

Connection:

For the bipolar use a special bipolar cable is required (see Chapter 7).


On the one hand, the cable is connected to the needle electrode terminal (1) and on the other hand to the neutral electrode terminal (3).

On the end of the cable two perfusion needles are connected.


Pump:

Furthermore an additional perfusion pump is required. This should have similar features as the HFIT-T®-Pump (pumping rate, syringe size, etc.).

Connect the outgoing hose from the external pump to the perfusion needle, which is connected to the neutral electrode connection.

Schema: Monopolar

Schema: Bipolar

Ensure, during bipolar application that both pumps are set to the same pumping rate.

Please note that the HF current during bipolar application has to travel a shorter distance and requires less power to overcome the tissue resistance.


8.3 Start of treatment

Activate the perfusion before the incision so that the needle electrode is vented. Place the needle always under perfusion.

8.3.1 Activation of the flow rate

- Actuate the flow rate key (23) in the main menu.
- A window will appear with the prompt to confirm the flow rate.
- Confirm by pressing the button "YES"
- If you actuated the key "flow rate" by mistake, you can acknowledge with the "NO" button.


Activation of the flow rate via foot switch

If a foot switch is connected to the HF surgery device, the flow rate can be activated by pressing the green pedal. A confirmation is not required.

8.3.2 Activation of the HF power

- Actuate the power button (22) in the main menu.
- A window will appear with the prompt to confirm the power.
- Confirm by pressing the button "YES"
- If you actuated the key "power" by mistake, you can acknowledge with the "NO" button.
- After activating the HF power, the defined treatment time runs down and it plays a soothing melody.

Activation of the HF power via foot switch

If a foot switch is connected to the HF surgery device, the flow rate can be activated by pressing the blue pedal. A confirmation is not required.

8.4 End of the treatment

8.4.1 Switching off the flow rate and the HF power

When the defined treatment time is over, the HF power will be switched off automatically. Meanwhile the flow rate continues.

Remove the needle electrode under perfusion and switch off the flow rate by pressing the green pedal on the foot switch or the flow rate button (23) on the device.

By actuating the blue pedal on the foot switch or the “power” button on the device the HF power can be switched off manually.

It is recommended to remove the active electrode already a few seconds before the end of the treatment period. Thereby, the working channel is made desolated and the bleeding will be stopped. The thus triggered alarm is acknowledged by pressing the blue pedal on the footswitch or on the "power button" on the HF device.

8.5 Synchronization mode

After switching on the device the "synchronization mode" is always active. This is visualized by the symbol (20).

In the synchronization mode, the parameters power, tissue impedance and flow rate are continuously measured and adjusted during the treatment.

If the tissue impedance exceeds the value of $700 \Omega \pm 50 \Omega$, then the power will be reduced for a period of 3 s to 2 W and the flow rate will be increased to the maximum rate of 200ml/h for a period of 3s.

Thus, the heat propagation in the tissue is optimized and the exsiccation of the tissue is reduced.


It is possible to deactivate the “synchronization mode”.

You must proceed as follows:


Actuate the button “Power setting” (14) in the main menu.

The menu “Power setting” appears.

- Actuate the key “Sync P-Z-F”.

- The “Sync PZF” key is disappeared and the “synchronization mode” is not active.
- Actuate the “Exit” key and hold for a few seconds to return to the main menu.

It is recommended to always work in “synchronization mode” to avoid an exsiccation of the tissue and a carburization.

8.6 Automatic mode

In addition to the “synchronization mode”, the HEBU HFiT®-T can be operated in the “Automatic mode”. In this case, only the power is changed, while the flow rate corresponding to the power is adjusted automatically. Conversely, when these parameters are changed in the “automatic mode” the power is automatically adjusted to the flow rate.


This is visualized by the symbol (21).

The relations between these two parameters were created from experimental studies and they are listed in the following table:


Power [W]	Flow rate [ml / h]
5	40
10	47
15	54
20	60
25	67
30	74
35	80
40	90
45	100
50	107
55	114
60	120
65	126
70	133
75	139
80	146
85	152
90	159
95	166
100	172

To deactivate the “automatic mode” proceed as follows:
Actuate the button “Power setting” (14) in the main menu.
It appears the menu “power setting”.

- Actuate the key “Auto P-F”.

- The “Auto P-F” key is disappeared and the “automatic mode” is not active.
- Actuate the “Exit” key and hold for a few seconds to return to the main menu.

8.7 Real-time monitoring

During the entire treatment period the tissue impedance is visualized via a bar chart (28). The tissue impedance in Ohm, the HF voltage in volts, the HF current in milliampere, the HF power in Watt and, so far as one works with a thermal sensor, the temperature in °Celsius are shown at the bottom of the LCD-display.

The total introduced HF energy is also shown on the LCD display (26). The value is reset after switching on or directly via pressing the display of the “introduced energy” (26).

8.8 Setting recommendations

The following parameters serve as guidelines for treatment with the HEBU HFIT®-T system:

Tumor diameter	Active component - needle	Art.-No. needle	Power	Flow rate	Time
10-20 mm	15mm	HF 9946-15	30-40 W	90-120 ml/min	5-8 min
15-25 mm	20mm	HF 9946-20	40-50 W	110-140 ml/min	7-10 min
20-30 mm	25mm	HF 9946-25	50-60 W	130-160 ml/min	9-12 min
25-35 mm	30mm	HF 9946-30	60-70 W	150-180 ml/min	11-14 min

Please note that the settings have to be individually adapted to the circumstances and entirely depending on the doctor’s advice. The HEBU medical GmbH assumes no responsibility for these settings.

9 Accessories

HF 9945-15	HEBU HFiT®-T Needle Electrode, \varnothing x shaft length x electrode length: 2,0 x 150 x 15 mm
HF 9945-20	HEBU HFiT®-T Needle Electrode, \varnothing x shaft length x electrode length: 2,0 x 150 x 15 mm
HF 9946 15	HEBU HFiT®-T Needle Electrode, \varnothing x shaft length x electrode length: 2,0 x 150 x 15 mm
HF 9946-20	HEBU HFiT®-T Needle Electrode, \varnothing x shaft length x electrode length: 2,0 x 150 x 15 mm
HF 9946-25	HEBU HFiT®-T Needle Electrode, \varnothing x shaft length x electrode length: 2,0 x 150 x 15 mm
HF 9946-30	HEBU HFiT®-T Needle Electrode, \varnothing x shaft length x electrode length: 2,0 x 150 x 15 mm

HF 9518-01	Perfusion Pump for 50ml syringe, 0,1 – 200 ml/h
HF 9418-02	RS232 Interface cable
HF 9529	Trolley
HF 9537-02	Double foot switch for HEBU HFiT®-T-system, class AP, 6 PIN plug, cable length 4,5 m
HF 9564 H	Single-use security electrode neutral with double seperated contact surface, Hydrogel, Conductive area 110 cm ² , 50 pcs per carton
HF 9949-00	Temperature probe, Type K, \varnothing 0,5 mm, with cable, length 3 m, sterilisable
HF 9949-01	Y-connector for insertion of the Temperature probe into the HEBU HFiT®-T-Needle and for connection of the perfusion hose, 5 pcs per carton
HF 9949-02	Connection hose with Luer-Lock-Connection, length 1,9 m, 25 pcs per carton
HF 9949-03	Syringe 50 ml, 50 pcs per carton

HF 9526-05	Potential equalization cable, 0,5 m
HF 9526-10	Potential equalization cable, 1,0 m
HF 9526-30	Potential equalization cable, 3,0 m
HF 9569-30	Cable for disposable return plates, 3 m
HF 9569-50	Cable for disposable return plates, 5 m
HF 9571-30	Monopolar cable for HEBU HFiT®-T Needles, sterilisable, 3 m
HF 9572-30	Bipolar cable for HEBU HFiT®-T Needles, sterilisable, 3 m
HF 9979-01	Mains cable, e.g. for Germany, France, Austria, Skandinavia, 2,5 m
HF 9979-02	Mains cable for Switzerland, 2,5 m
HF 9979-03	Mains cable e.g. for Taiwan, Korea, USA, Canada, 2,5 m
HF 9979-04	Mains cable for Italy 2,5 m
HF 9979-05	Mains cable for Great Britain, 2,5 m
HF 9979-06	Mains cable for Argentina and Australia, 2,5 m

10 Operating errors


Error	Cause of error	Remedial actions
"Timeout to reach power"	<ul style="list-style-type: none"> The voltage of the power supply is reached within < 700 ms, which is necessary for the set power on the display. 	<ul style="list-style-type: none"> Actuate the "Power button"(22) or the blue pedal of the foot switch.
"ERROR output POWER"	<ul style="list-style-type: none"> The output measured by the monitoring module is 20% greater than the set output. 	<ul style="list-style-type: none"> Actuate the "Power button"(22) or the blue pedal of the foot switch.
"NE Error "	<ul style="list-style-type: none"> The HF power has been activated before the neutral electrode was connected or detected. The neutral electrode is defective. 	<ul style="list-style-type: none"> Actuate the "Power button"(22) or the blue pedal of the foot switch. Check the correct connection of the neutral electrode. Check the neutral electrode for integrity. Replace as necessary.
"Patient Detect Error"	<ul style="list-style-type: none"> The HF power has been activated before the tissue of the patient was recognized. During the treatment, the contact between needle and patient has been interrupted. The HF cable is damaged. 	<ul style="list-style-type: none"> Actuate the "Power button"(22) or the blue pedal of the foot switch. Ensure that the needle electrode is positioned correctly. Check the HF cable of the needle electrode for integrity. Replace as necessary.
"Thermocouple Error"	<ul style="list-style-type: none"> At a set temperature, the HF power was activated without connecting a temperature sensor. The temperature sensor is defective. 	<ul style="list-style-type: none"> Actuate the "Power button"(22) or the blue pedal of the foot switch. Check the correct connection of the temperature sensor. Check the cable of the temperature sensor for integrity. For treatment without temperature sensor, switch off the temperature control.
"Power Supply Error - No 1"	<ul style="list-style-type: none"> Missing voltage at the power factor correction block. 	<ul style="list-style-type: none"> Switch off the HF device and switch on again.
"Power Supply Error - No 2"	<ul style="list-style-type: none"> The necessary voltage from the block for regulating the power voltage for the output generator in idle running is missing. 	<ul style="list-style-type: none"> Switch off the HF device and switch on again.
"Foot power control Error"	<ul style="list-style-type: none"> The necessary voltage from the block is missing for correction of the output factor. 	<ul style="list-style-type: none"> Actuate the "Power button"(22) or the blue pedal of the foot switch.

Error	Cause of error	Remedial actions
		<ul style="list-style-type: none"> ▪ Before switching off, ensure that the foot switch isn't activated. ▪ Check the foot switch for integrity and function.
"Foot fluid control Error"	<ul style="list-style-type: none"> ▪ When switching on the device, the green pedal is activated or short circuited. 	<ul style="list-style-type: none"> ▪ Actuate the "Power button"(22) or the blue pedal of the foot switch. ▪ Before switching off, ensure that the foot switch isn't activated. ▪ Check the foot switch for integrity and function.
"PILOTc communication Error"	<ul style="list-style-type: none"> ▪ Activating of the HF power without connection to the perfusion pump. ▪ Connection to the perfusion pump is interrupted. ▪ Perfusion pump isn't switched on. 	<ul style="list-style-type: none"> ▪ Switch off the HF device. ▪ Check the connection between pump and HF device. ▪ Switch on the HF device and the pump. Connect the pump as described in chapter 5.3.
"PILOTc operation Error"	<ul style="list-style-type: none"> ▪ Perfusion pump is defective. 	<ul style="list-style-type: none"> ▪ Switch off the HF device. ▪ Check the pump for function with help of the operating manual.
"ERROR in CPU ROM"	<ul style="list-style-type: none"> ▪ Incorrect checksum of the ROM memory from control processor. 	<ul style="list-style-type: none"> ▪ Switch off the HF device and switch on again.
"ERROR in CPU RAM"	<ul style="list-style-type: none"> ▪ Incorrect checksum of the RAM memory from control processor. 	<ul style="list-style-type: none"> ▪ Switch off the HF device and switch on again.

11 Performance diagram

The following diagram shows the half and the full power in dependence of the resistance.

12 Protective measures and warnings

Please observe also the instructions provided on working with the neutral electrode in chapter 6.1, the operating instructions with the accessories and the applicable safety measures!

	<p>Before any procedure, check the cables, accessories (handles, electrodes, forceps) and foot switches for any signs of visible irregularity. Do not use cables or instruments with broken or defective insulation. These can result in burns if they come into contact with the skin!</p>
---	--

The unit should only be operated at a mains socket with protective earth contact. When performing electrosurgical procedures, the risk of burns must be kept as small as possible. For this, the following safety measures must be observed:

	All plug-in connections may only be carried out with compatible accessories.
	Use the greatest care when working with the neutral conductor!
	Protect the patient from contact with earthed metallic elements (OP table, infusion stand etc.) during use.
	Avoid contact between supply cables and patients, other cables or instruments! Ensure that there are no loops in the cables of handles and electrodes or fixture to the OP table.
	Avoid small-area skin-to-skin contact in the patient (e.g. a finger making contact with the thigh).
	Avoid contact between HF instruments and non-insulated instruments.
	Prevent the penetration of fluids between the patient's body and the neutral electrode.
	The output power should not be set higher than is absolutely necessary for the procedure.
	Do not touch the patient during activation of the coagulation or cutting current. Should this be necessary, a piece of fabric can be used for insulation purposes.
	In patients connected to monitoring devices (ECG), the monitoring electrodes must be attached at the greatest possible distance to the application areas of the electrosurgical electrodes.
	The use of highly flammable narcotic substances and oxidizing gases (N ₂ O or O ₂) should be avoided. Should this be necessary, take particular precautions.

	Exclusively non-flammable disinfectants must be used. If alcohol-based solutions are required, wait before performing surgery until these have completely evaporated. Bear in mind that disinfectants can flow down and collect in recesses of the body. A spark from the active electrode can cause combustion of the fluid in normal use.
	Sparks from the active electrode can set bandages (textiles), metabolic gases and combustible fluids alight.
	In the case of procedures involving high-frequency current applied to body parts with a minimal cross-section, the bipolar technique should be used in order to prevent the risk of unintentional coagulation.
	In the case of procedures on body parts with a close link to the remainder of the body, the bipolar technique should be used for safety reasons.
	Surgical gloves do not afford the surgeon sufficient protection from burns!
 	During surgery, pacemakers are at risk of destruction. We recommend consulting the treating cardiologist and/or an authorized representative of the pacemaker manufacturer. As a rule, in such cases bipolar application technology is advisable. Electrosurgical procedures are prohibited on outpatients with pacemakers.
	A marked drop in power at a normal setting can be an indication of incorrect application of the neutral electrode.
	High-frequency creep currents can cause burns in areas distant from the application site through conductive connections.
	To remove the cable from the unit, always hold the plug! Never pull on the cable, as this can result in damage to the cable insulators, burns on the skin or also fire.
	Instruments which are temporarily not required during the procedure must be placed away from the patient on the instrument table or a similar surface.
	Reusable electrosurgical instruments which are connected to the unit are generally not delivered in a sterile condition!

Electromagnetic safety

In the development of electrosurgical units, HEBUmedical attaches particular importance to the stringent electromagnetic emission regulations. Consequently, the solutions we supply have a correspondingly minimal emission level for the required output. Measurements conform the high level of electromagnetic safety offered by HEBUmedical electrosurgical units.

Under typical working conditions based on an 8-hour day, a field occurs at a distance of 5 to 15 cm from the cables. At a distance of 20-40 cm, the value drops well below the prescribed limit.

Electromagnetic fields occur primarily around the cable.
The electrosurgical unit itself does not constitute a significant emitting element.

13 Review of the technical status

Before any procedure, check the cables, accessories (handles, electrodes, forceps) and foot switches for any signs of visible irregularity. Do not use cables or instruments with broken or defective insulation. These can result in burns if they come into contact with the skin!

13.1 Mechanical faults

If damage occurs to switches, connecting sockets, houses or the touch screen, or if the unit is dropped on the floor, please contact your authorized dealer for advice on further use.
The dealer can arrange for a detailed technical check to be performed.

The manufacturer accepts no liability for calibration work or repairs to the unit carried out by the user.

13.2 Recurring safety inspections

A safety inspection of the electrosurgical system is required once a year. The manufacturer only permits the use of units for which an up-to-date inspection by suitably trained personnel can be verified.

All information relating to the inspections (date, inspecting staff member, result etc.) must be recorded in the unit's log book.

13.3 Servicing and manufacturer address

Should you require the instructions for use in paper form, please use the contact details below. The instructions for use in paper form will be made available to you within seven calendar days of receipt of the request.

Alternatively, you can print out the electronic instructions for use yourself.

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen / Germany
Tel. +49 7461 94 71 - 0
Fax +49 7461 94 71 - 22
eMail: service@HEBUmedical.de
Web: www.HEBUmedical.de

14 System maintenance work

14.1 Cleaning

The electrosurgical unit can be wiped using a cloth dampened with standard disinfectant, and so kept clean using the simplest possible method. Avoid allowing cleaning fluid and disinfectant to penetrate into the unit.

14.2 Fuses

The fuse sockets are located on the back of the unit. The unit works with two fuses of 6.3 A each.

Exchanging the fuses:

- Disconnect the unit from the mains
- Release the fuse from the retainer using a flat blade screwdriver.
- Exchange the fuse

After an exchange of fuses, the unit must be checked by suitably qualified personnel.

14.3 Transport

Prevent mechanical damage and moisture.

If the unit has been stored for a long period in a cold environment, wait before removing the protective packaging until the unit has reached room temperature.

During transportation, the standard safety measures additionally apply.

Only ever transport the unit in its original packaging.

14.4 Environmental directives

From the implementation date of the European directive 2002/96/EU as national legislation, the following regulations apply:

- Electrical and electronic units must not be disposed of with normal household waste.
- The consumer is required by law to return electrical and electronic devices at the end of their service life to designated public collection points or to the place of purchase.

Disposal details are set out by the applicable local legislation in the relevant country. The symbol on the product, the instructions for use or the packaging indicates the relevant applicable legislation. By sending end-of-life devices for reuse, material recycling or other forms of utilization you are making a major contribution to environmental protection.

Please also send the unit packaging for environmentally responsible recycling after the end of the product's life.

15 Warranty card

Warranty card

We grant on all electrical medical products

We grant a warranty in accordance with the statutory and country-specific regulations on all electrical medical products (verification required in the form of an invoice or delivery note). The warranty period is two years. In the event of material or production errors, we provide free repair of the unit. Any attempt to carry out independent repairs will render the warranty cover null and void.

Damage caused by incorrect handling, overloading or normal wear and tear is not covered under the terms of the warranty.

In the event of a warranty claim and resulting repair, please enclose a copy of the invoice or delivery note (or batch number) with the instrument.

Manufacturer and Aftersales Service Address:

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen
Tel. +49 7461 9471 -0
Fax +49 7461 9471 -22
service@HEBUmedical.

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**Our vast range of products covers over 10 000 various instruments.
They can be found in Internet or request our catalog.**

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