HF 9523

HBS Argon

**Argon Modul / Argon Module / Module à l'argon /
Módulo de argón / Modulo Argon**

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Seite 3

DEUTSCH

Operating Manual

Page 27

ENGLISH

Mode d'emploi

Page 49

FRANÇAIS

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Página 71

ESPAÑOL

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Pagina 95

ITALIANO


HF 9523

Inhalt

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HF 9523

Contents

1	Symbol descriptions	28
2	Introduction	29
3	Scope of delivery	30
4	Scope	30
4.1	Inteded Use	30
4.2	Contraindications	31
5	Structure and basic components of the system	31
5.1	Connection with the HF- unit	33
5.2	Main menu of HBS Argon	33
5.3	Overview of symbols for operating modes and settings	34
6	Specifications	35
7	Accessories	36
7.1	Accessories	36
7.2	Foot switch	37
8	Safety functions	37
9	Commissioning the unit	38
9.1	Before starting work	38
9.2	Installation of the unit	38
9.3	Connection of Argon	39
9.3.1	Gas supply of the gas transmission system	40
9.3.2	Connection of the Handle	40
9.3.3	Setting of the Argon gas flow	40
9.3.4	Saving and selecting the program	41
9.4	Operating errors	41
10	Guideline value for power settings	42
11	Protective measures and warnings	43
11.1	Notes on the Argon- supported Coagulate and Cutting	43
11.2	Special notes for minimal-invasive surgery	44
11.3	Use of foot switches	44
12	Review of the technical status	45
12.1	Mechanical faults	45
12.2	Recurring safety inspections	45
12.3	Servicing and manufacturer address	45
13	System maintenance work	46
13.1	Cleaning	46
13.2	Repair	46
13.3	Fuses	46
13.4	Transport	47
13.5	Environmental directives	47
	Warranty card	48

1 Symbol descriptions

Symbol	Definition
	CE-labelling
	Attention
	Manufacturer
	Lot-description
	Reference code
	Medical device / FDA Prescription device
	Medical device
	Non sterile
	Keep away from sunlight
	Dry storage required
	Temperature
	(Electronic) instruction for use
	No intervention on patients with cardiac pacemakers.
	WEEE-labelling
	Green dot

2 Introduction

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

Nowadays the high-frequency surgery is the most commonly used method among the surgical techniques. In combination with argon gas, this technology can be used even more effectively.

The HBS Argon is an argon aggregate, which is designed with the high-frequency device to cut and coagulate biological tissue with high-frequency current under argon protection.

When the applicator is activated, the argon gas flows through the electrode. The high-frequency current flows through the plasma beam to the tissue.

A plasma consists of free electrons, positive ions and neutral atoms (or molecules). The conductivity of the plasma is determined by the amount of electrons and ions and their ability for moving free in the surrounding gas of ions and neutral atoms.

A plasma is generated out of the strong current field in immediate proximity of the electrode during the gas ionization.

The application of argon gas during cutting ensures cooling of the surrounding tissue and prevents access of oxygen. This leads to the prevention of burns, smell and tissue burning, which are subject to interference, by ensuring safe and reliable hemostasis.

The argon – plasma – coagulation is applied in the open laparoscopic thorax surgery, especially for nursing surface bleeding for 10 years. The thermic effect of the ionized plasma is used.

The advantages of argon – plasma – coagulation to the previous coagulation methods are the following:

- noncontact coagulation
- less blood loss, especially in the case of parenchyma
- possibility to avoid already coagulated areas
- limited penetration – up to 3 mm.
- effect as protective gas against carbonation

The most important performances of the HBS Argon are:

- LCD-Touchscreen-Display with intuitive, high-resolution user menu.
- Independent regulation of the output of the argon gas for both operating modes – cut und coagulating with their operating modes.
- Automatic switching of the gas supply from empty to full cylinder.
- controlled Argon supply with the pressure and output, set by the standards equal to the realistic clinical conditions.
- Acoustic alarm during the Activation of diverse operating modes as well as an alarm in case of errors.
- Storage of 8 programs, 4 of them with preset parameters (argon gas flow for cutting and coagulating) and 4 user programs.

Advantages of the HBS Argon:

- Independent regulation of argon gas flow rate for both operating modes – cutting and coagulating with their operating modes.
- Automatic switching of the gas supply from empty to full cylinder.
- Controlled argon supply with the pressure and output, set by the standards equal to the realistic clinical conditions.

- Acoustic alarm when various operating modes are activated, as well as alarm sound when errors occur.

The Operator of HBS Argon can use the maximum of all advantages of the device, if he considers this user’s guide. It will help you to work professional with the device and to choose the correct parameters.

3 Scope of delivery

HF 9523	HBS Argon
HF9979-01	Power cable
HF9412-03	Connecting cable
	Instruction for use

4 Scope

The device is developed for the application in the following disciplines:

- General surgery
- Traumatology
- Cardiac and Thoracic Surgery
- Neurosurgery
- Vascular surgery
- Graft surgery
- Abdominal surgery
- Urology
- ENT, jaw- and facial surgery
- Pediatric surgery
- Laparoscopic surgery
- Gynecology
- Orthopedy
- Gastroenterology

4.1 Inteded Use

An electrosurgical unit is a generator that converts electronic energy into high frequency alternating current (RF) power. If these currents flow through biological tissue, heat is generated, which is used to cut and coagulate the tissue.

The device has monopolar and bipolar application modes.

In monopolar mode, RF current is applied to the tissue with an active electrode. The cutting or coagulation effect is caused by a high concentration of current, i. a high current density over the small area of the active electrode. As a result, an increase in temperature is generated and the water evaporates from the tissue immediately surrounding the active electrode. Depending on the intensity of the HF current, this leads to a standstill of bleeding or to a cut in the tissue. The HF current flows from the active electrode to the neutral electrode and is distributed on its large surface area. Thus, there reduces the density of the current and there are usually no unintended thermal effects at the

position of the neutral electrode. Through the neutral electrode, the current flows back to the electro-surgery device.

In bipolar mode, the current path is limited only to tissue located between the two electrodes of a bipolar instrument. The use of the neutral electrode is not required.

Optionally, the HF units can be modularly expanded with an argon unit. The argon device is an argon unit designed for cutting and coagulation of biological tissue with high-frequency currents under argon protection.

4.2 Contraindications

Because of the increased risk of an arterial embolism, the application of argon gas is not recommended in patients with aortic-pulmonary septal defect or such with arterial septal defect.

5 Structure and basic components of the system


The housing of the HBS Argon is made of metal and glass. As no openings are available for ventilation, the unit is simple to clean (cf. 13.1).

Front view of the HBS Argon

Located on the front of the unit are:

- (1) Socket for connection of the output of the High-Frequency device
- (2) Socket for monopolar Instrument to cut /coagulate (working with two- and three pole handle of the Type Martin and ValleyLab);
- (3) Argon gas connection
- (4) LCD Touchscreen Display
- (5) Main switch (on/off)

Rear view of the HBS Argon


Located on the back of the housing are:

- (6) Socket for power cable
- (7) System outlet for connecting an additional device
- (8) Connection of the earth lead
- (9) Socket for Foot switch
- (10) Socket for switching on the sensor line of cylinder 2
- (11) Socket for the connection of the gas supply hose of cylinder 2
- (12) Socket for switching on the sensor line of cylinder 1
- (13) Socket for the connection of the gas supply hose of cylinder 1

5.1 Connection with the HF- unit

5.2 Main menu of HBS Argon

Main menu of HBS Argon

- (14) Display of the gas flow rate in mode “CUT”. Button – Entrance to the side for setting the output
- (15) Display of the gas output in mode “COAG”. Button to open the side for setting the output
- (16) Program selection button
- (17) Display of the amount of gas in cylinder 1
- (18) Display of the amount of gas in cylinder 2
- (19) Button – approve the gas operation mode
- (20) Button – Activation of the “cleaning” function of the orifice and emptying the systems after closing the gas cylinders.

5.3 Overview of symbols for operating modes and settings

Symbol	Meaning
Control elements	
	Increase output
	Reduce output
	On/ off switch of the Argon mode
	Activates flushing of the argon probe
	„Exit“: saves settings and returns to the main menu
Special denominations	
	See user's guide
	The neutral electrode is isolated from the earth potential at high frequencies
	Symbol, which indicates a "CF" -Class, being resistant to the impacts of a cardiodefibrillator
	Entrance for HF-Energy, connection to the HV-exit of the electrosurgical device
	Exit Argon Accessories
Symbols on the backplate	
	Socket for foot switch
	Grounding
	Electrical gas cylinder - connection

6 Specifications

Mains connection	
Supply voltage	100-260 Vac, 50-60 Hz
Consumed power – max.	18 VA
Fuses	2 St. T 1,0 A (inerte)
NF leakage current	
In accordance with	EN 60601-1; IEC 60601
Leakage currents against Earth	< 500 μ A
Leakage currents against Housing	< 1000 μ A
Patients Leakage currents	< 10 μ A
Patients – auxiliary stray currents	< 100 μ A
HF leakage currents – measured in a system with a electrosurgical device	
In accordance with	EN 60601-2; EN 60601
Neutral electrode against Earth	< 100 mA
Aktive electrode against Earth	< 100mA
Mode with Argon	
Argon Minimum purity	99.998 %
Argon gas sources	<ul style="list-style-type: none"> • 1 or 2 pieces gas cylinders with maximum pressure 200bar and recommended filling volume 25 l • Gas plant with highest inlet pressure 2,5bar
Control of the output on the exit	0,3 – 6,0 l/min with Step 0,1 l/min regardless for both operating modes Cut and Coagulation
Exchange empty cylinder with full cylinder	At output on entrance of the System <9 l/min
Safety measures	
Type of the construction according to	EN 60601; IEC 60601
Protection class	I
Type	CF
Earthing terminal	yes
Switching-on the neutral electrode	Floating Output

- Monitoring for wrong dosing of the output and switching-off the hf – generator.
- Monitoring of the fault conditions and switching-off the delivered output and hf – energy

Miscellaneous	
Storage	  
Surrounding temperature	+10°C to +60°C
Humidity	10 - 85 %
Air pressure during operation	500 - 1060hPa
During operation	
Surrounding temperature	+10°C to +40°C
Humidity	30 - 75 %
Air pressure during operation	700 - 1060hPa
Dimensions and weight	
Dimensions (L/W/H)	370 x 380 x 100 mm
Weight	5,9 kg

7 Accessories

The HBS Argon Module is compatible with the HBS 200i and HBS 300i.

Apply for a copy of our extensive electrosurgery accessories catalogue today!

To avoid an incompatibility and an unsafe operation, only use accessories recommended by the manufacturer when working with the HBS Argon.

Note:

Please observe the attached instruction for use, cleaning and preparation for the accessories.

Check the connection cables, handles and the coagulation instruments for interference before any sterilization. In particular, check the electrode cables and accessories used for endoscopic operations for defective insulation.

7.1 Accessories

HF9412-03: Connection cable among HF-device and HBS ARGON Module, 0,3m

7.2 Foot switch

The following foot switches can be used:

**HF 9533-02: Foot switch,
cable length 4,5m**

**HF 9535-01: Double foot switch,
cable length 4,5m**

8 Safety functions

Microprocessor tracking of the Safety functions

In order to minimize the risks of high- frequency electrosurgery, the system is equipped with a number of safety functions. The microprocessor control and the sensor technology of the system allow the implementation of a number of monitoring functions. The functional states, as well as the states of the interference in the generator are signalized visually and tonally, whereat in case of necessity the HF – current is switched off immediately.

The HBS Argon has an automatic, self-diagnostic control.

Monitoring of the incorrect dose

By following the output energy and power, the HBS Argon module always makes a comparison of the acceptable voltage limits and of the values selected on the control panel. If unacceptable deviations are reached, the emitted high-frequency energy is interrupted and an error message is issued.

Protection against wrong settings of the system

The HBS Argon offers:

- Facilitating control because of the intuitive design of the menus on the LCD – display and highly sensitive virtual buttons; clear, easy-to-understand symbols for selecting the operating mode and setting the output power.

9 Commissioning the unit

The device can only be installed and used in rooms for medical purpose according to its specification.

Only Persons who are trained by the manufacturer or by the suppliers and who use the device correctly according to the instructions, are allowed to work with the device.

Do not put vessels with liquids on or next to the device.

9.1 Before starting work

Make sure that the surface on which the device is placed is at least as large as the device and is not inclined by more than 10° from the horizon. HEBUmedical recommends placing the device on a trolley which is provided by the company.

Check, that the cooling openings are not blocked on the lower device side and the air can circulate freely.

9.2 Installation of the unit

Mains connection

- Connect mains cable
- Connect the device to the mains only with the delivered mains cable or with a cable of the same quality, which has the required sign for passed test from the national legislation.
- Connect the device with the yellow – green earthing cable on the earth socket (8), which is located besides the mains supply socket.
- Ensure that there are no other connections and switch it on with the On /Off switch in the lower right corner of the front of the device to check if it lights up. Then switch the device off again.

Foot switch

The foot switches used are of 2 types:

- Double- foot switch: The yellow section is used to activate cutting and the blue section control coagulation.
- Single-foot switch: blue section – for coagulation.

The connection of both foot switches occurs by the socket (9) for the foot switch on the backplate of the device.

Auto-test when switching on

Caution: Make sure that all connections are connected correctly and that the active electrodes do not have contact with the patient, a third person or with a conductive material, before each switch-on.

Switch-on the device. You can hear an original message on the LCD Panel the welcome page is displayed and the auto-test is carried out.

Error in the system

Do not use the device in any case, if a error message appears at a given time and the LCD monitor glows in a red colour!

9.3 Connection of Argon


Use Argon gas with a purity of at least 99,998 % or higher. At a low purity, it is recommended to use a membrane- gas filter in order not to endanger the safety of the patient.

The HBS Argon Module can work with gas from a cylinder as well as with gas from a gas transmission system

Use of gas cylinders


- Secure the gas cylinder (A) with the provided fixture (B)
- Do NOT remove the Identification label
- Mount an original HEBU pressure relief valve HEBU (C) at the gas cylinder (article number: HF 9530-02)

Caution!

The gas cylinders can only be connected to the HBS Argon with original valves and connecting hoses.

- Connect the hose (A) of gas cylinder 1 to connection point 1 (21) and the gas pressure measurement line (B) from the valve to terminal 22.
- Connect the hose (C) of gas cylinder 2 to connection point 2 (24) and the gas pressure measurement line 2 (D) to the corresponding terminal 23.

When the pressure in the active gas cylinder drops to below 5% the HBS Argon automatically switches to the full gas cylinder. If both cylinders are empty, the message "Argon level LOW" appears on activation at the display.

9.3.1 Gas supply of the gas transmission system

Connect the hose of the gas transmission system with one of both gas sockets.


9.3.2 Connection of the Handle

Connect the monopolar output (25) of the HF-device to the input of the Argon module (26).

Connect the coupling of the handle to the HF output (27) and the hose connector (of the Handle) to the gas supply socket (28).

9.3.3 Setting of the Argon gas flow

The gas flow can be set in an interval of 0,3 L/min to 6,0 L/min in steps of 0,1 l/min.


The gas quantity in the cylinder is shown on the display by filling in the image.

Different gas flows can be set for each of the cutting CUT and coagulate COAG functions.


To set the gas flow, touch the correspondent display- CUT FLOW or COAG FLOW. Then the menu for setting the gas flow opens. Use the arrow keys to set the desired gas flow.

Select the "Exit" field to save and return to the main menu.

Press the flow display under „CUT FLOW" – The menu for selecting quantity in "CUT" mode opens:

Press the flow display under „COAG FLOW" –the menu for selecting gas the gas quantity in the "COAG" mode opens:

By pressing the "Argon Flow" -Button, the gas operation mode is permitted. The display changes the colour of the push button from red to green.

After finishing the process, it is advising to stop the gas supply to the module and to clean the system and the orifice with the Button "C".

9.3.4 Saving and selecting the program

By touching the program name on the main side opens a program selection menu.

Four factory programs (e.g. Gastro F3,2) are set, whose parameters can be changed. The factory settings are retained when the device is switched off and on again.

For example, four client-programs are entered, whose parameters can be set. After shut down, the last set values are saved.

Program: Gastro Φ 2.3
Program: Gastro Φ 3.2
Program: Minimal invasive surgery
Program: Open surgery
Program: User 1
Program: User 2
Program: User 3
Program: User 4

9.4 Operating errors

Following Error messages can occur:

Error	Localization	Explanation	Cause
11: "EEROR in CPU ROM"	CPU_Failure	Wrong check sum of the ROM – Memory of the control processor	Controller
12: "EEROR in CPU RAM"		Wrong check sum of the RAM – Memory of the control processor	
13: „ARGON level LOW"	ARGON_Block	The pressure in the cylinders is below the required minimum or there is an error in the pressure measurement system	Argon module
14: „ARGON no FLOW"		The flow volume set at the display cannot be reached	

10 Guideline value for power settings

When performing the power settings, the user should generally take into account the cutting speed, geometry of the electrode and properties of the tissue.

Monopolar cutting:

	Recommended settings
Smooth cuts with thin needle electrode	5 - 40 W
Cutting with knife electrode, lancet or wire loop	50 - 175 W
Cutting in wet environment (TUR)	80 - 150 W

Monopolar coagulation:

	Recommended settings
Coagulation of large surfaces	60 W
Coagulation of small surfaces	40 W

Argon Coagulation:

	Recommended settings
With flexible probe	20-40 W 1,5-2,5 L/min
With inflexible electrode	25-50 W 2,5-4,5 L/min

Bipolar application:

	Recommended settings
Coagulation	30 – 60 W
Cutting with sloughing	40 – 60 W
Cutting in wet environment (TUR)	50 – 75 W
Vessel sealing, BiSeal®	50 – 120 W

The recommended settings serve as guideline values only and must be adjusted to the relevant situation. In the event of unfavourable conditions, these values may deviate!

11 Protective measures and warnings

Before any procedure, check the cables, accessories (handles, electrodes, forceps) and foot switches for any signs of visible irregularity. Do not use cables or instruments with broken or defective insulation. These can result in burns if they come into contact with the skin!

The unit should only be operated at a mains socket with protective earth contact. When performing electrosurgical procedures, the risk of burns must be kept as small as possible. For this, the following safety measures must be observed:

The unit is not allowed to be used in areas with explosive risk.

The reparation of the unit is only allowed to be carried out by HEBU or by HEBU authorized agent.

It is necessary to ensure, that the unit is used only by persons knowing the user's guide.

The unit and the connecting cables should be checked before each use.

Interferences must be eliminated in time.

Make sure, that you are familiar with all features and safety of the electrosurgical device connected with the HBS Argon Module.

11.1 Notes on the Argon- supported Coagulate and Cutting

High purity gas (99.998% or higher) must be used. The use of argon gas with a lower purity can be dangerous for the patients.

The studies show that the argon gas that has reached the blood is excreted from the lungs during respiration. For young children and patients with respiratory diseases, the exhalation of the argon gas may be slower.

The HBS Argon should only be used with Argon gas. The use of other gases can lead to a risk of explosion.

It is required, that the gas cylinder and the connection component are tested and certified and are equal to the standards and the regional requirements. The gas cylinder should not be operated under pressure that is higher than specified.

	Only use self-adhesive neutral electrodes.
	Unsecured gas cylinders and / or defective connections on the gas transmission lines can lead to an accident.
	Check that the accessories have the required functionality. The defective connections can lead to sparking, damage of the accessories or to ignition.


11.2 Special notes for minimal-invasive surgery

	During minimal-invasive operations, the recommended insufflation pressure should not be exceeded. The additional use of argon gas increases the risk of increasing the intra-abdominal pressure. To minimize the risk, the gas supply should not exceed 5L/min.
	During minimally invasive surgery, an accumulation of gas in the abdominal cavity can lead to gas embolism.
	Even if the HF-current is switched off, the tip of the electrode can still be hot, which can lead to burns.
	An activation by accident or movement of the activated electrode outside of the field of view can lead to injuries of the patient.

11.3 Use of foot switches

For safe operation of the foot switches, observe the following:

- Ensure that only Class AP foot switches are used in the environment described as „medical“.

(1) Medical environment

Medical environment 1 includes the area where flammable mixtures with air can cause an explosion with a normal ignition energy. The anaesthetic gases and the vapours of the disinfectants are heavier than the air and fall to the floor.

There is no risk of explosion outside this area because of the stronger dilution of these substances.

12 Review of the technical status

Before any procedure, check the cables, accessories (handles, electrodes, forceps) and foot switches for any signs of visible irregularity. Do not use cables or instruments with broken or defective insulation. These can result in burns if they come into contact with the skin!

12.1 Mechanical faults

If damage occurs to switches, connecting sockets, houses or the touchscreen, or if the unit is dropped on the floor, please contact your authorized dealer for advice on further use.

The dealer can arrange for a detailed technical check to be performed

The manufacturer accepts no liability for calibration work or repairs to the unit carried out by the user.

12.2 Recurring safety inspections

A safety inspection of the electrosurgical system is required once a year. The manufacturer only permits the use of units for which an up-to-date inspection by suitably trained personnel can be verified.

All information relating to the inspections (date, inspecting staff member, result etc.) must be recorded in the unit's log book.

12.3 Servicing and manufacturer address

Please contact your supplier in case of a problem that is insoluble for the operator or if the service is necessary. When installing the device, it informs you about the possibilities for technical maintenance and service at the specific point.

If you have trouble with handling the device at your place, please contact the manufacturer of HBS Argon

Should you require the instructions for use in paper form, please use the contact details below. The instructions for use in paper form will be made available to you within seven calendar days of receipt of the request.

Alternatively, you can print out the electronic instructions for use yourself.

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen / Germany

Tel. +49 7461 94 71 - 0

Fax +49 7461 94 71 - 22

eMail: service@HEBUmedical.de

Web: www.HEBUmedical.de

13 System maintenance work

13.1 Cleaning

	<p>When cleaning the device:</p> <p>Always disconnect the appliance from the mains before cleaning.</p> <p>Before connecting the device to the mains, ensure that all flammable cleaning agents have completely evaporated.</p>
---	--

Observe the local hygiene instructions:

- Clean all surfaces of the device (the front side included) with cleaning and disinfecting agents without alcohol content.
- Follow the instructions of the manufacturer for the disinfecting agents.
- Wipe the agents off with a sponge or cloth moistened in clean water.
- Dry the device with a clean, lint-free cloth.
- Do not clean or disinfect the surfaces with flammable or explosive products. If their use cannot be avoided, flammable or explosive products must be vaporised completely before switching on the device.
- Propanol and ethanol affect the surfaces.
- Do not immerse the device in liquids.
- Do not sterilize the device with steam, ethyl oxide or other agents.

13.2 Repair

	<p>In case of service:</p> <p>The repair of the device should only be carried out by HEBUmedical or by HEBUmedical trained and authorized staff.</p>
---	--

Do not open the device housing. The HBS Argon does not contain any user-replaceable parts. The system may only be repaired by a qualified service engineer, who is specially trained to repair the HBS Argon.

Interventions of other persons can lead to cancellation of the guarantee and pose a risk for health.

13.3 Fuses

The fuses are located on the back of the device in socket (6), where the mains cable is connected.

The fuses must be replaced as follows:

- Use a small screwdriver to remove the fuse holder.
- Replace defective fuses with the type described in chapter 5.
- Put the fuse holder back in place.

After the exchange of the fuses, the device must be checked by qualified personnel.

13.4 Transport

Prevent mechanical damage and moisture.

If the unit has been stored for a long period in a cold environment, wait before removing the protective packaging until the unit has reached room temperature.

During transportation, the standard safety measures additionally apply.

Only ever transport the unit in its original packaging.

13.5 Environmental directives

From the implementation date of the European directive 2002/96/EU as national legislation, the following regulations apply:

- Electrical and electronic units must not be disposed of with normal household waste.
- The consumer is required by law to return electrical and electronic devices at the end of their service life to designated public collection points or to the place of purchase.

Disposal details are set out by the applicable local legislation in the relevant country. The symbol on the product, the instructions for use or the packaging indicates the relevant applicable legislation. By sending end-of-life devices for reuse, material recycling or other forms of utilization, you are making a major contribution to environmental protection.

Please also send the unit packaging for environmentally responsible recycling after the end of the product's life.

Warranty card

On all electrical medical products we grant

We grant a warranty in accordance with the statutory and country-specific regulations on all electrical medical products (verification required in the form of an invoice or delivery note). The warranty period is two years. In the event of material or production errors, we provide free repair of the unit. Any attempt to carry out independent repairs will render the warranty cover null and void.

Damage caused by incorrect handling, overloading or normal wear and tear is not covered under the terms of the warranty.

In the event of a warranty claim and resulting repair, please enclose a copy of the invoice or delivery note (or batch number) with the instrument

Manufacturer and Aftersales Service Address:

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen
Tel. +49 7461 9471 -0
Fax +49 7461 9471 -22
service@HEBUmedical.

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