HF 9504

HRF 100 x-touch

Elektrochirurgie-Gerät / Electrosurgical Unit / Aparato de cirugía electrodoméstica / Apparecchie di elettrochirurgia

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Seite 3

DEUTSCH

Operating Manual

Page 35

ENGLISH

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Página 67

ESPAÑOL

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Página 97

ITALIANO

HF 9504

Inhalt

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
HF 9505

Contents

1	Symbol descriptions	36
2	Introduction	37
2.1	Monopolar operation	38
2.2	Bipolar operation	38
3	Scope of delivery	38
4	Scope	39
4.1	Inteded Use	39
4.2	Contraindications.....	39
5	Structure and basic components of the system	40
5.1	Main menu of HRF 100 x-touch	42
5.2	Overview of symbols for operating modes and settings	43
6	Technical data	45
7	Accessories	46
7.1	Accessories for the monopolar operating mode	46
7.2	Accessories for the bipolar operating mode.....	47
7.3	Foot switch	47
8	Commissioning the unit	48
8.1	Before starting work	48
8.2	Installation of the device.....	48
8.3	Switch on.....	49
8.4	Foot switch settings.....	49
8.5	Handles with buttons	50
8.6	Neutral electrode settings and attachment.....	50
8.7	Setting the operating mode and power level.....	53
8.8	Saving and accessing the program.....	55
8.9	Saving the program against changes.....	56
8.10	Switching off	56
8.11	Security functions	56
8.12	Error handling	57
9	Guideline value for power settings	58
10	Protective measures and warnings	59
11	Review of the technical status	62
11.1	Mechanical faults.....	62
11.2	Recurring safety inspections	62
11.3	Servicing and manufacturer address	62
12	System maintenance work	63
12.1	Cleaning	63
12.2	Fuses.....	63
12.3	Transport.....	63
12.4	Environmental directives	64
	Warranty card	65

1 Symbol descriptions

Symbol	Definition
	CE-labelling
	Attention
	Manufacturer
	Lot-description
	Reference code
	Medical device / FDA Prescription device
	Medical device
	Non sterile
	Keep away from sunlight
	Dry storage required
	Temperature
	(Electronic) instruction for use
	No intervention on patients with cardiac pacemakers.
	WEEE-labelling
	Green dot

2 Introduction

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

An electrosurgical unit is a generator used to transform electronic energy into high-frequency alternating current (HF current).

If this current flows through human tissue, it creates heat which is used for cutting and for coagulation of the tissue. Electrosurgery is one of the most important energy methods in surgery.

Of all surgical disciplines, it remains the most effective combination tool between tissue cutting and hemostasis. The HRF 100 x-touch improves the surgical results thanks to 4.0 MHz radio wave technology.

The high frequency of 4.0 MHz minimizes temperature losses and thus does not damage the cells in the vicinity. The clinical advantages include the reduction of post-operative discomfort, the minimization of necrosis, the maximum recognition of the histological specimens, improvement of the healing process with excellent aesthetic results.

The HRF 100 x-touch is the best user-friendly electrosurgical unit for outpatient and clinical applications. The maximum output power is 100 W, both in monopolar and bipolar operation. The display, the selection of the modes and the output power values are carried out on a color screen.

The HRF 100 x-touch is a versatile and compact electrosurgical device that meets a high safety standard and combines economy with operational flexibility.

It offers the following application modes:

For monopolar surgical use

- Pure cut
- Blend cut
- Soft/Contact coagulation (no sparks/desiccate mode)
- Forced coagulation

For bipolar use:

- Bipolar coagulation (mode "Bipolar Forceps")
- Bipolar coagulation (mode "Bipolar Forceps AUTO Start")

The most important functions of HRF 100 x-touch are:

- Better control and safety, thanks to the integrated microprocessor technology in the modules for activating high frequency energy (HF energy) via the handles and foot switches.
- High level of security thanks to the microprocessor module for monitoring the neutral electrode and the interior design for operation with separate neutral electrodes
- Autotest after switching on
- One output for monopolar operation "Mono" and one output for bipolar operation "Bi".
- Polyphonic tone signaling with the possibility of regulating the tone strength when activating the various operating modes, as well as information tone without the possibility of switching off in the event of errors
- Light indication that uses the panel's RGB display for different operating modes or error states.
- Saving of 9 basic programs for monopolar and bipolar work (entered power and corresponding operating model)

- HEBU generators have an automatic output dosage by impedance with voltage regulation and arcing regulation and an automatic power adjustment with dynamic control and stop in case of anomalies.

2.1 Monopolar operation

In the monopolar mode, HF current is applied to the tissue using an active electrode. The cutting or coagulation effect is caused by a high concentration of current, i.e. a high current density over the small surface of the active electrode. This generates an increase in temperature, causing the water from the tissue immediately surrounding the active electrode to vaporize. Depending on the intensity of the HF current, this causes bleeding to stop or generates a cut in the tissue.

The HF current flows from the active electrode to the neutral electrode and is distributed over its large surface. This reduces the density of the current in this area and generally prevents unwanted thermal effects at the position of the neutral electrode. The current flows back to the electrosurgical unit through the neutral electrode

2.2 Bipolar operation

In the bipolar mode, the current path is limited only to tissue located between the two electrodes of a bipolar instrument.

Use of the neutral electrode is not necessary. Bipolar application does not entail any risk of the current flowing through the patient's body. Consequently, there is no danger of burns due to a poorly or incorrectly positioned neutral electrode.

Bipolar coagulation is therefore safer than the monopolar method and is particularly recommended for patients with pacemakers or for operations on organs with a small diameter.

3 Scope of delivery

HF 9504	HRF 100 x-touch
HF 9979-XX	Power cable
	Operating instructions in five languages
	Device manual in five languages
	Service manual in five languages

4 Scope

The use of electrosurgical equipment is predominantly in the following areas:

- General surgery
- Vascular Surgery
- ENT
- Gynecology
- Neurosurgery
- Pediatric surgery
- Laparoscopic surgery
- Plastic surgery
- skin surgery

4.1 Inteded Use

An electrosurgical unit is a generator that converts electronic energy into high frequency alternating current (RF) power. If these currents flow through biological tissue, heat is generated, which is used to cut and coagulate the tissue.

The device has monopolar and bipolar application modes.

In monopolar mode, RF current is applied to the tissue with an active electrode. The cutting or coagulation effect is caused by a high concentration of current, i. a high current density over the small area of the active electrode. As a result, an increase in temperature is generated and the water evaporates from the tissue immediately surrounding the active electrode. Depending on the intensity of the HF current, this leads to a standstill of bleeding or to a cut in the tissue. The HF current flows from the active electrode to the neutral electrode and is distributed on its large surface area. Thus, there reduces the density of the current and there are usually no unintended thermal effects at the position of the neutral electrode. Through the neutral electrode, the current flows back to the electrosurgery device.

In bipolar mode, the current path is limited only to tissue located between the two electrodes of a bipolar instrument. The use of the neutral electrode is not required.

Optionally, the HF units can be modularly expanded with an argon unit. The argon device is an argon unit designed for cutting and coagulation of biological tissue with high-frequency currents under argon protection.

4.2 Contraindications

HF surgical devices should not be used with simulators such as bone stimulators. It also discourages the use of electrosurgical devices when pacemakers, such as cardiac pacemakers, are in use.

5 Structure and basic components of the system

Front view of HRF 100 x-touch

Rear view HRF 100 x-touch

The housing of the HRF 100 x-touch is made of metal and glass. As no openings are available for ventilation, the unit can be cleaned easily (cf. 12)


Located on the front of the unit are:

- **(1)** Terminal for a neutral electrode
- Connections for
 - **(2)** one bipolar instrument: 2-pin standard, 5-pin plug or MARTIN plug
 - **(3)** one monopolar instrument: 3-pin standard or MARTIN plug
- **(4)** LCD touch screen display
- **(5)** Main switch (on/off)

Located on the back of the housing are:

- **(6)** Volume controller
- **(7)** Rating plate with manufacturer's information regarding serial number, mains voltage, output and fuses
- **(8)** Connections for foot switch for monopolar or bipolar mode (mono / bi)
- **(9)** Earthing pins
- **(10)** Terminal for a power supply cable (100-260 V, 50-60 Hz)

5.1 Main menu of HRF 100 x-touch

- **(11)** Taste - Programmwahl
- **(12)** Input for program selection, name of the program
- **(13)** Display / key for activating or deactivating the value saving - performance and modes when changing in the corresponding program
- **(14)** Display / key of the set output power (adjustable in ten levels) for monopolar cutting for output MONOPOLAR
- **(15)** Display / key for the type of monopolar cutting for the MONOPOLAR output
- **(16)** Display / key of the set output power (adjustable in ten levels) for monopolar coagulation for output MONOPOLAR
- **(17)** Key for the type of monopolar coagulation for the MONOPOLAR output
- **(18)** Key for the set maximum power (adjustable in ten levels) for bipolar coagulation for the BIPOLAR output
- **(19)** Bipolar coagulation button for BIPOLAR output
- **(20)** Display of the neutral electrode type
- **(21)** button - selection of the foot switch functions

5.2 Overview of symbols for operating modes and settings

Symbol	Meaning
	Increase output
	Reduce output
	Pure Cut
	Blend Cut
	Soft/Contact coagulation (no sparks/desiccate mode)
	Forced coagulation
	Bipolar coagulation (Mode "Bipolar Forceps")
	Bipolar coagulation (Mode "Bipolar Forceps" AUTO Start")
	„Exit“: saves the setting and return to the main menu

Symbol	Status	Meaning
NEUTRAL TYPE		Undivided neutral electrode detected
		Divided neutral electrode detected
	 alarm	Neutral electrode not detected or not connected
	MONO	The foot switch is set to MONOPOLAR
	BI	The foot switch is set to BIPOLAR

Special terms:

The device is resistant to the following influences

Symbol	Meaning
	See operating instructions
	Power switch: on / off
	The neutral electrode is isolated from the earth potential at high frequencies
	Symbol class "CF", resistant to the influence of cardio defibrillator

Designations on the back side

Symbol	Meaning
	Volume control
	Grounding
	Pedal socket

6 Technical data

Mains connection	
Supply voltage	100-260 Volt, 50 - 60 Hz.
Power consumption	160 VA (100 Volt, PF = 1,0) 150 VA (230 Volt, PF = 0,99)
Continuous operation	Interval 10s / 30 s (activation / pause)
Nominal frequency	4 MHz
Output power in monopolar operation	
Pure cut	100 W at 400 Ω
Blend cut	100 W at 400 Ω
Soft/Contact coagulation (no sparks/desiccate mode)	100 W at 75 Ω
Forced coagulation	100 W at 400 Ω
Output power in bipolar operation	
Bipolar coagulation (Mode „Bipolar Forceps“)	100 W at 50 Ω
Bipolare Koagulation (Modus „Bipolare Pinzette AUTO Start“)	100 W at 50 Ω
Safety measures	
Type	CF
Protection class	I
Equipotential connection	Yes
Safety	In accordance with DIN EN ISO 60601-1
Electromagnetic compatibility	In accordance with DIN EN ISO 60601-1-2
Safety of high-frequency surgical units	In accordance with DIN EN ISO 60601-2-2
Miscellaneous 	
Operating temperature	+10°C to +40°C
Humidity during operation	30-75%
Air pressure during operating	700-1060hPa
Storage temperature	-10°C to +60°C
Storage humidity	10-85%
Air pressure in storage	500-1060hPa
Dimensions (L/W/H)	370 x 370 x 145 mm
Weight	8,5 kg

7 Accessories

Request a copy of our extensive
electrosurgery accessories catalogue today!


In order to avoid incompatibility and unsafe operation, only use the accessories recommended by the manufacturer when working with the HRF 100 x-touch.

7.1 Accessories for the monopolar operating mode


Handles for monopolar electrodes can be connected to the unit using a 2-pin standard or MARTIN plug or using the 4 mm plug (banana plug) with only one contact.

3-pin standard plug

MARTIN monopolar
plug

Monopolar 4mm 1-pin plug
(banana plug)

Neutral electrode

Valleylab plug for the neutral electrode

The HRF 100 x-touch detects neutral electrodes with undivided and divided contact surface. A comparison of the electrode detected by the machine and the used electrode provides a direct indication of correct adhesion of the neutral electrode. Both single-use and reusable neutral electrodes can be used.

HEBUMedical recommends the use of the single-use safety
neutral electrode with double divided contact surface
(HF9564H) for maximum safety!

7.2 Accessories for the bipolar operating mode

Bipolar instruments are connected to the electrosurgical unit using a 2-pin standard, 5-pin or Martin plug.

A 5-pin plug is used for automatic detection of instruments

MARTIN bipolar plug

2-pin-plug for bipolar instruments

5-pin-plug for automatic detection of instruments

7.3 Foot switch

The following foot switches can be used:

HF 9533-02: Foot switch, cable length 4,5m

HF 9535-01: Double foot switch, cable length 4,5m

HF: 9537-01: Wireless double foot switch, Incl. transmitter

8 Commissioning the unit

8.1 Before starting work

According to its specification, the device can only be installed and used in rooms for medical purposes.

Only people who have been trained by the manufacturer or the supplier and who operate the device correctly in accordance with the instructions for use may work with the device.

Do not place any containers with liquids on or near the device.

Make sure that the surface on which the device is placed is at least as large as the device and not inclined by more than 10 ° to the horizon. HEBUmedical recommends putting the device on a trolley provided by the company.

Check that the cooling openings on the lower side of the device are not blocked and that air can circulate freely

8.2 Installation of the device

Mains connection

- Connect the power cord
- Only connect the device to the mains with the power cord supplied or with a cable of the same quality that has the mark required by national legislation for a passed test.
- Connect the device with the yellow - green earth cable to the earth socket (8), which is located next to the mains supply socket.
- Make sure that there are no other connections and switch on with the on / off switch in the lower right corner of the front of the device to check whether it lights up. Then switch off the device again

Autotest at switching on:

Caution: Before switching on, make sure that all connections are correctly connected and that the active electrodes are not in contact with the patient, with a third person or with a conductive material.


Turn on a device. You will hear an original message, the LCD panel will display the start page and an autotest will be carried out.

Errors in the system

Under no circumstances should you use the device if an error message appears at any given time!

Volume control

The volume of the activation and message signals can be adjusted using the controller on the rear of the device. The speakers are also on the back. The sounds cannot be turned off.

8.3 Switch on

When using an electrosurgical unit which produces HF current, always be aware of two fundamental principles:

Current flows through everything located between the active and the neutral electrode (monopolar)

The system is switched on by actuating the main switch (ON/OFF). However, high-frequency current only flows once a button is activated at the handle or foot switch.

Consequently, accessories can also be connected to the unit when it is switched on.

However, extreme care should always be taken when doing so!

Always ensure that the system is not inadvertently activated during installation by pressing the foot switch or the button in the handle.

8.4 Foot switch settings

After connecting the foot switch at the back, the required mode must be set at the touchscreen. To set the mode, touch the field indicating the foot switch symbol. This changes the text which appears above the symbol.

8.5 Handles with buttons

Activating the handle with buttons takes place independently of the foot switch settings. These cannot be activated using the foot switch.

Make the necessary connections with the device switched off!

8.6 Neutral electrode settings and attachment

In the monopolar operating mode, a neutral electrode is required.

Do not use undivided but only separate neutral electrodes.

- Place the neutral electrode on the patient.
- Connect the cable of the neutral electrode to the electrode.
- Connect the cable of the neutral electrode to the socket (3) on the front of the device.

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If the displayed neutral electrode type does not agree with the used neutral electrode type, this is indicative of poor adhesion of the neutral electrode on the patient, a defective neutral electrode or a defective cable!

Please note the following information when attaching the neutral electrode

Application information for single-use neutral electrodes

- When using single-use products, always pay attention to the use-by date. Do not use electrodes which are over this date!
- Neutral electrodes designed for single use must not be used a second time!
- The electrode must not be attached to surfaces with a high amount of hair growth. If necessary, shave the relevant area before attaching.
- The neutral electrode must be attached in such a way that its entire surface adheres to the patient's body.
- After every change of the patient's position, check that the neutral electrode is adhering correctly and also check the cable routing.
- As soon as the electrode has been attached, it must not be pulled off and attached again. Should the electrode have to be repositioned, use a new neutral electrode.
- The neutral electrode must not come into contact with fluids and must not be coiled round.
- Do not apply any additional electrically conductive gel to the neutral electrode.
- When performing surgery on small children, special electrodes suitable for children or babies must be used (see manufacturer's instructions)
- Never under any circumstances use damaged products!
- When removing the neutral electrode, never pull on the cable!
When pulling off single-use electrodes, ensure that the patient's skin is not damaged. Avoid any abrupt movements

Application information for reusable neutral electrodes


- The neutral electrode must be wiped with disinfectant before repeated use. To disinfect reusable rubber neutral electrodes, we recommend the products Incidin perfekt, Minutil and Incidur F from Ecolab.
- The neutral electrode must be attached in such a way that its entire surface adheres to the patient's body. A rubber band provided for this purpose can be used to ensure better adhesion.
- When attaching the neutral electrode, ensure that the long side is pointing towards the operation field.
- After each change in the patient's position, check the correct adhesion of the electrode and the connecting cable!
- The electrode must not be attached to surfaces with a high amount of hair growth. If necessary, shave the relevant area before attaching.
- The electrode must not be wet or coiled round.
- Fluids between the patient's skin and the neutral electrode must be avoided.
- Do not apply any additional electrically conductive gel to the neutral electrode.
- When performing surgery on small children, special electrodes suitable for children or babies must be used (see manufacturer's instructions)
- When removing the neutral electrode, never pull on the cable!
- Never under any circumstances carry out repairs on the neutral electrode!

Please note that rubber electrodes lose their conductive capability if active substances are eroded due to cleaning of the material. This type of electrode considerably increases the risk of burns. Consequently


ensure that not only the unit but also the reusable neutral electrode is included in regular maintenance inspections.

Areas for attachment of the neutral electrode


- Attach the neutral electrode to clean, dry skin.
- The electrode must be attached close to the operation area, but with a minimum direct distance from the operating field of 20 cm.
- Attach over tissue with a good supply of blood without dips or curves in the skin, for instance on the upper arm or thigh. (cf. picture below)
- Do not attach the neutral electrode to projecting surfaces.
- Do not attach above bones, scars, cuts or scratches.
- The electrode must not be attached to surfaces with a high amount of hair growth. If necessary, shave the relevant area before attaching.
- Do not attach above particularly adipose areas such as the abdomen or buttocks.
- The neutral electrode must not be attached above implants.
- Electrodes of monitoring systems must be placed at a distance of at least 20 cm from the operation area and from the position of the attached neutral electrode
- Correct attachment of the neutral electrode on adult patients


- Correct attachment of the neutral electrode on children

- Attachment on patients with pacemakers
Note! Consult the treating cardiologist!

RIGHT


WRONG!

8.7 Setting the operating mode and power level

■ Operating mode


HRF100 x-touch has two different operating modes for monopolar cutting, two for monopolar coagulation and two operating modes for bipolar applications (see chapter 4.2 Overview of the symbols for operating modes and settings). The selection of all modes is based on the same principle. The field with the illustration of the operating modes must be touched. The selection menu then opens.

The required operating modes must be selected by touching. To save and return to the main menu, select the "Exit" field.


- (22) button - cutting type I "cutting with minimal tissue injury
- (23) button - cutting type II "cutting with medium Tissue injury
- (25) button - Save the selected mode and return to the work page

Menu: Selection of operating modes (MONO - CUT)

- (25) button – Contact coagulation
- (26) button – forced coagulation
- (27) button – Save the selected mode

Menu: Power setting (MONO- CUT)


- (28) button – Bipolar forceps
- (29) button – Bipolar forceps with automatic start
- (30) button – Save the selected mode and return to the work page

Menu: Power setting (BIPOLAR – COAG)

■ Power setting

The power selection menu is identical for all modes.

To set the required power level, touch the power display. The menu for power level selection appears.

Use the arrow keys **(32/33)**, in order to set the required power setting **(31)**.

To save and return to the main menu, select the "Exit" field. **(34)**.

Menu: Power setting

All changes are saved automatically. Even after turning off and then again turn the device on, all the things will remain. If you want to save certain settings from being changed, you will find further information in chapter 7.8 "Saving and accessing the program" and chapter 8.9 „ Saving the program against changes"

8.8 Saving and accessing the program

Touching the "Program" field causes a selection of saved programs to appear.

To select the required program touch the relevant line. The program you are located in is highlighted in grey.

In the main menu, the number of the accessed program and the name are displayed


Nr		⊗
1	DR. MULLER	✓
2	Name 2	
3	Name 3	
4	Name 4	
5	Name 5	
6	Name 6	
7	Name 7	
8	Name 8	

PROGRAM1 DR MULLER

To change the program name, touch the field with the name.

In our example, this is: DR MULLER. A window opens with a keypad. The maximum symbol length is 15 characters.

To save the name, press "Enter", and to quit without saving, press "Esc".

Every change of the settings is automatically saved in the program you are currently in.

In order to save individual settings in a new program, initially select a free storage location. Settings and the name of the program can then be changed at will.

8.9 Saving the program against changes

By touching the lock on the display the current program will be protected against changes. Nevertheless it is possible to change settings during surgery. Once the device is switched off, the settings will be restored at the time of activation the lock,

8.10 Switching off

After terminating the operation, switch off the system at the main switch and pull out the plug. After switching off the system, all electrodes and forceps must be disconnected from the cable and the electrode cable from the unit.

8.11 Security functions

Microprocessor tracking of the safety functions

In order to minimize the risks associated with high frequency electrosurgery, the system is equipped with a number of safety functions. The system's microprocessor control and sensor technology allow a number of monitoring functions to be implemented. The functional states as well as the states of the fault in the generator are signaled visually and tonally, whereby the HF current is switched off immediately if necessary.

The HRF 100 has an automatic, self-diagnosing control

Monitoring the wrong dosage

By following the output energy and power, the HBS argon module constantly compares the permissible voltage limits and the values selected on the control panel. If impermissible deviations are reached, the radiated HF energy is interrupted and an error message is output.

8.12 Error handling

In the event of an error, a message appears on the display.

The following error messages can occur:

Error	No.	Cause of failure	Error management
Err 1 Total Failure	1	"Power timeout error"	The supply voltage was not reached. Check the controller and the power pack, replace if necessary
	10	"Power Error"	The difference between the display and the measured values of the monitoring is greater than 20%. Check the controller, the interface and the power pack, replace if necessary.
Err 2 Modules Failure	2	"NE electrode error"	Check the connection on the device and the electrode on the patient, if necessary replace the NE and / or connection cable.
	6	"Mono - Cut hand /FS1 Error"	Check CUT connections, foot switch settings, check settings.
	7	"Mono – Coag (Bi) hand /FS1 Error"	Check COAG connections, foot switch settings, check settings.
Err 03 PS Failure	4	"Power Supply 1 Error"	The block voltage required to correct the power factor is missing. Check the power connection, check the power supply unit, replace if necessary.
	5	"Power Supply 2 Error"	There is no voltage required from the block for regulating the power voltage for the power generator when idling. Check the power connection, check the power supply unit, replace if necessary..
Err 5 CPU Failure	11	"EEROR in CPU ROM"	Incorrect check sum of the ROM memory from the control processor. Check controller, replace if necessary
	12	"EEROR in CPU RAM"	Incorrect check sum of the RAM memory from the control processor. Check controller, replace if necessary.

If a different numerical code is displayed or the device has other recognizable defects, please contact service@hebumedical.de.

9 Guideline value for power settings

When performing the power settings, the user should generally take into account the cutting speed, geometry of the electrode and properties of the tissue.

Monopolar cutting:

	Recommended settings
Smooth cuts with thin needle electrode	Step 1 – 4
Cutting with knife electrode, lancet or wire loop	Step 2 – 6
Cutting with wire loop	Step 2 – 6
Cutting in fatty tissue with lancet or knife electrode	Step 3 – 8

Monopolar coagulation:

	Recommended settings
Coagulation of large surfaces	Step 6
Coagulation of small surfaces	Step 4

Bipolar application:

	Recommended settings
Coagulation	Step 1 – 3

The recommended settings serve as guideline values only and must be adjusted to the relevant situation.
In the event of unfavourable conditions, these values may deviate!

10 Protective measures and warnings

Please observe also the instructions provided on working with the neutral electrode in chapter 8.6 the operating instructions with the accessories and the applicable safety measures!

Before any procedure, check the cables, accessories (handles, electrodes, forceps) and foot switches for any signs of visible irregularity. Do not use cables or instruments with broken or defective insulation. These can result in burns if they come into contact with the skin!

Before any procedure, check the cables, accessories (handles, electrodes, forceps) and foot switches for any signs of visible irregularity. Do not use cables or instruments with broken or defective insulation. These can result in burns if they come into contact with the skin!

	The device must not be used in areas where there is a risk of explosion.
	The device may only be repaired by HEBU or a representative authorized by HEBU.
	It must be ensured that the device is only operated by people who are familiar with the operating instructions.
	The device, including the connecting cables, should be checked before each use.
	Faults must be eliminated in good time
	All plug-in connections may only be carried out with compatible accessories.
	Use the greatest care when working with the neutral conductor!
	Protect the patient from contact with earthed metallic elements (OP table, infusion stand etc.) during use.
	Avoid contact between supply cables and patients, other cables or instruments! Ensure that there are no loops in the cables of handles and electrodes or fixture to the OP table.

	Avoid small-area skin-to-skin contact in the patient (e.g. a finger making contact with the thigh).
	Avoid contact between HF instruments and non-insulated instruments..
	Avoid penetration of liquids between the patient's body and the neutral electrode.
	The output power should not be set higher than is absolutely necessary for the procedure.
	Do not touch the patient during activation of the coagulation or cutting current. Should this be necessary, a piece of fabric can be used for insulation purposes.
	In patients connected to monitoring devices (ECG), the monitoring electrodes must be attached at the greatest possible distance to the application areas of the electrosurgical electrodes.
	The use of highly flammable narcotic substances and oxidizing gases (N20 or O2) should be avoided. Should this be necessary, take particular precautions..
	Exclusively non-flammable disinfectants must be used. If alcohol-based solutions are required, wait before performing surgery until these have completely evaporated. Bear in mind that disinfectants can flow down and collect in recesses of the body. A spark from the active electrode can cause combustion of the fluid in normal use.
	Sparks from the active electrode can set bandages (textiles), metabolic gases and combustible fluids alight.
	In the case of procedures involving high-frequency current applied to body parts with a minimal cross-section, the bipolar technique should be used in order to prevent the risk of unintentional coagulation.
	In the case of procedures on body parts with a close link to the remainder of the body, the bipolar technique should be used for safety reasons.
	The frequently used indirect coagulation mode should only be used with well insulated forceps. We recommend forceps with specially insulated grip areas (ask your HEBUmedical contact)
	Surgical gloves do not afford the surgeon sufficient protection from burns!
 	During surgery, pacemakers are at risk of destruction. We recommend consulting the treating cardiologist and/or an authorized representative of the pacemaker manufacturer. As a rule, in such cases bipolar application technology is advisable. Electrosurgical procedures are prohibited on outpatients with pacemakers.

	<p>A marked drop in power at a normal setting can be an indication of incorrect application of the neutral electrode.</p>
	<p>Soiled electrodes can bring about a drop in the quality of the unit. This applies particularly to gentle and bipolar coagulation. To prevent the tissue adhering to the active electrode (monopolar or bipolar), it is advisable to wait for a moment after coagulation of the tissue and to interrupt the current flow to allow the active tip to cool before separating the electrode from the tissue.</p>
	<p>High-frequency creep currents can cause burns in areas distant from the application site through conductive connections.</p>
	<p>To remove the cable from the unit, always hold the plug! Never pull on the cable, as this can result in damage to the cable insulators, burns on the skin or also fire.</p>
	<p>Instruments which are temporarily not required during the procedure must be placed away from the patient on the instrument table or a similar surface.</p>
	<p>Reusable electrosurgical instruments which are connected to the unit are generally not delivered in a sterile condition!</p>
	<p>In a “medical” environment (environment in which flammable mixtures with air can cause an explosion with normal ignition energy) only use class AP footswitches</p>

■ **Electromagnetic safety**

In the development of electrosurgical units, HEBUmedical attaches particular importance to the stringent electromagnetic emission regulations. Consequently, the solutions we supply have a correspondingly minimal emission level for the required output. Measurements conform the high level of electromagnetic safety offered by HEBUmedical electrosurgical units.

Under typical working conditions based on an 8-hour day, a field occurs at a distance of 5 to 15 cm from the cables. At a distance of 20-40 cm, the value drops well below the prescribed limit.

Electromagnetic fields occur primarily around the cable.
 The electrosurgical unit itself does not constitute a significant emitting element.

11 Review of the technical status

Before any procedure, check the cables, accessories (handles, electrodes, forceps) and foot switches for any signs of visible irregularity. Do not use cables or instruments with broken or defective insulation. These can result in burns if they come into contact with the skin!

11.1 Mechanical faults

If damage occurs to switches, connecting sockets, houses or the touchscreen, or if the unit is dropped on the floor, please contact your authorized dealer for advice on further use.

The dealer can arrange for a detailed technical check to be performed.

The manufacturer accepts no liability for calibration work or repairs to the unit carried out by the user.

11.2 Recurring safety inspections

A safety inspection of the electrosurgical system is required once a year. The manufacturer only permits the use of units for which an up-to-date inspection by suitably trained personnel can be verified.

All information relating to the inspections (date, inspecting staff member, result etc.) must be recorded in the unit's log book.

11.3 Servicing and manufacturer address

Should you require the instructions for use in paper form, please use the contact details below. The instructions for use in paper form will be made available to you within seven calendar days of receipt of the request.

Alternatively, you can print out the electronic instructions for use yourself.

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen / Germany

Tel. +49 7461 94 71 - 0

Fax +49 7461 94 71 - 22

eMail: service@HEBUmedical.de

Web: www.HEBUmedical.de

12 System maintenance work

12.1 Cleaning

When cleaning the device:

Always disconnect the device from the mains before cleaning.

Before connecting the device to the mains, make sure that all flammable cleaning agents have completely evaporated.

Follow the local hygiene instructions.

- Clean all surfaces of the device (including the front) with cleaning agents and disinfectants without alcohol.
- Follow the manufacturer's instructions for the disinfectants.
- Wipe off the agent with a sponge or cloth dampened in clean water.
- Dry the device with a clean, lint-free cloth.
- Do not clean or disinfect surfaces with flammable or explosive products. If their use cannot be avoided, flammable or explosive products must be completely vaporized before switching on the device.
- Propanol and ethanol attack surfaces.
- Do not immerse the device in liquid.
- Do not sterilize the device with steam, ethyloxide or other agents

12.2 Fuses

The fuse sockets are located on the back of the unit. The unit works with two fuses of 6.3 A each.

Exchanging the fuses:

- Disconnect the unit from the mains
- Release the fuse from the retainer using a flat blade screwdriver
- Exchange the fuse

After an exchange of fuses, the unit must be checked by suitably qualified personnel.

12.3 Transport

Prevent mechanical damage and moisture.

If the unit has been stored for a long period in a cold environment, wait before removing the protective packaging until the unit has reached room temperature.

During transportation, the standard safety measures additionally apply.

Only ever transport the unit in its original packaging.

12.4 Environmental directives

From the implementation date of the European directive 2002/96/EU as national legislation, the following regulations apply:

- Electrical and electronic units must not be disposed of with normal household waste.
- The consumer is required by law to return electrical and electronic devices at the end of their service life to designated public collection points or to the place of purchase.

Disposal details are set out by the applicable local legislation in the relevant country. The symbol on the product, the instructions for use or the packaging indicates the relevant applicable legislation. By sending end-of-life devices for reuse, material recycling or other forms of utilization you are making a major contribution to environmental protection.

Please also send the unit packaging for environmentally responsible recycling after the end of the product's life.

Warranty card

We grant on all electrical medical products

We grant a warranty in accordance with the statutory and country-specific regulations on all electrical medical products (verification required in the form of an invoice or delivery note). The warranty period is two years. In the event of material or production errors, we provide free repair of the unit. Any attempt to carry out independent repairs will render the warranty cover null and void.

Damage caused by incorrect handling, overloading or normal wear and tear is not covered under the terms of the warranty.

In the event of a warranty claim and resulting repair, please enclose a copy of the invoice or delivery note (or batch number) with the instrument.

Manufacturer and Aftersales Service Address:

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen
Tel. +49 7461 9471 - -0
Fax +49 7461 9471 -22
service@HEBUmedical.de

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