HF 9928-10

HF 9928-11

HF 9928-12

HF 9928-13

APC-Sonde / APC probe / Sonde APC / Sonda APC / Sonda APC

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- HF 9928-10**
- HF 9928-11**
- HF 9928-12**
- HF 9928-13**

Contents

1	Symbol descriptions	16
2	Introduction	17
3	Storage.....	17
4	Scope	17
4.1	Intended Use	18
4.2	Contraindications	18
4.3	Maximum electrical capacity	18
5	Application remarks	19
5.1	Connecting the probe	20
6	Warnings	21
7	Safety remarks	23
8	Preparation and Sterilization.....	24
9	Inspection.....	25
9.1	Functional check.....	25
10	Warranty/Repair	25
11	Servicing and manufacturer address	26

1 Symbol descriptions

Symbol	Definition
	CE-labelling
	Attention
	Manufacturer
	Lot-description
	Reference code
	Medical device / FDA Prescription device
	Medical device
	Keep away from sunlight
	Dry storage required
	(Electronic) instruction for use
	No intervention on patients with cardiac pacemakers.
	Temperature
	sterile

Do not reuse

2 Introduction

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following. In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

3 Storage

Store the flexible APC probes and their connecting cable in a cool and dry location. Protect from the effects of sunlight. Medical products should never be stored or operated in dusty environments on principle.


4 Scope

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

APC probes are used in the following areas:

- General medicine
- Gynecology
- liver surgery
- Traumatology
- Cardiac and Thoracic Surgery
- ENT

Example of products:

HEBU APC probe [1] HF 9928-10

Connecting cable [2] HF 9928-00

4.1 Intended Use

HEBU disposable APC Probe: The flexible probes for argon plasma coagulation (APC probes) are used for argon gas-assisted coagulation and are used in oesophago-, gastro-, broncho-, colono- and rectoscopically assisted interventions in the gastrointestinal tract or in the tracheo-bronchial space. The APC probes can be used in conjunction with electrosurgical equipment (high frequency generators) and an argon unit with an output voltage of max.4 kVs. They are connected to the unit with the reusable and steam sterilizable cable. The product is supplied sterile and is intended for single use. The activation of the gas supply as well as the HF-current takes place via the footswitch of the high frequency generator.

4.2 Contraindications

There are no known contraindications.

4.3 Maximum electrical capacity

Maximum electrical capacity: 4kVs.

5 Application remarks

- HEBU APC probes are generally operated using a spray coagulation current mode. Select this current mode and always operate using the blue pedal!
- To "ignite" the argon plasma, approach the tissue with the distal probe opening to a distance of around 1-5 mm.
During the argon-supported coagulation process, the distance can be increased to appr. 2 - 2.5 mm before the plasma is extinguished.

Flexible endoscopy and laparoscopy applications:

The diameter of the endoscope working channel should be at least 0.3 mm bigger than the specified probe diameter.

The maximum flow rate of 2 l/min must not be exceeded!

During endoscopic application, flow rates of 0.3 to 0.8 l/min are recommended at a spray output rating of appr. 30 W.

For the individual probe sizes, this means:

HEBU article no.	Diameter	Length	Flow rate
HF 9928-10	1,5 mm	1,5 m	0,3 - 0,4 l/min
HF 9928-12	2,3 mm	2,2 m	0,4 - 0,6 l/min
HF 9928-11		1,0m	
HF 9928-13	3,2 mm	2,2 m	0,5 - 0,8 l/min

In order to ensure that a critical intra-abdominal pressure level is not exceeded, the trocar valve must be opened during laparoscopic applications. If an insufflator pressure alarm gives off a signal, the argon application must be stopped until the intra-abdominal pressure has dropped back to below the critical level.

Whichever method is used to monitor intra-abdominal pressure by an insufflator, independent and continuous monitoring of intra-abdominal pressure by the OP team is a requirement. The use of a pressure release-controlled insufflator is recommended.

5.1 Connecting the probe

Initially, ensure that the electrical surgical device and the (integrated) argon unit are ready for operation. Ensure that the gas cylinder valve is open and that there is an adequate supply of gas (> 30 bar).

1. Remove the connecting cable (p. 3, [2]) and HEBU APC probe (page 3, [1]) from the sterilization container under sterile conditions.
2. Carry out a visual inspection of the instrument for possible damage to the insulation. Do not use damaged instruments! HEBUmedical warns users not to attempt to repair or patch up damaged areas. The instrument is operated at high coagulation voltage levels!
3. Check whether the working channel of the endoscope is suitable for the required probe diameter (3.2 mm, 2.3 mm or 1.5 mm).
4. Connect the three-pin plug of the connecting cable to the HEBU electrical surgery device.
5. Insert the required probe with plug into the coupling at the end of the connecting cable. Ensure that the plug is pushed into the coupling as far as it will go.
6. Establish the argon gas connection to the electrical surgery device / the argon unit. To do this, connect the female Luer lock connector of the connecting cable to the male Luer lock gas outlet of the electrical surgery device / the argon beamer. If required, an intermediate sterile filter can be connected.
7. Activate the function “**C**” at the electrical surgery device / the argon unit in order to flood the gas lumina of the connecting cable and the probe with argon gas.
8. Actuate the blue foot pedal and check the system function (activate the argon beamer and generator).

The probe is ready for operation.

6 Warnings

The products may only be used in combination with HEBU accessories and only by clinically instructed and qualified personnel. Electrodes must not be twisted, as this may cause damage to the material and the electrode might break. Power output of the HF unit: Always select the lowest possible output of the HF unit for any procedure. If the coagulation capacity of the electrode is below normal, do not increase the output of the HF unit without a thorough preliminary check. In particular, the following must be checked: correct attachment of all HF cables and plugs, correct activation of current by the finger - or footswitches, undamaged insulation of the HF cables and instrument, and check that the distal end of the electrode is dean and undamaged. Before operating, check that the signal transmission of the monitors works without interference (e.g. without noise); this check is made by pressing the yellow button (CUT) and the blue button (COAG).

Always grasp the plug when pulling the HF cable. **Never pull on the HF cable**, as this could result in damage. Cables with defective/brittle insulation or a break in the electrical cable (due to severe kinking or crushing of the cable) can cause burns to the user/patient or even cause fire. Active HF instruments and HF handles must not be placed on the patient to prevent personal injury in the event of accidental HF activation.

To clean the active electrodes, the "Auto Start" mode (only for bipolar application) must be deactivated or the instrument must be disconnected from the generator. The tip area may still be hot immediately after activation, i.e. there is a risk of burns.

Endogenous burns are those caused by high current density in the patient's tissue. Possible reasons are: The patient is inadvertently positioned in contact with electrically conductive parts. Direct contact between skin areas and HF cables may lead to electrical capacitance which in turn may cause burns.

Exogenous burns are those caused by heat of ignited fluids or gases. They may also be caused by explosions. Possible reasons include: Ignition of skin cleaning agents and disinfectants, ignition or narcotic gases etc.

Cardiac pacemakers may be damaged by electrosurgical current. Consult a cardiologist before operating. Never use electrosurgical current on outpatients fitted with pacemakers.

7 Safety remarks

	Do not reuse products which are intended for a single use only!
	Never place the gas outlet nozzle of the HEBU APC probe directly on the tissue during gas activation! Danger of gas emphysema!
	Do not guide the gas outlet nozzle vertically over the tissue during activation. An angle of around 45° to the tissue has proven successful.
	Do not blast argon into the vessel system. Danger of embolism!
	Always hold the applicator end within the visual range of the endoscope during activation! Never activate the coagulation stream without visual inspection!
	Endogenous gas in the gastrointestinal tract can be combustible or explosive. Before any electrosurgical application, clear any such gases in suspected areas by flushing with argon gas (e.g. by activating the beamer on the HF output setting "0" or using the purge function "C/Argon Flow" of the beamer).
	Argon plasma can ignite combustible materials (e.g. the plastics used in the endoscope or tubes). This can result in a hazard, particularly in oxygen-enriched environments!
	The distal end of the probe must project by at least 10 mm from the working channel of the endoscope. There are marking rings provided at the distal probe end for orientation purposes.
	Do not introduce oxygen or other combustible gases or fluids before and during application in the tracheo-bronchial system. During protracted beamer applications, carry out alternating ventilation!
	In the event of laparoscopic applications, the argon gas flow brings about an increase in intra-abdominal pressure. The use of electronic pressure release-controlled CO ₂ insufflators is recommended. During activation, monitor intra-abdominal pressure and where necessary release pressure through the trocar!
	For endoscopic applications with video endoscopes, do not direct the argon gas jet towards the camera chip.
	Protect the probe from mechanical overloading. Ensure that the hose does not kink!
	As high HF voltage levels are required for gas ionization, pay particular attention to ensuring that the electrical insulation of cables and applicators is undamaged.

	Electrode wires which project out from the flexible tube at the distal end can result in mechanical damage to tissue
	As argon gas-supported application is a monopolar high-frequency form of surgical application, the safety rules for monopolar high-frequency surgery must be observed!
	Due to the manufacturing method used, quality 4.8 argon in gas cylinders is germ free. In order to guarantee freedom from germs, the cylinder valve and the gas outlet at the gas supply device must be disinfected prior to connection.
	Observe manufacturer's instructions for electrical load capacity!
	Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

8 Preparation and Sterilization

Flexible APC probe, single use are sterile packed and must not be cleaned and sterilized before use. They must be disposed of after single use.

Optionally, it is also possible to purchase the flexible APC probe, single use non-sterile.

The non-sterile HFIT®-T needle electrodes must be cleaned and sterilized before use as follows:

Rinse APC probe thoroughly with demineralized water.

	On using the fractionated vacuum method, sterilization must be performed with at least 134°C (USA 132° C) with a minimum dwell period of 3 minutes . Vacuum drying must then be carried out for at least 20 minutes.
---	---

Do not sterilize flexible APC probe, single use in hot air sterilizers.

This product is a disposable product and must **NOT** be reused!

9 Inspection

The products must be visually inspected for protein residues and other contamination, paying particular attention to cavities, blocks, inclusions, pipes, and other inaccessible areas..

Carry out the functional check mentioned above. Products with stains, which are blunt, bent, no longer function or which are otherwise damaged must be segregated!

To help identify faulty products that need to be sorted out, we recommend the brochure "Instrument Reprocessing" from the Working Group "Instrumenten Aufbereitung". This includes in particular Chapter 8 " Checks and Care" and Chapter 12 " Surface Changes: Deposits, Discoloration, Corrosion, Aging, Swelling and Stress Cracks".

9.1 Functional check

A newly purchased product must be subjected to a thorough visual and function check after its delivery and before each use. If a new product has externally visible defects or does not operate according to the instructions, contact the manufacturer or his responsible representative immediately.

When used as intended, these products are more or less subject to wear, depending on the intensity of use. This kind of wear is caused by technical factors and cannot be avoided. Therefore carry out a thorough visual examination and a function test prior to each application of these products. Carry out all plug connections only with the compatible articles and plug in the plug completely.

Products must be checked for irregularities. Paying attention to cracks, fractures and the occurrence of corrosion. Carry out a function check in accordance with the intended application of the product.

Defective products must not be used!

10 Warranty/Repair

Our products are made of high-quality materials and are carefully checked before delivery. However, even when used as intended, they are subject to wear to a greater or lesser extent depending on the intensity of use.

This wear is due to technical reasons and is unavoidable.

However, if faults occur that are not related to wear, please contact our customer service department.

Defective products may no longer be used and must have undergone the entire preparatory process before being returned.

11 Servicing and manufacturer address

Should you require the instructions for use in paper form, please use the contact details below. The instructions for use in paper form will be made available to you within seven calendar days of receipt of the request.

Alternatively, you can print out the electronic instructions for use yourself.

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HEBUmedical GmbH