HF 9530-02

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HF 9530-02

Contents

- 1 SYMBOL DESCRIPTIONS14
- 2 INTRODUCTION15
- 3 SCOPE16
 - 3.1 INTENDED USE16
 - 3.2 CONTRAINDICATION16
- 4 WARNING16
- 5 PRODUCT DESCRIPTION.....17
- 6 COMMISSIONING AND USE17
 - 6.1 CONNECTION TO GAS CYLINDER.....17
 - 6.2 CONNECTION TO HF-DEVICE.....18
 - 6.3 USE19
 - 6.3.1 Opening19
 - 6.3.2 Closing.....19
 - 6.3.3 Adjusting the mechanical style pressure switch adjustment screw19
- 7 CLEANING19
- 8 MAINTENANCE.....19
 - 8.1 CHANGING GASKET RING20
 - 8.2 CHANGING FILTER SCREW20
 - 8.3 FUNCTIONAL CHECK.....20
- 9 STORAGE21
- 10 WARRANTY / REPAIR21
- 11 TECHNICAL DATA.....21
 - 11.1 PRESSURE REDUCER21
 - 11.2 MECHANICAL STYLE PRESSURE SWITCH.....21
- 12 SERVICING AND MANUFACTURER ADDRESS.....22

1 Symbol descriptions

Please read these instructions carefully and keep them for future reference. Improper handling can cause personal injury or damage the product.

Symbol	Definition	Danger	Consequence
	DANGER!	Immediate danger for people	Death or serious injuries
	WARNING!	Possible danger for people or objects	Health damages or severe material damages
	ATTENTION!	Possible danger to objects	Material damages

Additional symbol descriptions:

Symbol	Definition
	CE-labelling
	Attention
	Manufacturer
	Validated Parameters
	Lot-description
	Reference code
	Medical device / FDA Prescription device

	<p>Medical device</p>
	<p>Non sterile</p>
	<p>Keep away from sunlight</p>
	<p>Dry storage required</p>
 <p>Hinweis auf eIFU</p>	<p>(Electronic) instruction for use</p>
	<p>No intervention on patients with cardiac pacemakers.</p>
	<p>Temperature</p>
	<p>WEEE-labelling</p>
	<p>Green dot</p>

2 Introduction

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

You are now the owner of a high-quality product whose use and correct handling are described in the following. In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

3 Scope

This product is designed to power HBS-Xtouch devices with argon gas. It reduces the gas pressure of the gas source to the required operating pressure and maintains it constant.

3.1 Intended Use

Argon pressure reducer: This product is designed to power Argon gas HBS x-touch devices. It reduces the incoming gas pressure of the gas source to the required operating pressure and maintains it constant.

3.2 Contraindication

No contraindications are known.

4 Warning

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

	Incorrect use of the pressure reducer can cause serious damage.
	Avoid mechanical overstressing of the device beyond the structural design, this can lead to breakage and loss of function!
	Run the pressure reducer only at the intended pressure and with the intended gas. The use of improper gas or at improper pressure might cause risks of explosion.
	Protect the pressure reducer from dust, oil and other sources of dirt.
	Do not use the pressure reducer if it is not in perfect working conditions.
	When you draw gas, the cylinder must be placed upright and protected from falling.
	Danger of explosion. At startup and also during use, hands and tools must be free from oil and grease.
	Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

5 Product description

- High pressure gauge (1)
- Mechanical style pressure switch (2)
- Pressure adjusting screw (3)
- Relief valve (4)
- Outlet connection (5)
- Body (6)
- Gasket ring (7)
- Coupling nut (8)
- Inlet connection (9)
- Filter screw (10)
- Adjusting screw mechanical style pressure switch (11)


Figure 1

Figure 2

6 Commissioning and use

6.1 Connection to gas cylinder

The pressure reducer is constructed for gas cylinder connections in line with DIN 477.

- At every gas cylinder change, initial check the gasket ring (7) on intactness.
- Before connecting the pressure reducer to the gas cylinder, make sure the openings of the gas cylinder and the pressure reducer are free from dirt. Otherwise blow out with compressed-air.
- Turn the pressure adjusting screw (6) anti-clockwise as far as it will go. Do the same with the valve of the gas cylinder.
- Put the inlet connection (9) of the pressure reducer on the opening of the gas cylinder. Align the high pressure gauge (1) upright and close the coupling nut (8) manually. If necessary again align the high pressure gauge.

- Tighten the coupling nut by using a 30 mm flat wrench.


Stay away from heat sources and easily inflammable materials during the commissioning.

Do not direct the opening of the gas cylinder on yourself or other persons.

6.2 Connection to HF-Device

Ports of the HF-Device

- gas pressure measure line 2 (12)
- gas connection 2 (13)
- gas pressure measure line 1 (14)
- gas connection 1 (15)

Figure 3

- Plug in the mechanical style pressure switch (3) leaving cable in gas measure line 1 or 2 (14 or 12).
- Turn the lock clockwise until it is engaged. You can hear that at the click-sound.
- Plug in the outlet connection (5) leaving tube with its adapter at the end to gas connection 1 or 2 (15 or 13). The lock engages by itself (Figure 4 and 5).


Figure 4

Figure 5

6.3 Use

6.3.1 Opening

- Slowly open the gas cylinder valve. The high pressure gauge (1) will show you the cylinder pressure.
- Slowly open the pressure adjusting screw (3) as far as it will go.

Opening the gas cylinder valve too quickly might cause gauges malfunction.

6.3.2 Closing

- Close the gas cylinder valve.
- Release the gas out of the pressure reducer.
- Turn the pressure adjusting screw (3) anti-clockwise till it is completely closed.

6.3.3 Adjusting the mechanical style pressure switch adjustment screw

The mechanical style pressure switch fulfills the purpose to provide the control system of the HF-device with information, whether the required pressure is on or not. For that you can set a minimum pressure with the aid of the adjustment screw. Pay attention that the minimum pressure is always set under 4 bar, otherwise the gas supply will be inhibited.

To set the minimum pressure, perform the following instructions:

- Turn the adjustment screw (11) clockwise with the aid of a slotted screwdriver, to increase the pressure.
- Turn the adjustment screw (11) anti clockwise with the aid of a slotted screwdriver, to decrease the pressure.
- The adjusted pressure can be seen on the scale.

7 Cleaning

In case you need to clean the pressure reducer use a dry, dust-free cloth to wipe the item. For stubborn soiling use a damp cloth.

Avoid water inside the pressure reducer. It might cause risks of explosion.

8 Maintenance

Hands and tools must be free from oil and grease during maintenance

8.1 Changing gasket ring

At a high level of abrasion perform the following instructions:

- Lever the gasket ring with the aid of a slotted screwdriver and take it off.
- Pull the fresh gasket ring without any aid of tools on its place.

8.2 Changing filter screw

At a high level of soiling, but at least every 2nd year perform the following instructions:

- Release the filter screw with the aid of a 6 mm allen wrench.
- Hold the pressure reducer so that the inlet connection (9) is facing to the ground. This way it is prevented that dirt enters the inside.
- Turn out manually the filter screw.
- Turn in manually the fresh filter screw.
- Tighten the filter screw with the aid of a 6 mm allen wrench.

8.3 Functional check

A newly purchased product must be subjected to a thorough visual and function check after its delivery and before each use. If a new product has externally visible defects or does not operate according to the instructions, contact the manufacturer or his responsible representative immediately.

When used as intended, these products are more or less subject to wear, depending on the intensity of use. This kind of wear is caused by technical factors and cannot be avoided. Therefore carry out a thorough visual examination and a function test prior to each application of these products. Carry out all plug connections only with the compatible articles and plug in the plug completely.

Products must be checked for irregularities. Paying attention to cracks, fractures and the occurrence of corrosion. Carry out a function check in accordance with the intended application of the product.

Defective products must not be used and must have undergone the complete preparatory treatment process again before being returned.

9 Storage

The pressure reducer must be stored dry and dust-free. The surrounding air must be non condensing.

It is recommended to keep the wrapping, so the item can be storage in it when required.

10 Warranty / Repair

Our products are manufactured from high-grade materials and carefully checked prior to dispatch. However, even if used properly in accordance with their intended purpose they are subject to a greater or lesser degree of wear depending on their intensity of use.

This wear is technically induced and unavoidable.

Should faults occur independently of wear, please contact our customer services. Defective products should no longer be used.

They must undergo the complete preparatory treatment process before being returned.

11 Technical data

11.1 Pressure reducer

Inlet pressure P1 [bar]	max. 230
Outlet pressure P2 [bar]	4
Delivery flow Q1 [m ³ /h]	< 1
Operating temperature [°C]	-20 bis +60

11.2 Mechanical style pressure switch

Model	IS 1000		
Proof pressure [bar]	10		
Max. pressure [bar]	7		
Operating pressure range [bar]	1 bis 6		
Hysteresis [bar]	0,8		
Ambient and gas temperature [°C]	-5 bis +60		
Max. contact capacity	2 VAC AC / 2 W DC		
Voltage	24 V AC/DC	48 V AC/DC	100 V AC/DC
Max. operating current	50 mA	40 mA	20 mA

12 Servicing and manufacturer address

Should you require the instructions for use in paper form, please use the contact details below. The instructions for use in paper form will be made available to you within seven calendar days of receipt of the request.

Alternatively, you can print out the electronic instructions for use yourself.

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen / Germany
Tel. +49 7461 94 71 - 0
Fax +49 7461 94 71 - 22
eMail: service@HEBUmedical.de
Web: www.HEBUmedical.de

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**Our vast range of products covers over 10 000 various instruments.
They can be found in Internet or request our catalog.**

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