*HF 9581-20*  
*HF 9581-37*

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**HF 9581-20**  
**HF 9581-37**

# Contents

1	Symbol descriptions.....	22
2	Introduction.....	23
3	Scope .....	23
3.1	Intended Use .....	23
3.2	Contraindication.....	23
3.3	Maximum electrical capacity .....	23
4	Application.....	24
4.1	Functional check in the sterile operation area .....	24
4.2	Operating.....	24
5	Assembly and disassembly.....	25
5.1	Assembly of blade.....	25
5.2	Cleaning .....	26
5.3	Blade assembly .....	26
6	Warnings .....	27
7	Safety Remarks .....	28
8	Scope of delivery .....	30
9	Preparation.....	30
9.1	Reutilization restrictions .....	30
9.2	Information on instrument preparation.....	30
9.3	Preparation at the place of use .....	31
9.4	Ultrasound bath (optional).....	31
9.5	Manual cleaning.....	31
9.6	Mechanical cleaning .....	31
9.7	Drying .....	34
10	Maintenance, inspection .....	34
10.1	Functional check.....	34
11	Sterilization.....	35
11.1	Packaging.....	36
12	Lifetime.....	36
13	Storage.....	36
14	Warranty/Repair.....	36
15	Servicing and manufacturer address.....	36

# 1 Symbol descriptions

Symbol	Definition
	CE-labelling
	Attention
	Validated Parameters
	Manufacturer
	Lot-description
	Reference code
	Medical device / FDA Prescription device
	Medical device
	Non sterile
	Keep away from sunlight
	Dry storage required
	(Electronic) instruction for use
	No intervention on patients with cardiac pacemakers.
	Temperature

## 2 Introduction

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following. In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

## 3 Scope

The use of HF bipolar clamps is used in HF surgery for coagulation and tissue cutting. In addition, for the permanent closure of veins, arteries and tissue bundles with subsequent separation.

### 3.1 Intended Use

A monopolar or bipolar product covered with electrical insulation down to the working end. It is used in open surgery to coagulate vascular or tissue bundles for hemostasis. It is connected either to an electrosurgical diathermy unit or to an electrocautery generator that delivers the bipolar electrical current. It is typically designed as a scissors-type instrument with a long, curved, pincer-like mouth. It is made of metal and synthetic material in different sizes and shapes. It is a reusable product.

### 3.2 Contraindication

The RF bipolar clamps should not be used in conjunction with the central nervous or circulatory system. Tube ligation is not a safe application.


### 3.3 Maximum electrical capacity

Maximum electrical capacity: 0,3 kVp

## 4 Application

The vessel sealing instrument BiSeal® is intended to be used exclusively with units which are equipped with the BiSeal®-function from HEBU medical GmbH. Attention must be paid to the directions given in the manual instructions as delivered with the unit. The BiSeal® vessel sealing instrument is used for minimally invasive surgery, in particular for laparoscopy. The instrument can be inserted through the 10 mm trocar. Its purpose is to permanently seal veins, arteries and bundles of tissue, and thereafter implement the separation by cutting, without having to change instruments. Hereby vessels with a diameter of up to 7 mm can be permanently sealed. A cable connects the instrument to the HF-generator's bipolar sockets.

The vessel sealing instrument BiSeal® is available in working lengths of 200 mm (HF 9581-20) and 370 mm (HF 9581-37).

Due to interaction of mechanical pressure ( $p$ ) and high-frequency current ( $I$ ) a permanently and reliably sealed area of veins, arteries or tissue bundles is created. Hereby it is of significant importance that the pressure is being kept constant during the entire application period.

Instruments for coagulation may only be used by persons who received special training or instructions in the use of such. Vessels with a diameter of up to 7 mm can be permanently sealed. For large vessels it is recommended to make a double coagulation, one next to the other. After sealing, it can be cut in the middle.

### 4.1 Functional check in the sterile operation area

The mechanical function must be checked before begin of operation, and by laparoscopic applications before inserting into the trocar: Opening, closing, locking, unlocking, moving the blade.

The electric function of the instrument should be checked before each application, e.g. on a sterile swab soaked with NaCl-solution.

### 4.2 Operating

After correct assembly and connection to one of the bipolar sockets on the HF-generator, set to BiSeal® mode.

Then grasp the tissue which is to undergo treatment. Afterwards close the handle. The sealing current is activated by use of the foot switch.

The system ends the process by an automatic stop-function (automatic stop of flow of the high-frequency current, audible confirmation).

Following this, the sealed tissue can be cut with the integrated blade.


Before beginning the sealing process, make sure that no ureter, nerve cords, intestinal loops or any other structures are grasped unintentionally!

For large vessels it is recommended to make a double coagulation, one next to the other. After sealing, it can be cut in the middle.

## 5 Assembly and disassembly

### 5.1 Assembly of blade

1. Pull down safety latch and push halfway in


2. Remove blade with needle holder and dispose according to safety standards.

### Attention!


**Never** remove blade with fingers – danger of injury!!!

## 5.2 Cleaning

3. Clean upper and lower parts of blade duct

## 5.3 Blade assembly

4. Insert new blade until it is caught in. Only then remove protective cap.

## 6 Warnings

The products may only be used in combination with HEBU accessories and only by clinically instructed and qualified personnel. Electrodes must not be twisted, as this may cause damage to the material and the electrode might break. Power output of the HF unit: Always select the lowest possible output of the HF unit for any procedure. If the coagulation capacity of the electrode is below normal, do not increase the output of the HF unit without a thorough preliminary check. In particular, the following must be checked: correct attachment of all HF cables and plugs, correct activation of current by the finger - or footswitches, undamaged insulation of the HF cables and instrument, and check that the distal end of the electrode is dean and undamaged. Before operating, check that the signal transmission of the monitors works without interference (e.g. without noise); this check is made by pressing the yellow button (CUT) and the blue button (COAG).

To pull the HF cable, always touch the plug. Never pull on the HF cable, as this could result in damage. Cables with defective/brittle insulation or a break in the electrical cable (due to severe kinking or crushing of the cable) can cause burns to the user/patient or even cause fire. Active HF instruments and HF handles must not be placed on the patient to prevent personal injury in the event of accidental HF activation.

To clean the active electrodes, the "Auto Start" mode (only for bipolar application) must be deactivated or the instrument must be disconnected from the generator. The tip area may still be hot immediately after activation, i.e. there is a risk of burns.

**Endogenous burns** are those caused by high current density in the patient's tissue. Possible reasons are: The patient is inadvertently positioned in contact with electrically conductive parts. Direct contact between skin areas and HF cables may lead to electrical capacitance which in turn may cause burns.

**Exogenous burns** are those caused by heat of ignited fluids or gases. They may also be caused by explosions. Possible reasons include: Ignition of skin cleaning agents and disinfectants, ignition or narcotic gases etc.

**Cardiac pacemakers** may be damaged by electrosurgical current. Consult a cardiologist before operating. Never use electrosurgical current on outpatients fitted with pacemakers.

## 7 Safety Remarks

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

Wrong application and disregard of cautionary instructions can cause serious incidents. Therefore please inform yourself with the assistance of this user manual, as well as the instructions of the electrosurgical unit, about all functions and the safe handling of this instrument. Improper handling or maintenance, as well as inappropriate use can lead to premature wear and tear of the instrument.

	All instruments must be completely cleaned, disinfected and sterilized, before being put into use and before every other further application.
	It is very important to inspect each surgical instrument for visible damage and wear, e.g. cracks, fractures or fissures in the insulation before every use. In particular areas such as blades, points, locks, ratchets, latches and any moving parts, insulation and ceramic parts are to be carefully examined.
	Do not use damaged instruments.
	Do not use in the presence of flammable or explosive substances.
	The instrument should not be placed on the patient.
	Bear in mind the limitations of this treatment method! Do not attempt to seal vessels which are wider than the maximum diameter of 7 mm!
	Use the instrument solely according to its intended purpose.
	Coagulation of fallopian tubes with the BiSeal® instrument is <b>not</b> a safe method of sterilization!
	A functional check must be done before implementation.
	The blade carrier must be exchanged after each application. Never use the same blade carrier twice!
	When submerging the jaw parts into electrically conducting liquids, no safe vessel sealing is possible, due to the current being conducted by the liquid (shunt circuit)!
	Even after ending the coagulation process, the outer surfaces of the jaw parts can be very hot and cause unwanted burns. Any touching of tissue immediately after the end of the sealing process is to be avoided!

	Do not activate the blade carrier before the sealing of the grasped vessel is completed – danger of bleeding! Only after the electrosurgery unit has confirmed the proper vessel sealing process, the cut can be made.
	Between the branches, in the blade duct and in the joints, rest substances might get deposited, which can cause a deterioration of function. (Electric shunt circuit, blade carrier gets stuck). In such cases the instrument must be cleaned during applications by use of delivered cleaning tools, a brush or a swab, as well as sterilized water.
	During vessel sealing the blade must be retracted into the jaw parts, otherwise it can cause a short circuit and no sealing function can be effected.
	Start with a low power output and increase the output according to requirement. If the sealing function is poor or not at all effective, look for cause. Do not simply turn the output power up to maximum.
	The connection cable should be aligned in such a way that it does not touch the patient, nor has any contact with other cables. Protect the cable from mechanical damage (running over, squeezing, bending).
	Avoid mechanical overstressing of the instrument beyond the structural design, this can lead to breakage and deformation!
	Avoid the simultaneous use of HF surgery and electrically grounded electrodes and sensors of diagnostic devices!
	If the products are used on patients with transmissible spongiform encephalopathy or HIV infection, we decline any responsibility for their reuse.
	Observe manufacturer's instructions for electrical load capacity!
	Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

#### Additional directions for laparoscopic applications

	Only do vessel sealing when the jaw parts of the instrument are visible for the operator. Unintentional activation or motion of instrument can lead to injury of patient!
	The electrodes inside the jaw parts carrying current must have a safety distance of at least 2 cm to metallic end of trocar, when activated.
	Insert and withdraw the instrument into/from the trocar with closed jaw parts only!

## 8 Scope of delivery

The instrument must be checked for correct contents as well as for being undamaged immediately upon delivery. In case of damages occurred during transport, these are to be reported right away.

HF 9581-20 or HF 9581-37	Vessel sealing instrument BiSeal® 200 mm or 370 mm
HB 2310	Needle holder – for save removal of the blades
HF 9581-01	Disposable replacement blades
-	Manual Instruction

## 9 Preparation

The person in charge of preparatory treatment is responsible for ensuring that the treatment is duly carried out using the relevant equipment, materials and personnel in the treatment facility and so achieves the desired result. This necessitates validation and routine monitoring of the process used. We urge you to take note of the national regulations dealing with instrument preparation.

### 9.1 Reutilization restrictions

Frequent repeat preparatory treatment has minimal effects on the product. The end of the product life is normally determined by wear and damage due to use

### 9.2 Information on instrument preparation

- Use cleaning and/or disinfection agents with a pH-value within 9-10. Please observe manufacturer instructions regarding dosage, exposure time and renewal of solutions.
- Do **not** use hard brushes (e.g. metal brushes or metal sponges) or coarse abrasive cleaners.
- Never leave instruments in cleaning or disinfection agents for longer than the specified time.
- Only used demineralized water for rinsing.
- Rinse and dry carefully through channels and pipes.
- Sensitive instruments must be cleaned in a storage or clamping fixture.
- Observe manufacturer instructions of cleaning – and sterilizing equipment.

### 9.3 Preparation at the place of use

Directly after using remove coarse dirt of products and rinse out the working cannulas. Do not use fixing agents or hot water (> 40°C), as this results in residues becoming fixed and can affect the success of the subsequent cleaning operation

Dismantle and/or open products as far as possible. Within short time after use clean the products for reducing a drying of the residues.

This enables an easier cleaning. If products come into contact with corroding medicines or cleaning agents, wash these up with water immediately after use.

Longer drying times, e.g. for dry disposal are not validated and therefore not recommended.

 The drying time during validation was 1 hour.

### 9.4 Ultrasound bath (optional)

All products must be opened, dismantled and any cavities rinsed through.

Place products in the screen basket in such a way that overlaps and contact between products are avoided. Add cleaning agent to the water and adjust the temperature of the solution in line with the cleaning agent manufacturer's instructions.

The cleaning in the ultrasound bath should be at 35-40 kHz, 5 minutes at least.

 To validate cleaning in an ultrasonic bath, the test items were ultrasonically treated in Neodisher Mediclean forte 0,5 % for 5 minutes.

Subsequently rinse products including all cavities before cleaning and disinfection.

Medicine products which are possessing a bad ultrasound transmission, e.g. soft material are not usable for the ultrasound bath.

### 9.5 Manual cleaning

Since mechanical processes can be standardized, reproduced and therefore validated, mechanical cleaning/disinfection should be preferred to manual processes.

Manual cleaning and disinfection process is not validated and therefore needs to be validated additionally by the end user.

### 9.6 Mechanical cleaning

On the basis of international standards (EN ISO 15883) and national directives, only validated machine cleaning and disinfection methods may be used. For the mechanical cleaning we recommend a standard program.



Only completely demineralized water should be used for cleaning, neutralisation and rinsing, in accordance with the „Guidance Compiled by the DGKH (Germany Society for Hospital Hygiene), DGSV (German Society for Sterile Supply) and AKI (Working Group on Instrument Reprocessing) for the Validation and Routine Monitoring of Automated Cleaning and Thermal Disinfection Processes for Medical Devices as well as Advice on Selecting Washer-Disinfectors“ (which refers to DIN EN ISO 15883-1 Point 6.4.2)

Flexible (complex) products with invisible surfaces have to be pre-cleaned manually before mechanical cleaning.


We recommend for all push shafts, -and pipe shafts instruments and instruments whose surfaces are on top of each others during the cleaning (f.e. bone forceps and gouge forceps) a manually pre-cleaning for an optimal cleaning result without residues

Observe the following by loading:


- Place the dismantled/opened instruments securely in the tray. Products with openings and gaps have to be faced down with the opened side so that they can be cleaned and no water of the cleaning process is collecting inside them. If available use balanced devices for rinsing



- Place the instruments with joints in an opened position into the cleaning, and disinfection machine.

- If needed use an adapter for cleaning

- Do not overload trays, avoid creating any overlaps.

Preliminary rinsing (cold, if applicable fully demineralized water without additives) is followed by chemical. The chemical cleaning should take place at **40°C -60°C** for at **least 5 minutes**.

We recommend products with a **pH-value within 9-10**, e.g. Neodisher MediClean forte from Dr. Weigert. The cleaning agents used should be selected depending on the material and properties of the instruments and in accordance to national regulations: If there is a high chloride concentration in the water, pitting and tension crack corrosion can occur on the products. The occurrence of this type of corrosion is minimized by using alkaline cleaning agents and demineralized water. By adding an acid-based neutralization medium, the rinsing off of alkaline cleaning agent residues is facilitated during the first intermediate rinsing process (warm or cold water).

In order to prevent the formation of deposits, it is advisable to use neutral cleaners where the water quality is unfavourable. After the second intermediate rinsing process, thermal disinfection takes place. The thermal disinfection should take place at temperatures of between **80 and 95°C**, with an **exposure time as outlined in EN ISO 15883**.

After the finishing of the program take the good out of the machine because corrosion can arise if the instrument remains in the machine.

<b>V Parameters used for the validation of preparation</b>	
<b>Pre-rinsing</b>	1 minute with cold tap water
<b>Cleaning</b>	Temperature: 55 °C
	Soaking Time: 5 minutes (worst case)
	Neodischer Mediclean forte 0,4% (worst case)
<b>Neutralization</b>	Temperature: cold DI water
	Soaking Time: 2 minutes
	Neodisher Z 0,1%
<b>Post-rinsing</b>	2 minutes with cold DI water
<b>Disinfection</b>	Temperature: 90 °C (A <sub>0</sub> 3000)
	Soaking Time: 5 minutes

## 9.7 Drying

Ensure adequate drying by the cleaning and disinfection device or using other suitable measures.

**V** Drying was omitted in the validation (worst case condition).

## 10 Maintenance, inspection

After cooling to room temperature, the products must be visually inspected for protein residues and other contamination, paying particular attention to cavities, blocks, inclusions, pipes, and other inaccessible areas. Instruments which are not free of residues must be returned for a complete retreatment process. To ensure that surgical instruments can be used according to their intended purpose after reprocessing, it is necessary to carry out a functional test after cleaning, disinfection and drying of the visual inspection and care measures. Carry out functional checks described in point 10.1.

Products with stains, which are blunt, bent, no longer function or which are otherwise damaged must be segregated!

To help identify faulty products that need to be sorted out, we recommend the brochure "Instrument Reprocessing" from the Working Group "Instrumenten Aufbereitung". This includes in particular Chapter 8 "Checks and Care" and Chapter 12 "Surface Changes: Deposits, Discoloration, Corrosion, Aging, Swelling and Stress Cracks".

### 10.1 Functional check

A newly purchased product must be subjected to a thorough visual and function check after its delivery and before each use. If a new product has externally visible defects or does not operate according to the instructions, contact the manufacturer or his responsible representative immediately.

When used as intended, these products are more or less subject to wear, depending on the intensity of use. This kind of wear is caused by technical factors and cannot be avoided. Therefore carry out a thorough visual examination and a function test prior to each application of these products. Carry out all plug connections only with the compatible articles and plug in the plug completely.

Products must be checked for irregularities. Paying attention to cracks, fractures and the occurrence of corrosion.

If there are joints, the product should be oiled with a care product before the functional test. We recommend a medical white oil based on paraffin oil.

Carry out a function check in accordance with the intended application of the product.

Essential tests for instruments with serration, jaw and ratchet are among others:

- Correct opening and closing (smooth movement, complete)
- Intact serration (bent, broken off)
- The jaws must close properly
- No throwing the jaws over (cross bite)

Defective products must not be used and must have undergone the complete preparatory treatment process again before being returned

## 11 Sterilization

Prior to sterilization, products must undergo cleaning and disinfection, be rinsed off without residue using demineralized water and subsequently dried. HEBUmedical recommends using a validated steam sterilization process (e.g. sterilizer in compliance with EN 285 and validated in accordance with DIN EN ISO 17665-1).

On using the fractionated vacuum method, **sterilization** must be performed with at least **134°C (USA 132° C)** with a **minimum dwell period of 3 minutes**. Vacuum drying must then be carried out for at least 20 minutes.

<b>V Parameters used for the validation of steam sterilization</b>	
<b>Prevacuum</b>	3 times
<b>Sterilization temperature</b>	132 °C
<b>Sterilization time</b>	1,5 minutes (half cycle method)
<b>Drying time</b>	20 minutes

The vapour must be free of ingredients, recommended limiting values of feed water and vapour condensate are determined through EN 285.

Other sterilization processes are compatible but not validated from HEBUmedical.

When loading, observe the recommended total weight. After the sterilization, check the sterile product packaging for damage, and inspect the sterilization indicators

The vessel sealing instruments BiSeal® must not be sterilized in hot air sterilizers.

## 11.1 Packaging

Compliant packaging of products for sterilization in line with ISO 11607. Packaging used must be suitable for the instruments and protect them from microbiological contamination during storage. The seal must not be under tension. HEBUmedical recommends container or hospital common sterilization paper/film packagings for sterilized packaging.

 During validation the instruments were packaged in hospital common sterilization packagings (paper/film packagings) and steam sterilized.

## 12 Lifetime

The steam sterilization procedure was validated by laboratory tests. The products were sterile validated at a pre-vacuum of at least 5min duration and a temperature of 134°C for a lifetime of 50 cycles. You can continue to use the instruments at your own responsibility over this cycle value if the tests described in chapter 7 have been successfully completed.

## 13 Storage

Storage of the sterilized instruments in a dry, clean and dust-free environment. Products must be protected from direct sunlight in storage.

## 14 Warranty/Repair

Our products are made of high-quality materials and are carefully checked before delivery. However, even when used as intended, they are subject to wear to a greater or lesser extent depending on the intensity of use.

This wear is due to technical reasons and is unavoidable.

However, if faults occur that are not related to wear, please contact our customer service department. Defective products may no longer be used and must have undergone the entire preparatory process before being returned.

## 15 Servicing and manufacturer address

Should you require the instructions for use in paper form, please use the contact details below. The instructions for use in paper form will be made available to you within seven calendar days of receipt of the request.

Alternatively, you can print out the electronic instructions for use yourself.

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### **HEBUmedical GmbH**