HEBU fx

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Fax +49 7461 94 71 - 22
eMail: service@HEBUmedical.de
Web: www.HEBUmedical.de

Contents

1	Symbol descriptions	22
2	Introduction	22
3	Scope	23
3.1	Intended Use	23
3.2	Contraindications.....	23
4	Warning.....	23
5	Product description	24
5.1	Rail clamping base	24
5.2	Universal clamp.....	24
5.3	Insulated rail clamping base.....	25
5.4	Columns	26
5.5	Articulated arms	27
5.6	Combined instrument holder	28
5.7	Endoscope holder	28
5.8	Sets	29
6	Installation and use	29
6.1	Assembly	29
6.2	Disassembly	30
6.3	Articulated arms	30
6.4	Combined instrument holder	31
6.5	Endoscope holder	31
7	Preparation.....	32
7.1	Reutilization restrictions	32
7.2	Information on instrument preparation	32
7.3	Preparation at the place of use	32
7.4	Ultrasound bath (optional)	33
7.5	Manual cleaning	33
7.6	Mechanical cleaning.....	33
7.7	Drying	35
8	Maintenance, inspection	35
8.1	Functional check	35
9	Sterilization.....	36
9.1	Packaging.....	36
10	Lifetime.....	37
11	Storage.....	37
12	Warranty / Repair	37
13	Manufacturer and servicing address	37

1 Symbol descriptions

Symbol	Definition
	CE-labelling
	Attention
	Validated Parameters
	Manufacturer
	Lot-description
	Reference code
	Medical device / FDA Prescription device
	Medical device
	Non sterile
	Keep away from sunlight
	Dry storage required
	(Electronic) instruction for use

2 Introduction

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them. By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following. In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

3 Scope

The HEBU-fx Articulated arm system is used in almost every area, mainly in abdominal surgery. The procedure may only be carried out by suitably trained and qualified specialist personnel.

3.1 Intended Use

A product for temporarily holding and positioning instruments, endoscopes and other products during a surgical procedure. It consists of a base, adjustable articulated arms and height-adjustable columns to adjust the instruments to the desired position. It is a reusable product.

3.2 Contraindications


No contraindications are known.

4 Warning

	Medical products are delivered in a non-sterile condition and must be cleaned, disinfected and sterilized prior to their initial use.
	The use of faulty instruments is in principle forbidden and they have to go through the whole cleaning process before return.
	Please take into consideration that through higher power a bigger damage of the tissue can result: f.e. on forceps: the power at the end of the jaw is higher than at the tip of the jaw
	Please observe the additional information enclosed with the products.
	Remove all protective sleeves and films prior to first using or preparation for use.
	The safe combination of different products or of products with implants must be reviewed prior to clinical application by the user.
	Avoid improper throwing or dropping of instruments
	Avoid mechanical overstressing of the instrument beyond the structural design, this can lead to breakage and deformation!
	A visual inspection of the instrument for damage and contamination must be carried out before each use!
	To prevent all contact corrosion, instruments with damaged surfaces must be separated immediately.
	If the products are used on patients with transmissible spongiform encephalopathy or HIV infection, we decline any responsibility for their reuse.
	After ophthalmical use, please pay attention to water quality during treatment (according to the specifications of AAMI TIR34 and the recommendations of the Rober-Koch-Institute on preparation of medical devices)!
	Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

5 Product description

5.1 Rail clamping base


Rail clamping base (HB 3328-30) made of aluminum with double clamping and with height adjustment for columns $\varnothing 16\text{mm}$. Not electrically insulate. Clamping range:

-Rail: 25x10 – 35x8mm

Can be steam-sterilized up to 134°C.

5.2 Universal clamp

Universal clamp (HB 3328-31) made of aluminum with height adjustment for columns $\varnothing 16\text{ mm}$. Not electrically insulated. Clamping range:

-Rail 25x10 – 35x8mm

-Round tube: $\varnothing 10 – 40\text{mm}$

-Table tops/flat section: 0 – 40mm

Can be steam-sterilized up to 134 °C.

5.3 Insulated rail clamping base

Insulated rail clamping base (HB 3328-32) for high-frequency and insulated uses, manufactured from robust, tested materials, front-handling, gentle and secure clamping of high-adjustable columns $\varnothing 16\text{mm}$. Tension resistance up to 1500V AC. Clamping range:

-Rail: 25x10 – 35x8 mm

Can be steam-sterilized up to 134°C


DANGER:

Verify regularly that the insulating plastic parts are free of cracks, porosity, and damages. Damaged plastic parts can reduce the penetration resistance. Damaged Insulated rail clamping bases must not be used again. For further information on the handling see product information for Insulated rail clamping bases.


5.4 Columns

The columns (see following table) serve to adjust the height of articulated arm and as fixed connections between the rail clamping base and the articulated arm.

Picture	HB-No.	Description	Dimensions [mm]	Connection threading top/bottom
	HB 3328-20	Column	Ø16x300	M12x1 / M8
	HB 3328-21	Column	Ø16x400	M12x1 / M8
	HB 3328-22	L-column	Ø16, 450x100	M12x1 / M8
	HB 3328-23	U-column	Ø16, 226x330x116	M12x1 / M8
	HB 3328-25	Stop lid	Ø20	M8
	HB 3328-26	End piece	Ø16	M8

5.5 Articulated arms

The HEBU articulated arms can be positioned exactly in any desired position. All three joints can be fixed and quickly released with the central handle. The purely mechanical central clamping unit allows easy and quick positioning of endoscopes and other instruments. All articulated arms listed below can be steam-sterilized up to 134°C.

Technical data

HB-No.	Description	Length L	Connection threading bottom	Connection threading top
HB 3328-10	Articulated arm	280mm	M8	M12x1
HB 3328-11	Articulated arm	400mm	M8	M12x1
HB 3328-12	Articulated arm	550 mm	M8	M12x1

5.6 Combined instrument holder

The combined instrument holder (HB 3328-01) made of stainless steel permits the locking of round parts of \varnothing 4-18 mm or cubical parts 0-18 mm. The combined instrument holder can be steam-sterilized up to 134 °C. Connection threading M8.

Compatibility with other products:

All standard instruments and wound-hook holders with the above listed dimensions can be held with the HEBU Combined Instrument holder.

5.7 Endoscope holder

The various optic holders (see table below) allow form-fittingly clamping of endoscopes or other optics. The clamping jaws and the wing-handle are made of plastic, the other parts of stainless steel. The optic holder can be steam-sterilized up to 134 °C.

HB-Nummer	Description	Clamping- \varnothing	Connection threading
HB 3328-02	Optic holder	\varnothing 5 und \varnothing 10	M8
HB 3328-03	Optic holder	\varnothing 4	M8

Amendment for sales in the United States:

US Federal Law restricts the endoscope holder system to sale by or on the order of a physician.

Indications for Use: HEBU Holding Systems consist of a table-mounted endoscope holder system intended for use by surgeons to hold endoscopes and arthroscopes with a diameter of 4mm to 10mm during general diagnostic and therapeutic procedures. The device is also intended for use by qualified surgeons for holding endoscopes during diagnostic and therapeutic neurologic procedures.

5.8 Sets

A complete set includes:

- a rail clamping base
- a column (I, L, or U-shape),
- one or two articulated arms (double-arm system)
- an optic or combined instrument holder for each articulated arm

Individual HEBU sets will be assembled on customer demand.

Before delivery to the clients, the column and the combined instrument holder or optic holder are normally rigidly threaded to the articulated arm and additionally fastened with pins.

The complete set can be steam-sterilized up to 134 °C.

6 Installation and use

For a “non-grounded” patient positioning (complying to the safety norms EN IEC 60601-1 and 2-2, VDE 0750 Part 2 and VDE 0753) an insulated rail clamping base must be used. Special care should be taken when these products are used together with high-frequency applications. Contacts between the product and the high-frequency device must be avoided.

6.1 Assembly

The assembly of the set can be done easily and without additional tools. All handling should be done with as little force as possible and as much force as necessary!

1. Attach the rail clamping base to the rail of the surgical table and secure it by turning the handle clockwise (Pos. 1, see figure section 3.1 and 3.2)
2. If necessary, remove the stop lid (HB 3328-25, see section 3.4) from the column by turning counterclockwise. Place the column with one hand into the column opening of the rail clamping base and secure the column with the other hand by turning the handle (Pos. 2) clockwise. Then screw the stop lid into the column clockwise and tighten.
3. The desired height and orientation of the column including the articulated arm can be adjusted at any time by loosening the handle (Pos. 2)

Important:

Hold the column with one hand and loosen the handle with the other hand. When the desired position is reached, tighten the handle.

DANGER:

If the rail clamping base is not correctly tightened, it may become loose and may cause injuries to the patient.


6.2 Disassembly

The disassembly of the set can be performed also without additional tools and occurs in the reverse order of the assembly:

1. If necessary, remove the stop lid (HB 3328-25) from the column by turning counterclockwise.
Important:
Hold the articulated arm with one hand and loosen the handle (Pos. 2) with the other hand.
2. Pull the articulated arm together with the column out of the opening of the rail clamping base. Thread the stop lid (HB 3328-25) lightly into the columns to prevent loss.
3. Remove the rail clamping base from the rail of the surgical table after loosening the handle (Pos. 1).

6.3 Articulated arms

Holding-Releasing-Positioning-Fixing-Controlling:

1. Hold the articulated arm with one hand in the anterior segment and manipulate the central clamping device with the other hand
2. To release turn the central handle (Z) counterclockwise as far as necessary
3. Move the articulated arm into the desired position
4. To secure, turn the central handle (Z) clockwise
5. Check that the articulated arm is functioning properly

DANGER:

- If the articulated arm is not tightened properly it can become loose and move, potentially causing injuries.
- Do not put too much weight on the articulated arm.
- The articulated arm can transmit electrical current and heat to the patient. Avoid the contact between the articulated arm and any sources of electrical current or heat.

WARNING:

- The fixation of the articulated arm is based on the principle of friction. Changing the position without losing the clamping mechanism can cause damage and will shorten the lifespan of the articulated arm. The articulated arm can be adjusted with little force. If the central clamping handle (Z) has been fully loosened it has to be turned in the clockwise direction!
- Store the articulated arms so that they are not damaged.

6.4 Combined instrument holder

By turning the handle counterclockwise the instrument holder is opened and the instrument or wound hook can be placed between the clamping jaw. By turning the handle clockwise the device is secured.

6.5 Endoscope holder

By turning the handle counterclockwise the endoscope holder is opened and the endoscope can be inserted from the front. By turning the handle clockwise the endoscope is secured.

WARNING:

- The endoscope can be damaged if the handle is tightened too strongly.
- Do not try to change the position of the endoscope by turning the camera. This may damage the endoscope. Always open the central handle of the articulated arm for repositioning.

Amendment for sales in the United States:

US Federal Law restricts the endoscope holder system to sale by or on the order of a physician.

Indications for Use: HEBU Holding System consists of a table-mounted endoscope holder system intended for use by surgeons to hold endoscopes and arthroscopes with a diameter of 4mm to 10mm during general diagnostic and therapeutic procedures. The device is also intended for use by qualified surgeons for holding endoscopes during diagnostic and therapeutic neurologic procedures.

7 Preparation

The person in charge of preparatory treatment is responsible for ensuring that the treatment is duly carried out using the relevant equipment, materials and personnel in the treatment facility and so achieves the desired result. This necessitates validation and routine monitoring of the process used. We urge you to take note of the national regulations dealing with instrument preparation.

The validated parameters refer to reusable surgical instruments. The validated parameters should also be observed for the other products described, unless a different procedure is explicitly described.

7.1 Reutilization restrictions

Frequent repeat preparatory treatment has minimal effects on the product. The end of the product life is normally determined by wear and damage due to use

**WARNING:**

The articulated arm must not be immersed in liquids.

7.2 Information on instrument preparation

- Use cleaning and/or disinfection agents with a pH-value within 9-10. Please observe manufacturer instructions regarding dosage, exposure time and renewal of solutions.
- Do **not** use hard brushes (e.g. metal brushes or metal sponges) or coarse abrasive cleaners.
- Never leave instruments in cleaning or disinfection agents for longer than the specified time.
- Only used demineralized water for rinsing.
- Rinse and dry carefully through channels and pipes.
- Sensitive instruments must be cleaned in a storage or clamping fixture.
- Observe manufacturer instructions of cleaning – and sterilizing equipment.
- Products made of aluminum alloys will be damaged by basic (pH>9) cleaning agents and solvents.

**WARNING:**

The central handle of the articulated arm must be tightened during pretreatment and cleaning/disinfection!

7.3 Preparation at the place of use

Directly after using remove coarse dirt of the instruments and rinse out the working cannulas. Do not use fixing agents or hot water (> 40°C), as this results in residues becoming fixed and can affect the success of the subsequent cleaning operation

Dismantle and/or open instruments as far as possible. Within short time after use the instruments clean the instruments for reducing a drying of the residues.

This enables an easier cleaning. If instruments come into contact with corroding medicines or cleaning agents, wash these up with water immediately after use.

Longer drying times, e.g. for dry disposal are not validated and therefore not recommended.

The drying time during validation was 1 hour.

7.4 Ultrasound bath (optional)

All instruments must be opened, dismantled and any cavities rinsed through.

Place instruments in the screen basket in such a way that overlaps and contact between instruments are avoided. Add cleaning agent to the water and adjust the temperature of the solution in line with the cleaning agent manufacturer's instructions.

The cleaning in the ultrasound bath should be at 35-40 kHz, 5 minutes at least.

To validate cleaning in an ultrasonic bath, the test items were ultrasonically treated in Neodisher Mediclean forte 0,5 % for 5 minutes.

Subsequently rinse instruments including all cavities before cleaning and disinfection.

Medicine products which are possessing a bad ultrasound transmission, e.g. soft materials are not usable for the ultrasound bath.

7.5 Manual cleaning

Since mechanical processes can be standardized, reproduced and therefore validated, mechanical cleaning/disinfection should be preferred to manual processes.

Manual cleaning and disinfection process is not validated and therefore needs to be validated additionally by the end user.

7.6 Mechanical cleaning

Based on international standards (EN ISO 15883) and national directives, only validated machine cleaning and disinfection methods may be used. For the mechanical cleaning we recommend a standard programme for surgical instruments, f.e. instruments from Miele.

Only completely demineralized water should be used for cleaning, neutralisation and rinsing, in accordance with the „Guidance Compiled by the DGKH (Germany Society for Hospital Hygiene), DGSV (German Society for Sterile Supply) and AKI (Working Group on Instrument Reprocessing) for the Validation and Routine Monitoring of Automated Cleaning and Thermal Disinfection Processes for Medical Devices as well as Advice on Selecting Washer-Disinfectors“ (which refers to DIN EN ISO 15883-1 Point 6.4.2)

Flexible (complex) instruments with invisible surfaces must be pre-cleaned manually before mechanical cleaning.



We recommend for all push shafts, -and pipe shafts instruments and instruments whose surfaces are on top of each others during the cleaning (f.e. bone forceps and gouge forceps) a manually pre-cleaning for an optimal cleaning result without residues.

Observe the following by loading:

- Place the dismantled/opened instruments securely in the tray.
- Instruments with openings and gaps have to be faced down with the opened side so that they can be cleaned and no water of the cleaning process is collecting inside them. If available use balanced devices for rinsing

- Place the instruments with joints in an opened position into the cleaning, -and disinfection machine.
- Do not overload trays, avoid creating any overlaps.

Preliminary rinsing (cold, if applicable fully demineralized water without additives) is followed by chemical. The chemical cleaning should take place at **40°C -60°C** for at **least 5 minutes**.

We recommend products with a **pH-value within 9-10**, e.g. Neodisher MediClean forte from Dr. Weigert. The cleaning agents used should be selected depending on the material and properties of the instruments and in accordance to national regulations: If there is a high chloride concentration in the water, pitting and tension crack corrosion can occur on the instruments. The occurrence of this type of corrosion is minimized by using alkaline cleaning agents and demineralized water. By adding an acid-based neutralization medium, the rinsing off of alkaline cleaning agent residues is facilitated during the first intermediate rinsing process (warm or cold water).

In order to prevent the formation of deposits, it is advisable to use neutral cleaners where the water quality is unfavourable. After the second intermediate rinsing process, thermal disinfection takes place.

The thermal disinfection should take place at temperatures of between **80 and 95°C**, with an **exposure time as outlined in EN ISO 15883**.

After the finishing of the programme take the good out of the machine because corrosion can arise if the instrument remains in the machine.

V Parameters used for the validation of preparation	
Pre-rinsing	1 minute with cold tap water
Cleaning	Temperature: 55 °C
	Soaking Time: 5 minutes (worst case)
	Neodischer Mediclean forte 0,4% (worst case)
Neutralization	Temperature: cold DI water
	Soaking Time: 2 minutes
	Neodisher Z 0,1%
Post-rinsing	2 minutes with cold DI water
Disinfection	Temperature: 90 °C (A ₀ 3000)
	Soaking Time: 5 minutes

7.7 Drying

Ensure adequate drying by the cleaning and disinfection device or using other suitable measures.

V Drying was omitted in the validation (worst case condition).

8 Maintenance, inspection

After cooling to room temperature, the instruments must be visually inspected for protein residues and other contamination, paying particular attention to cavities, blocks, inclusions, pipes, and other inaccessible areas. Instruments which are not free of residues must be returned for a complete retreatment process.

Carry out the functional check mentioned above. Instruments with stains, which are blunt, bent, no longer function or which are otherwise damaged must be segregated!

To help identify faulty instruments that need to be sorted out, we recommend the brochure "Instrument Reprocessing" from the Working Group "Instrumenten Aufbereitung". This includes Chapter 8 "Checks and Care" and Chapter 12 "Surface Changes: Deposits, Discoloration, Corrosion, Aging, Swelling and Stress Cracks".

8.1 Functional check

A newly purchased product must be subjected to a thorough visual and function check after its delivery and before each use.

Products must be checked for irregularities. Paying attention to cracks, fractures and the occurrence of corrosion.

If there are joints, the instruments should be oiled with a care product before the functional test. We recommend a medical white oil based on paraffin oil.

Check instruments with joints for ease of movement. Carry out a function check in accordance with the intended application of the instrument.

Defective products must not be used and must have undergone the complete preparatory treatment process again before being returned.

9 Sterilization

Prior to sterilization, products must undergo cleaning and disinfection, be rinsed off without residue using demineralized water and subsequently dried. HEBUmedical recommends using a validated steam sterilization process (e.g. sterilizer in compliance with EN 285 and validated in accordance with DIN EN ISO 17665-1).

NOTE:

The central tension handle of the articulated arm must be open during sterilization!

WARNING:

By releasing the central handle, the articulated arm may suddenly slide downwards and can cause injury or damage.

To avoid this put the articulated arm down, hold it with one hand and with the other hand release the central handle.

The validated parameters refer to reusable surgical instruments. The validated parameters should also be observed for the other products described, unless a different procedure is explicitly described.

On using the fractionated vacuum method, **sterilization** must be performed with at least **134°C (USA 132° C)** with a **minimum dwell period of 3 minutes**. Vacuum drying must then be carried out for at least 20 minutes.

 Parameters used for the validation of steam sterilization	
Prevacuum	3 times
Sterilization temperature	132 °C
Sterilization time	1,5 minutes (half cycle method)
Drying time	20 minutes

The vapour must be free of ingredients, recommended limiting values of feed water and vapour condensate are determined through EN 285.

Other sterilization processes are compatible but not validated from HEBUmedical.

When loading, observe the recommended total weight. After the sterilization, check the sterile product packaging for damage, and inspect the sterilization indicators

9.1 Packaging

Compliant packaging of products for sterilization in line with ISO 11607. Packaging used must be suitable for the instruments and protect them from microbiological contamination during storage. The seal must not be under tension. HEBUmedical recommends container or hospital common sterilization paper/film packagings for sterilized packaging.

 During validation the instruments were packaged in hospital common sterilization packagings (paper/film packagings) and steam sterilized.

10 Lifetime

The steam sterilization procedure was validated by laboratory tests. The products were sterile validated at a pre-vacuum of at least 5min duration and a temperature of 134°C for a lifetime of 50 cycles.

You can continue to use the instruments at your own responsibility over this cycle value if the tests described in chapter 7 have been successfully completed.

11 Storage

Store products in a dry, clean and dust-free environment at moderate temperatures from 5°C to 40°C. Protect from the effects of the sun's rays and artificial light.

12 Warranty / Repair

Our products are manufactured from high-grade materials and carefully checked prior to dispatch. However, even if used properly in accordance with their intended purpose they are subject to a greater or lesser degree of wear depending on their intensity of use.

This wear is technically induced and unavoidable.

Should faults occur independently of wear, please contact our customer services. Defective products should no longer be used.

They must undergo the complete preparatory treatment process before being returned.

13 Manufacturer and servicing address

Should you require the instructions for use in paper form, please use the contact details below. The instructions for use in paper form will be made available to you within seven calendar days of receipt of the request.

Alternatively, you can print out the electronic instructions for use yourself.

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen / Germany
Tel. +49 7461 94 71 - 0
Fax +49 7461 94 71 - 22
eMail: service@HEBUmedical.de
Web: www.HEBUmedical.de

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**Our vast range of products covers over 10 000 various instruments.
They can be found in Internet or request our catalog.**

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