Einfache Schneidinstrumente /

Simple cutting instruments / Instruments de coupe simples / Instrumentos de corte simples / Strumenti di taglio semplici / Прости инструменти за рязане / Enkle skæreinstrumenter / Lihtsad lõikamisinstrumendid / Yksinkertaiset leikkausvälineet / Απλά όργανα κοπής / Jednostavni instrumenti za rezanje / Vienkārši griešanas instrumenti / Paprasti pļovimo instrumentai / Eenvoudige snijinstrumenten / Proste instrumenty do cięcia / Instrumentos de corte simples / Instrumente simple de tăiere / Enkla skärinstrument / Jednoduché strihové nástroje / Preprosti rezalni instrumenti / Jednoduché řezací nástroje / Egyszerű vágóeszközök

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Seite 5

Operating Manual

Page 23

Mode d'emploi

Page 41

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Página 59

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Página 77

DEUTSCH

ENGLISH

FRANÇAIS

ESPAÑOL

ITALIANO

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Table of contents

1	Symbol descriptions	24
2	Introduction	25
3	Intended use	25
	3.1 Intended purpose	25
	3.2 Indication	28
	3.3 Contraindications	29
	3.4 Patient target group	29
4	Warnings	30
5	Handling	31
6	Preparation	31
	6.1 Information on instrument preparation	31
	6.2 Preparation at the place of use	32
	6.3 Ultrasound bath (optional)	32
	6.4 Manual cleaning	32
	6.5 Mechanical cleaning	32
	6.6 Drying	34
7	Maintenance, inspection	35
	7.1 Functional check	35
8	Sterilization	36
	8.1 Packaging	36
9	Lifetime	36
10	Storage	37
11	Warranty / Repair	37
12	Disposal	37
13	Service and manufacturer address	38

1 Symbol descriptions

Symbols	Definition
	CE-labelling
	Attention
	Validated Parameters
	Manufacturer
	Lot-description
	Reference code
	Medical device / FDA Prescription device
	Medical device
	Non sterile
	Keep away from sunlight
	Dry storage required
	(Electronic) instruction for use

2 Introduction

By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following.

In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

3 Intended use

3.1 Intended purpose

Dermal Punches	A hand-operated surgical instrument used for removing blemishes or for sampling skin tissue. It usually consists of a handle with tips at the end of work in various configurations; some peaks include a closed cutting zone, edged by a cutting blade. It is a reusable instrument intended for transient use.
Dermatome	A surgical hand instrument for cutting thin skin slices for transplantation or excision of small skin lesions. This requires a special cutting blade, which is used for this purpose in the instrument. It is a reusable instrument intended for transient use.
Hand/ Foot instruments	An instrument for podiatrists to treat deformities and diseases of the foot. It is a non-invasive, reusable instrument.
Surgical scalpels; Knife; phalange knife; resection knife; dissection knife; microvascular knife; stick instruments, knives; Mucosal knife	A surgical instrument with a handle and a scalpel blade (not interchangeable) used by doctors to cut or prepare tissues with downward pressure. The instrument is made of stainless steel. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.
Scalpel Handles	An interchangeable component of a scalpel that works as a handle and picks up a matching blade. The instrument is made of stainless steel. It is a non-invasive, reusable instrument.
Amputating Knives	A heavy-duty surgical instrument with a handle and a single-sided or double-sided cutting blade of various shapes and sizes for use in limb amputations. It is a reusable instrument intended for transient use.
Enucleators	A surgical hand instrument for the dissection of the tonsils during an ENT procedure. It is usually spoon-like or rounded at the working ends. It is manufactured in various designs, e.g. with a shaft, which merges at the distal end in a rounded, bent or angled tip, or even double-ended. The working end may be tapered, flat, sharp or dull, straight or angled, serrated or smooth. It is made of stainless steel and is available in different sizes and shapes. It is a reusable instrument intended for transient use.

Tissue dissectors; Sepa-rator	A surgical hand instrument, usually in spoon shape or rounded at the end of work. It is used to separate soft tissue or body structures in interventions of general or plastic surgery. It usually has a handle that continues into a shaft that has a tip at the distal end. The tip may be pointed or flat, sharp or dull, angled or straight on the shaft. The instrument is made of high quality stainless steel and available in various shapes and dimensions. It is a reusable instrument intended for transient use. Es ist nicht zur Verwendung in direktem Kontakt mit dem Herz, dem zentralen Kreislaufsystem oder dem zentralen Nervensystem vorgesehen.
Micro Knives	A surgical instrument for cutting tissue that consists of a long slender metal handle that tapers towards the distal end. The internal radius of the hook has a cutting edge. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.
Plaster-of-Paris knives	A hand instrument with a short, flat, strong cutting blade with a sharp edge for cutting or shortening plaster casts. It is usually a one-piece product made of stainless steel. It is a non-invasive, reusable instrument.
Chisels and osteotomes; rhachiotome; flat chisel; Nasal osteotomes	A surgical, chisel-like hand instrument for cutting and / or shaping bones by bumping during orthopedic surgery. It is held by the surgeon who applies a manual force to the proximal end of the instrument via a surgical hammer. The distal end (the cutting or sharp edge) is sharp, often flat, sometimes curved (concave), and usually folded on both sides. It is made of stainless steel. It is a reusable instrument intended for transient use.
Cartilage knife	A surgical instrument for cutting, scraping or shaping cartilage tissue. It is a reusable instrument intended for transient use
Rhinoplastic and Nasal Knives; Tonsil Knives, Dissectors, Enucleators, Retractors; Tympanum Knives; Trigeminal Knives; Ear knives, dissectors and raspatory; Periosteal knife	A surgical hand instrument for cutting anatomical conditions during a general ENT procedure. It is made of stainless steel and has a handle in various designs, which merges into a narrow working end with a narrow, sharp blade. It is a reusable instrument intended for transient use.
Tendon stripper	A surgical instrument consisting of a handle and a semi-circular, domed working end that ends in a sharp cutting edge. Tendon strippers are used to cut out a piece of ligament, tendon or fascia to transplant them. It is a reusable instrument intended for transient use.
Elevators; Meniscus retractors	A hand-held surgical instrument used to lift, position, or lever open bone structures, other anatomical structures or surgical material during an orthopedic surgery. It is made of stainless steel and is available in a variety of sizes, styles and working ends. It is a reusable instrument intended for transient use.
Meniscotomy	A cutting, surgical instrument whose handle is T-shaped to be able to apply more force. The handle merges into a chisel-shaped edge. The instrument is used to cut crescent-shaped parts of the capsule or cartilage area. It is a reusable instrument intended for transient use.

Sternum chisel	A single-sided tapered blade surgical instrument for cutting and shaping bone during orthopedic surgery. The surgeon holds it in one hand while hitting the proximal end with a surgical gavel in the other hand. The distal end (the cutting edge) is sharp, flat or curved (concave), also called a hollow blade). It is made of high quality stainless steel and may consist of one piece, or have a synthetic handle. It is a reusable instrument intended for transient use.
Circumcision instruments	A surgical instrument for the controlled removal of the penis foreskin during circumcision. It is typically made of high quality stainless steel and may have various shapes, e.g. the bell clamp (bell-shaped clamp). The foreskin is pushed over the bell-shaped protective body and the glans are positioned therein. A screw mechanism is tightened so that the foreskin is compressed in an annular manner and fitted with e.g. a scalpel along the compression ring can be cut off. It is a reusable instrument intended for transient use.
Teflon cutting block	A product that serves as a base for cutting and shaping various materials. It is a non-invasive reusable product.
Septum Gouges	A surgical instrument with a single sided beveled blade for cutting and shaping nasal bones during ear, nose and throat (ENT) operations. The surgeon holds it in one hand while beating the proximal end of the instrument with a surgical gavel in the other hand. The distal end (the cutting edge) is sharp, but often flat but may be curved (concave). It is typically one piece and made of high quality stainless steel. It is a reusable instrument intended for transient use.
Adenoid currettes; Swivel Knives	A surgical hand instrument for cutting and excision of lymphoid tissue during an adenoidectomy. It is usually performed as a long, narrow instrument with a handle and a concave, spoon-like working part. The working part has a sharp edge for removing the tissue without injuring the surrounding muscles. It is usually made of stainless steel. It is a reusable instrument intended for transient use.
Adenotomes	A surgical instrument designed to excise hypertrophic lymphoid tissue in the nasopharynx (i.e., pharyngeal tonsils or polyps) during an adenoidectomy. It is usually a long, slender, manual hand instrument with cutting edges at the working end and a handle with a mechanism for actuating the blades. It is available in various designs and sizes and can be a one-piece instrument or require the insertion of blades into the handle. It is a reusable instrument intended for transient use.
Cataract Knives; Lancet Knife	An ophthalmic surgical instrument for placing precise cuts in ocular tissue during ophthalmic surgery on the eye and surrounding structures. It is usually designed as a one-piece instrument with a sharp, single-edged blade at the distal end and a handle at the proximal end. It can have different blade shapes, e.g. phaco (dull or sharp), crescent-shaped, straight, and micro-phaco (for cataract operations). It is usually made of stainless steel and some models may be equipped with diamond blades. It is a reusable instrument intended for transient use.
Foreign Body instruments	An ophthalmic instrument to remove foreign bodies from inside the eye or from the surface of the eyeball. It is typically made of high quality stainless steel and has a fine handle at the proximal end and a thin concave blade with a rounded nose at the distal end. The ophthalmologist uses it to remove the foreign body / object. It is a reusable, non-invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.

Tongue Holding Forceps	An instrument for grasping, holding or manipulating the tongue during a examinations or procedures. It usually has a self-retaining, scissor-like design with ring handles. The working end is different, e.g. straight, angled or curved with large, oval-fenced, cross-scored mouth surfaces for a better grip. Some models may have replaceable, grooved rubber inserts in the jaw part. The instrument is made of stainless steel. It is a reusable, invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.
Measuring probes	A product / instrument used for comparative measurement, e.g. to measure inner and outer diameters, lengths, depths or thicknesses. The instrument ist not calibrated. It is a reusable, invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.

3.2 Indication

Simple cutting instruments are used in the following areas:

Plaster-of-Paris knives	The instrument is used for treatments. It may only be used by trained and qualified professionals.
Foreign Body instruments	The instrument is used for treatments in ophthalmology. It may only be used by trained and qualified professionals.
Hand/ Foot instruments	The instrument is used for nail and foot care.
Cataract Knives; Lancet Knife	The instrument is used for surgical procedures in ophthalmology. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Circumcision instruments	The instrument is used for urological surgical procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Amputating Knives	The instrument is used in operations for the amputation of a limb. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Tendon stripper	The instrument is used in various surgical procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Tissue dissectors; Separator	The instrument is used for various procedures in general or plastic surgery. It is not intended for use in direct contact with the heart, central circulatory system or central nervous system. The procedure must be performed by trained and qualified personnel.
Dermal Punches	The instrument is used in plastic surgery. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Meniscotomy	The instrument is used during a menisectomy. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Tongue Holding Forceps; Measuring probes	The instrument is used for ENT examinations and treatments. It may only be used by trained and qualified professionals.
Scalpel Handles; Teflon cutting block	The instrument is used in various procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Simple cutting instruments

Dermatome; Surgical scalpels; The instrument is used in various surgical procedures. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Knives; Mucosal knife; Micro Knives

Chisels and osteotomes; rhammatome; flat chisel; Nasal osteotomes; Cartilage knife; The instrument is used in various surgical procedures in orthopedics. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Elevators; Meniscus retractors; Sternum chisel

Enucleators; Rhinoplastic and Nasal Knives; Tonsil Knives, Dissectors, Enucleators, Retractors; Tympanum Knives; The instrument is used during surgical ENT procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Trigeminal Knives; Ear knives, dissectors and raspatory; Periosteal knife; Septum Gouges; Adenoid currettes; Swivel Knives; Adenotomes

3.3 Contraindications

No contraindications are known.

3.4 Patient target group

The products are suitable for all patient groups.

4 Warnings

	Medical products are delivered in a non-sterile condition and must be cleaned, disinfected and sterilized prior to their initial use.
	The use of faulty instruments is in principle forbidden and they have to go through the whole cleaning process before return.
	Please take into consideration that through higher power a bigger damage of the tissue can result: f.e. on forceps: the power at the end of the jaw is higher than at the tip of the jaw
	Remove all protective sleeves and films prior to first using or preparation for use.
	The safe combination of different products or of products with implants must be reviewed prior to clinical application by the user.
	Avoid improper throwing or dropping of instruments.
	Avoid mechanical overstressing of the instrument beyond the structural design, this can lead to breakage and deformation!
	A visual inspection of the instrument for damage and contamination must be carried out before each use!
	To prevent all contact corrosion, instruments with damaged surfaces must be separated immediately.
	If the products are used on patients with transmissible spongiform encephalopathy or HIV infection, we decline any responsibility for their reuse.
	After ophthalmical use, please pay attention to water quality during treatment (according to the specifications of AAMI TIR34 and the recommendations of the Robert-Koch-Institute on preparation of medical devices)!
	Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

5 Handling

The type of treatment must be determined in each individual case by the surgeon in cooperation with the internist and the anaesthetist.

For operational use in various surgical disciplines must be done by appropriately trained and qualified personnel.

6 Preparation

The person in charge of preparatory treatment is responsible for ensuring that the treatment is duly carried out using the relevant equipment, materials and personnel in the treatment facility and so achieves the desired result. This necessitates validation and routine monitoring of the process used. We urge you to take note of the national regulations dealing with instrument preparation.

The validated parameters refer to reusable surgical instruments. The validated parameters should also be observed for the other products described, unless a different procedure is explicitly described.

6.1 Information on instrument preparation

- Use cleaning and/or disinfection agents with a pH-value within 9-10.
- Please observe manufacturer instructions regarding dosage, exposure time and renewal of solutions.
- Use suitable brushes that do not remove material, i.e. **no** hard brushes (such as metal brushes and metal sponges) or coarse abrasive cleaners.
- Never leave instruments in cleaning or disinfection agents for longer than the specified time.
- Rinse and dry carefully through channels and pipes.
- Sensitive instruments must be cleaned in a storage or clamping fixture.
- Observe manufacturer instructions of cleaning – and sterilizing equipment.

6.1.1 Water quality

For cleaning, neutralization and rinsing, we recommend the use of demineralized water in accordance with the "Guideline DGKH, DGSV, AKI for the validation and routine monitoring of automated cleaning and thermal disinfection processes for medical devices and on the principles of device selection" (the guideline refers to DIN EN ISO 15883-1 point 6.4.2).

According to the guideline, the following values are recommended:

- Conductivity: $\leq 15 \mu\text{s}/\text{cm}$ (deviating from the table in DIN EN 285)
- pH value: 5 - 7
- Total hardness: $\leq 0.02 \text{ mmol CaO}/\text{l}$
- Salt content: $\leq 10 \text{ mg}/\text{l}$
- Phosphate (as P_2O_5): $\leq 0.5 \text{ mg}/\text{l}$
- Silicate (as SiO_2): $\leq 1 \text{ mg}/\text{l}$

■ Chloride: ≤ 2 mg/l

6.2 Preparation at the place of use

Immediately after use, remove coarse dirt from the instruments with a suitable brush (see chapter 6.1) and rinse out the working cannulas. Do not use fixing agents or hot water ($> 40^{\circ}\text{C}$), as this results in residues becoming fixed and can affect the success of the subsequent cleaning operation

Dismantle and/or open instruments as far as possible. Within short time after use the instruments clean the instruments for reducing a drying of the residues.

This enables an easier cleaning. If instruments come into contact with corroding medicines or cleaning agents, wash these up with water immediately after use.

Longer drying times, e.g. for dry disposal are not validated and therefore not recommended.

V The drying time during validation was 1 hour.

6.3 Ultrasound bath (optional)

All instruments must be opened, dismantled and any cavities rinsed through.

Place instruments in the screen basket in such a way that overlaps and contact between instruments are avoided. Add cleaning agent to the water and adjust the temperature of the solution in line with the cleaning agent manufacturer's instructions.

The cleaning in the ultrasound bath should be at **35-40 kHz, 5 minutes** at least.

V To validate cleaning in an ultrasonic bath, the test items were ultrasonically treated in Neodisher Mediclean forte 0,5 % for 5 minutes.

Subsequently rinse instruments including all cavities before cleaning and disinfection.

The ultrasonic bath should not be used for medical devices where material could be removed in the ultrasonic bath, e.g. products made of soft plastic or coated instruments.

6.4 Manual cleaning


Since mechanical processes can be standardized, reproduced and therefore validated, mechanical cleaning/disinfection should be preferred to manual processes. Manual cleaning and disinfection process is not validated and therefore needs to be validated additionally by the end user.

6.5 Mechanical cleaning



Due to international standards (EN ISO 15883) and national guidelines, only validated mechanical cleaning and disinfection procedures should be used. For automated cleaning, we recommend a standard program for surgical instruments, e.g. instruments from Miele.

Observe the following by loading:

- Place the disassembled or opened instruments securely in the tray.
- Instruments with openings and recesses must be placed with the open side facing downwards so that they can be cleaned and no water from the cleaning process can collect in them.
- If available, use a coordinated rinsing device.

- Do not overload trays, avoid creating any overlaps

Preliminary rinsing is followed by chemical.

The chemical cleaning should take place at **40°C -60°C for at least 5 minutes.**

We recommend products with a **pH-value within 9-10**, e.g. Neodisher MediClean forte from Dr. Weigert. The cleaning agents used should be selected depending on the material and properties of the instruments and in accordance to national regulations: If there is a high chloride concentration in the water, pitting and tension crack corrosion can occur on the instruments.

The occurrence of this type of corrosion is minimized by using alkaline cleaning agents and demineralized water. By adding an acid-based neutralization medium, the rinsing off of alkaline cleaning agent residues is facilitated during the first intermediate rinsing process.

In order to prevent the formation of deposits, it is advisable to use neutral cleaners where the water quality is unfavourable.

After the second intermediate rinsing process, thermal disinfection takes place.

The thermal disinfection should be carried out with demineralized water that corresponds to microbiological drinking water quality (< 100 CFU/ml drinking water) at **80 to 95°C** and an expo-

sure time in accordance with **EN ISO 15883**.

After the finishing of the programme take the good out of the machine because corrosion can arise if the instrument remains in the machine.

V Parameters used for the validation of preparation	
Pre-rinsing	1 minute with cold tap water
Cleaning	Temperature: 55 °C
	Soaking Time: 5 minutes (worst case)
	Neodischer Mediclean forte 0,4% (worst case)
Neutralization	Temperature: cold DI water
	Soaking Time: 2 minutes
	Neodisher Z 0,1%
Post-rinsing	2 minutes with cold DI water
Disinfection	Temperature: 90 °C (A ₀ 3000)
	Soaking Time: 5 minutes

6.6 Drying

Ensure adequate drying by the cleaning and disinfection device or using other suitable measures.

V Drying was omitted in the validation (worst case condition).

7 Maintenance, inspection

After cooling to room temperature, the instruments must be visually inspected for protein residues and other contamination. Slits, barriers, locks, tubes and other areas that are difficult to access must be thoroughly inspected. Instruments that are not residue-free must be repeatedly subjected to the entire reprocessing process.

To ensure that surgical instruments can be used for their intended purpose after reprocessing, it is necessary to carry out a functional test after cleaning, disinfection and drying of the visual inspection and care measures. Carry out the functional tests described in chapter 7.1.

Instruments that are stained, blunt, bent, no longer functional, damaged in any other way or the instrument labeling is no longer recognizable (e.g. labeling, UDI code no longer legible) must be discarded!

To help identify faulty instruments that need to be sorted out, we recommend the brochure "Instrument Reprocessing" from the Working Group "Instrumenten Aufbereitung". This includes Chapter 8 "Checks and Care" and Chapter 12 "Surface Changes: Deposits, Discoloration, Corrosion, Aging, Swelling and Stress Cracks".

7.1 Functional check

A newly purchased product must be subjected to a thorough visual and function check after its delivery and before each use.

Products must be checked for irregularities. Paying attention to cracks, fractures and the occurrence of corrosion.

If there are joints, the instruments should be oiled with a care product before the functional test. We recommend a medical white oil based on paraffin oil.

Check instruments with joints for ease of movement.

The labeling must also be checked for legibility.

Carry out a function check in accordance with the intended application of the instrument.

Essential tests for instruments with cutting edge are among others:

- Smooth cutting edges

Defective products must not be used and must have undergone the complete preparatory treatment process again before being returned.

8 Sterilization

Prior to sterilization, products must undergo cleaning and disinfection, be rinsed off without residue using demineralized water and subsequently dried. HEBUmedical recommends using a validated steam sterilization process (e.g. sterilizer in compliance with EN 285 and validated in accordance with DIN EN ISO 17665-1).

The validated parameters refer to reusable surgical instruments. The validated parameters should also be observed for the other products described, unless a different procedure is explicitly described.

On using the fractionated vacuum method, **sterilization** must be performed with at least **134°C (USA 132°C)** with a **minimum dwell period of 3 minutes**. Vacuum drying must then be carried out for at least 20 minutes.

V Parameters used for the validation of steam sterilization	
Prevacuum	3 times
Sterilization temperature	132 °C
Sterilization time	1,5 minutes (half cycle method)
Drying time	Drying time

The vapour must be free of ingredients, recommended limiting values of feed water and vapour condensate are determined through EN 285.

Other sterilization processes are compatible but not validated from HEBUmedical.

When loading, observe the recommended total weight. After the sterilization, check the sterile product packaging for damage, and inspect the sterilization indicators.

8.1 Packaging

Compliant packaging of products for sterilization in line with ISO 11607. Packaging used must be suitable for the instruments and protect them from microbiological contamination during storage. The seal must not be under tension. HEBUmedical recommends container or hospital common sterilization paper/film packagings for sterilized packaging.

V During validation the instruments were packaged in hospital common sterilization packagings (paper/film packagings) and steam sterilized.

9 Lifetime

The steam sterilization procedure was validated by laboratory tests. The products were sterile validated at a pre-vacuum of at least 5min duration and a temperature of 134°C for a lifetime of 50 cycles.

You can continue to use the instruments at your own responsibility over this cycle value if the tests described in chapter 7 have been successfully completed.

10 Storage

Store products in a dry, clean and dust-free environment at moderate temperatures from 5°C to 40°C.

Protect from the effects of the sun's rays and artificial light.

11 Warranty / Repair

Our products are manufactured from high-grade materials and carefully checked prior to dispatch. However, even if used properly in accordance with their intended purpose they are subject to a greater or lesser degree of wear depending on their intensity of use.

This wear is technically induced and unavoidable.

Should faults occur independently of wear, please contact our customer services. Defective products should no longer be used.

They must undergo the complete preparatory treatment process before being returned.

12 Disposal

Before disposal, the product must undergo the complete reprocessing and sterilization process in accordance with the procedure described in these instructions for use.

Risk of infection from contaminated products!

- When disposing of the product, its components and their packaging, the national regulations must be observed.

Risk of injury from sharp-edged and/or pointed products!

- When disposing of the product, ensure that the packaging prevents injury from the product.

13 Service and manufacturer address

Should you require the instructions for use in paper form, please use the contact details below. The instructions for use in paper form will be made available to you within seven calendar days of receipt of the request.

Alternatively, you can print out the electronic instructions for use yourself.

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen / Germany
Tel. +49 7461 94 71 - 0
Fax +49 7461 94 71 - 22
eMail: service@HEBUmedical.de
Website: www.HEBUmedical.de

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