Scheren /

Scissors / Ciseaux / Tijeras / Forbici / Ножици / Saks / Käärid / Sakset / Ψαλίδι / Škare / Šķēres / Žirklēš / Schaar / Nożyce / Tesoura / Foarfece / Sax / Nožnice / Škarj / Nůžky / Olló

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Seite 5

Operating Manual

Page 21

Mode d'emploi

Page 37

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Página 53

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Página 71

DEUTSCH

ENGLISH

FRANÇAIS

ESPAÑOL

ITALIANO

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Table of contents

1	Symbol descriptions	22
2	Introduction	23
3	Intended use	23
	3.1 Intended purpose	23
	3.2 Indication	26
	3.3 Contraindications	26
	3.4 Patient target group	26
4	Warnings	27
5	Handling	28
6	Preparation	28
	6.1 Information on instrument preparation	28
	6.2 Preparation at the place of use	29
	6.3 Ultrasound bath (optional)	29
	6.4 Manual cleaning	29
	6.5 Mechanical cleaning	29
	6.6 Drying	31
7	Maintenance, inspection	32
	7.1 Functional check	32
8	Sterilization	33
	8.1 Packaging	33
9	Lifetime	33
10	Storage	34
11	Warranty / Repair	34
12	Disposal	34
13	Service and manufacturer address	35

1 Symbol descriptions

Symbols	Definition
	CE-labelling
	Attention
	Validated Parameters
	Manufacturer
	Lot-description
	Reference code
	Medical device / FDA Prescription device
	Medical device
	Non sterile
	Keep away from sunlight
	Dry storage required
	(Electronic) instruction for use

2 Introduction

By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following.

In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

3 Intended use

3.1 Intended purpose

<p>Cuticle- and Nail instruments; Curtnicle Nippers; Pedicure nippers; Nail Scissors; Nail Nippers</p>	<p>An instrument for cutting fingernails and toenails. It can be bent to match the natural curvature of human nails. The toenail models are usually made more robust to cut thick toenails. The instrument consists of two movable branches, with holes for fingers and thumbs in the handle, and cuts by closing the leaves over the nail. It is a non-invasive, reusable instrument</p>
<p>Operating Scissors; Delicate Scissors; Tenotomy Scissors, fine</p>	<p>A surgical instrument for cutting tissue. It consists of two movable branches, with holes for fingers and thumbs in the handle, and cuts by closing the leaves over the material to be separated. The working end can have different designs. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.</p>
<p>Scissors for plastic surgery</p>	<p>A surgical instrument used to cut tissue during cosmetic surgery. It consists of two moving leaves that are usually provided with finger and thumb grips. The distal end of the blades can have different shapes. It is a reusable instrument intended for transient use.</p>
<p>Incision Scissors; Cartilage Scissors; Dissecting Scissors; Scissors with tungsten carbide edges</p>	<p>A surgical instrument for cutting tissue by dissection and/or cutting. The scissors consist of two movable blades and handles with rings for fingers and thumbs and cut by shearing movements (the blades cut when the sharp blades meet when closed). The instrument is usually made of stainless steel and is available in a variety of sizes and styles. In some designs, the cutting edges may be reinforced with carbide. It is a reusable instrument and intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.</p>
<p>Perineum Scissors; Operating and Gynecological Scissors; Cephalotomy Scissors; Uterine Scissors</p>	<p>A surgical instrument for performing tissue sections in the area of the female reproductive system during gynaecological procedures. It consists of two movable blades, which are usually provided with finger and thumb ring handles. The distal end of the leaves can be designed differently. It is a reusable instrument intended for transient use.</p>

Umbilical Scissors	An obstetric surgical instrument for cutting the umbilical cord after birth. It consists of two connected branches with a thumb and finger ring grip which cut by a shearing motion (i.e., the leaves intersect when the sharp edges meet on closing). It is usually made of stainless steel and is available in different designs and sizes. It is a reusable instrument intended for transient use.
Lobectomy Scissors; Vascular Scissors	A surgical instrument for cutting tissue in the thoracic cavity. It consists of two movable blades, which are usually provided with finger and thumb ring handles. The distal end of the leaves can be designed differently. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.
Tissue Scissors; Anatomy Scissors	A surgical instrument for cutting tissue during a surgical procedure. It consists of two movable blades, usually with ring handles for the finger and thumb. The distal end of the blades may have different designs. It is a reusable instrument and intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.
Nasal Scissors (Tube shaft); Septum Scissors; Nasal Scissors	A surgical instrument for cutting tissue during an ENT procedure or during plastic surgery on the nose and its adjacent areas. The scissors consists of two movable blades, which are usually provided with finger and thumb ring handles. The distal end of the leaves can be designed differently. It is a reusable instrument intended for transient use.
Tonsil Scissors	A surgical instrument for cutting tonsillar tissue. It consists of two movable scissor blades, usually with holes for fingers and thumbs in the handles. The distal end of the scissor blades can have different designs. It is a reusable instrument intended for transient use.
Ligature Scissors	A surgical instrument for cutting or removing sutures. It consists of two leaves, usually with finger and thumb ring grips. The distal end of the leaves can be designed differently. It is a pair of scissors for general use, with some surgical scissors having a particular use. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.
Wire cutting scissors; Wire cutting scissors with TC	An instrument for external cutting of wires in various procedures. . The cutting blades will be more robust than ordinary scissors to tolerate the forces required to cut wire. It is a non-invasive, reusable instrument
Nail Splitting Scissors	A instrument for cutting and dividing a finger or toe nail, usually for removal. It consists of two swiveling branches with ring handles for thumb and fingers. The lower leaf is forced under the nail to be cut and the leaves are closed. These cutting leaves are usually very short, pointed and robust. It is a non-invasive, reusable instrument
Enterotomy Scissors	A surgical instrument for cutting tissue during a rectal procedure. It consists of two moving leaves, usually with holes for the fingers. The distal end of the leaves can be designed differently. It is a reusable instrument intended for transient use.
Bandage Scissors	An instrument for cutting dressings. It consists of two movable branches, with finger and thumb holes in the handle. The distal end can have different designs and is usually rounded for safety. The scissors are typically made of stainless steel and the blades can be reinforced with carbide. It is a non-invasive, reusable instrument.

Plaster tear-off pliers; Plaster scissors	A large, scissor-type, hand-held instrument whose leaves, distal to the end, merge into various mouth shapes suitable for cutting thick layers of plaster or synthetic material used for a plaster cast on various body parts. The tip of one of the distal blades of the scissor-type cutting mechanism is typically rounded at its lower edge to prevent injury to the patient when cutting the cast. The product is usually made of stainless steel and usually has a spring mechanism that helps to open the mouth. It is a non-invasive, reusable instrument.
Bone Shears; Rib Shears	A surgical tool to cut bone or cartilage tissue during orthopedic surgery. It consists of two movable blades, which are usually provided with finger and thumb ring handles. The distal end of the leaves can be designed differently. It is a reusable instrument intended for transient use.
Micro Scissors	A narrow knife with a tapering blade on a narrow handle, usually for cutting tendons. This general product group is sometimes used for fine scissors for cutting tendons, which are often provided with a spring to keep the scissors open. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.
Secondary Cataract Scissors; Delicate Eye Scissors	A surgical eye instrument used to cut either front or back segment tissue during eye surgery. Typically it has: 1) two movable blades, normally provided with a finger thumb-ring handle, that cut with a shearing function (i.e., cut the sharpened blades while they are being closed); or 2) an instrument tip, with blades at the distal end, which is inserted into a suitable handle by which the blades are operated. It is typically made of high quality stainless steel and plastic, and is available in various sizes. It is a reusable instrument intended for transient use.
Corneal Scissors; Iridectomy Scissors	A surgical eye instrument used to cut tissue during cataract surgery. It consists of two moving leaves, normally provided with a finger-thumb-ring handle, which cuts with a shearing function (i.e., cutting the sharpened leaves as they are closed). Typically, it is made in different designs and sizes, made of high quality stainless steel. It is a reusable instrument intended for transient use.
Iris Scissors	A surgical eye instrument used to cut tissue during eye surgery, involving the iris. It consists of two moving leaves, normally provided with a finger thumb-ring handle, which cut with a shearing function (i.e., cut the sharpened leaves while they are being closed). Typically, it is made in different designs and sizes, made of high quality stainless steel. It is a reusable instrument intended for transient use.

3.2 Indication

Scissors are used in the following areas:

Umbilical Scissors	The instrument is used during childbirth. It may only be used by trained and qualified professionals.
Perineum Scissors; Operating and Gynecological Scissors; Cephalotomy Scissors; Uterine Scissors	The instrument is used for surgical procedures in gynecology. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Bone Shears; Rib Shears	The instrument is used in various surgical procedures in orthopedics. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Scissors for plastic surgery	The instrument is used in plastic surgery. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Enterotomy Scissors	The instrument is used during rectal procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Cuticle- and Nail instruments; Cuticle Nippers; Pedicure nippers; Nail Scissors; Nail Nippers; Nail Splitting Scissors	The instrument is used for nail and foot care.
Nasal Scissors (Tube shaft); Septum Scissors; Nasal Scissors; Tonsil Scissors	The instrument is used during surgical ENT procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Wire cutting scissors; Wire cutting scissors with TC; Bandage Scissors; Plaster tear-off pliers; Plaster scissors	The instrument is used for treatments. It may only be used by trained and qualified professionals.
Secondary Cataract Scissors; Delicate Eye Scissors; Corneal Scissors; Iridectomy Scissors; Iris Scissors	The instrument is used for surgical procedures in ophthalmology. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Operating Scissors; Delicate Scissors; Tenotomy Scissors, fine; Incision Scissors; Cartilage Scissors; Dissecting Scissors; Scissors with tungsten carbide edges; Lobectomy Scissors; Vascular Scissors; Tissue Scissors; Anatomy Scissors; Ligature Scissors; Micro Scissors	The instrument is used in various surgical procedures. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system. The procedure may only be carried out by suitably trained and qualified specialist personnel.

3.3 Contraindications

No contraindications are known.

3.4 Patient target group

The products are suitable for all patient groups.

4 Warnings

	Medical products are delivered in a non-sterile condition and must be cleaned, disinfected and sterilized prior to their initial use.
	The use of faulty instruments is in principle forbidden and they have to go through the whole cleaning process before return.
	Please take into consideration that through higher power a bigger damage of the tissue can result: f.e. on forceps: the power at the end of the jaw is higher than at the tip of the jaw
	Remove all protective sleeves and films prior to first using or preparation for use.
	The safe combination of different products or of products with implants must be reviewed prior to clinical application by the user.
	Avoid improper throwing or dropping of instruments.
	Avoid mechanical overstressing of the instrument beyond the structural design, this can lead to breakage and deformation!
	A visual inspection of the instrument for damage and contamination must be carried out before each use!
	To prevent all contact corrosion, instruments with damaged surfaces must be separated immediately.
	If the products are used on patients with transmissible spongiform encephalopathy or HIV infection, we decline any responsibility for their reuse.
	After ophthalmical use, please pay attention to water quality during treatment (according to the specifications of AAMI TIR34 and the recommendations of the Rober-Koch-Institute on preparation of medical devices)!
	Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

5 Handling

The type of treatment must be determined in each individual case by the surgeon in cooperation with the internist and the anaesthetist.

For operational use in various surgical disciplines must be done by appropriately trained and qualified personnel.

6 Preparation

The person in charge of preparatory treatment is responsible for ensuring that the treatment is duly carried out using the relevant equipment, materials and personnel in the treatment facility and so achieves the desired result. This necessitates validation and routine monitoring of the process used. We urge you to take note of the national regulations dealing with instrument preparation.

The validated parameters refer to reusable surgical instruments. The validated parameters should also be observed for the other products described, unless a different procedure is explicitly described.

6.1 Information on instrument preparation

- Use cleaning and/or disinfection agents with a pH-value within 9-10.
- Please observe manufacturer instructions regarding dosage, exposure time and renewal of solutions.
- Use suitable brushes that do not remove material, i.e. **no** hard brushes (such as metal brushes and metal sponges) or coarse abrasive cleaners.
- Never leave instruments in cleaning or disinfection agents for longer than the specified time.
- Rinse and dry carefully through channels and pipes.
- Sensitive instruments must be cleaned in a storage or clamping fixture.
- Observe manufacturer instructions of cleaning – and sterilizing equipment.

6.1.1 Water quality

For cleaning, neutralization and rinsing, we recommend the use of demineralized water in accordance with the "Guideline DGKH, DGSV, AKI for the validation and routine monitoring of automated cleaning and thermal disinfection processes for medical devices and on the principles of device selection" (the guideline refers to DIN EN ISO 15883-1 point 6.4.2).

According to the guideline, the following values are recommended:

- Conductivity: $\leq 15 \mu\text{s}/\text{cm}$ (deviating from the table in DIN EN 285)
- pH value: 5 - 7
- Total hardness: $\leq 0.02 \text{ mmol CaO}/\text{l}$
- Salt content: $\leq 10 \text{ mg}/\text{l}$
- Phosphate (as P_2O_5): $\leq 0.5 \text{ mg}/\text{l}$
- Silicate (as SiO_2): $\leq 1 \text{ mg}/\text{l}$

■ Chloride: ≤ 2 mg/l

6.2 Preparation at the place of use

Immediately after use, remove coarse dirt from the instruments with a suitable brush (see chapter 6.1) and rinse out the working cannulas. Do not use fixing agents or hot water ($> 40^{\circ}\text{C}$), as this results in residues becoming fixed and can affect the success of the subsequent cleaning operation

Dismantle and/or open instruments as far as possible. Within short time after use the instruments clean the instruments for reducing a drying of the residues.

This enables an easier cleaning. If instruments come into contact with corroding medicines or cleaning agents, wash these up with water immediately after use.

Longer drying times, e.g. for dry disposal are not validated and therefore not recommended.

V The drying time during validation was 1 hour.

6.3 Ultrasound bath (optional)

All instruments must be opened, dismantled and any cavities rinsed through.

Place instruments in the screen basket in such a way that overlaps and contact between instruments are avoided. Add cleaning agent to the water and adjust the temperature of the solution in line with the cleaning agent manufacturer's instructions.

The cleaning in the ultrasound bath should be at **35-40 kHz, 5 minutes** at least.

V To validate cleaning in an ultrasonic bath, the test items were ultrasonically treated in Neodisher Mediclean forte 0,5 % for 5 minutes.

Subsequently rinse instruments including all cavities before cleaning and disinfection.

The ultrasonic bath should not be used for medical devices where material could be removed in the ultrasonic bath, e.g. products made of soft plastic or coated instruments.

6.4 Manual cleaning


Since mechanical processes can be standardized, reproduced and therefore validated, mechanical cleaning/disinfection should be preferred to manual processes. Manual cleaning and disinfection process is not validated and therefore needs to be validated additionally by the end user.

6.5 Mechanical cleaning



Due to international standards (EN ISO 15883) and national guidelines, only validated mechanical cleaning and disinfection procedures should be used. For automated cleaning, we recommend a standard program for surgical instruments, e.g. instruments from Miele.

Observe the following by loading:

- Place the disassembled or opened instruments securely in the tray.
- Instruments with openings and recesses must be placed with the open side facing downwards so that they can be cleaned and no water from the cleaning process can collect in them.
- If available, use a coordinated rinsing device.

- Do not overload trays, avoid creating any overlaps

Preliminary rinsing is followed by chemical.

The chemical cleaning should take place at **40°C -60°C for at least 5 minutes**.

We recommend products with a **pH-value within 9-10**, e.g. Neodisher MediClean forte from Dr. Weigert. The cleaning agents used should be selected depending on the material and properties of the instruments and in accordance to national regulations: If there is a high chloride concentration in the water, pitting and tension crack corrosion can occur on the instruments. The occurrence of this type of corrosion is minimized by using alkaline cleaning agents and demineralized water. By adding an acid-based neutralization medium, the rinsing off of alkaline cleaning agent residues is facilitated during the first intermediate rinsing process.

In order to prevent the formation of deposits, it is advisable to use neutral cleaners where the water quality is unfavourable.

After the second intermediate rinsing process, thermal disinfection takes place.

The thermal disinfection should be carried out with demineralized water that corresponds to microbiological drinking water quality (< 100 CFU/ml drinking water) at **80 to 95°C** and an expo-

sure time in accordance with **EN ISO 15883**.

After the finishing of the programme take the good out of the machine because corrosion can arise if the instrument remains in the machine.

V Parameters used for the validation of preparation	
Pre-rinsing	1 minute with cold tap water
Cleaning	Temperature: 55 °C
	Soaking Time: 5 minutes (worst case)
	Neodischer Mediclean forte 0,4% (worst case)
Neutralization	Temperature: cold DI water
	Soaking Time: 2 minutes
	Neodisher Z 0,1%
Post-rinsing	2 minutes with cold DI water
Disinfection	Temperature: 90 °C (A ₀ 3000)
	Soaking Time: 5 minutes

6.6 Drying

Ensure adequate drying by the cleaning and disinfection device or using other suitable measures.

V Drying was omitted in the validation (worst case condition).

7 Maintenance, inspection

After cooling to room temperature, the instruments must be visually inspected for protein residues and other contamination. Slits, barriers, locks, tubes and other areas that are difficult to access must be thoroughly inspected. Instruments that are not residue-free must be repeatedly subjected to the entire reprocessing process.

To ensure that surgical instruments can be used for their intended purpose after reprocessing, it is necessary to carry out a functional test after cleaning, disinfection and drying of the visual inspection and care measures. Carry out the functional tests described in chapter 7.1.

Instruments that are stained, blunt, bent, no longer functional, damaged in any other way or the instrument labeling is no longer recognizable (e.g. labeling, UDI code no longer legible) must be discarded!

To help identify faulty instruments that need to be sorted out, we recommend the brochure "Instrument Reprocessing" from the Working Group "Instrumenten Aufbereitung". This includes Chapter 8 "Checks and Care" and Chapter 12 "Surface Changes: Deposits, Discoloration, Corrosion, Aging, Swelling and Stress Cracks".

7.1 Functional check

A newly purchased product must be subjected to a thorough visual and function check after its delivery and before each use.

Products must be checked for irregularities. Paying attention to cracks, fractures and the occurrence of corrosion.

If there are joints, the instruments should be oiled with a care product before the functional test. We recommend a medical white oil based on paraffin oil.

Check instruments with joints for ease of movement.

The labeling must also be checked for legibility.

Carry out a function check in accordance with the intended application of the instrument.

Essential tests for Scissors are among others:

- Correct opening and closing (smooth movement)
- Smooth cutting edges
- Intact tips
- No "scraping" or grinding noise when closing

Defective products must not be used and must have undergone the complete preparatory treatment process again before being returned.

8 Sterilization

Prior to sterilization, products must undergo cleaning and disinfection, be rinsed off without residue using demineralized water and subsequently dried. HEBUmedical recommends using a validated steam sterilization process (e.g. sterilizer in compliance with EN 285 and validated in accordance with DIN EN ISO 17665-1).

The validated parameters refer to reusable surgical instruments. The validated parameters should also be observed for the other products described, unless a different procedure is explicitly described.

On using the fractionated vacuum method, **sterilization** must be performed with at least **134°C (USA 132°C)** with a **minimum dwell period of 3 minutes**. Vacuum drying must then be carried out for at least 20 minutes.

V Parameters used for the validation of steam sterilization	
Prevacuum	3 times
Sterilization temperature	132 °C
Sterilization time	1,5 minutes (half cycle method)
Drying time	Drying time

The vapour must be free of ingredients, recommended limiting values of feed water and vapour condensate are determined through EN 285.

Other sterilization processes are compatible but not validated from HEBUmedical.

When loading, observe the recommended total weight. After the sterilization, check the sterile product packaging for damage, and inspect the sterilization indicators.

8.1 Packaging

Compliant packaging of products for sterilization in line with ISO 11607. Packaging used must be suitable for the instruments and protect them from microbiological contamination during storage. The seal must not be under tension. HEBUmedical recommends container or hospital common sterilization paper/film packagings for sterilized packaging.

V During validation the instruments were packaged in hospital common sterilization packagings (paper/film packagings) and steam sterilized.

9 Lifetime

The steam sterilization procedure was validated by laboratory tests. The products were sterile validated at a pre-vacuum of at least 5min duration and a temperature of 134°C for a lifetime of 50 cycles.

You can continue to use the instruments at your own responsibility over this cycle value if the tests described in chapter 7 have been successfully completed.

10 Storage

Store products in a dry, clean and dust-free environment at moderate temperatures from 5°C to 40°C.

Protect from the effects of the sun's rays and artificial light.

11 Warranty / Repair

Our products are manufactured from high-grade materials and carefully checked prior to dispatch. However, even if used properly in accordance with their intended purpose they are subject to a greater or lesser degree of wear depending on their intensity of use.

This wear is technically induced and unavoidable.

Should faults occur independently of wear, please contact our customer services. Defective products should no longer be used.

They must undergo the complete preparatory treatment process before being returned.

12 Disposal

Before disposal, the product must undergo the complete reprocessing and sterilization process in accordance with the procedure described in these instructions for use.

Risk of infection from contaminated products!

- When disposing of the product, its components and their packaging, the national regulations must be observed.

Risk of injury from sharp-edged and/or pointed products!

- When disposing of the product, ensure that the packaging prevents injury from the product.

13 Service and manufacturer address

Should you require the instructions for use in paper form, please use the contact details below. The instructions for use in paper form will be made available to you within seven calendar days of receipt of the request.

Alternatively, you can print out the electronic instructions for use yourself.

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen / Germany
Tel. +49 7461 94 71 - 0
Fax +49 7461 94 71 - 22
eMail: service@HEBUmedical.de
Website: www.HEBUmedical.de

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