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Table of contents

1	Symbol descriptions	22
2	Introduction	23
3	Intended use	23
	3.1 Intended purpose	23
	3.2 Indication	25
	3.3 Contraindications	26
	3.4 Patient target group	26
4	Warnings	27
5	Handling	28
6	Preparation	28
	6.1 Information on instrument preparation	28
	6.2 Preparation at the place of use	29
	6.3 Ultrasound bath (optional)	29
	6.4 Manual cleaning	29
	6.5 Mechanical cleaning	29
	6.6 Drying	31
7	Maintenance, inspection	31
	7.1 Functional check	32
8	Sterilization	32
	8.1 Packaging	33
9	Lifetime	33
10	Storage	33
11	Warranty / Repair	33
12	Disposal	34
13	Service and manufacturer address	34

1 Symbol descriptions

Symbols	Definition
	CE-labelling
	Attention
	Validated Parameters
	Manufacturer
	Lot-description
	Reference code
	Medical device / FDA Prescription device
	Medical device
	Non sterile
	Keep away from sunlight
	Dry storage required
	(Electronic) instruction for use

2 Introduction

By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following.

In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

3 Intended use

3.1 Intended purpose

Hooks; Exploration- and Coagulation Hooks; Root Retractor; Spreaders; Tissue hook	A shaft-like surgical instrument, having a proximal end forming a handle which can be formed in various embodiments. The shaft at the distal end merges into one or more parallel hooks to hold tissue. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.
Skin Hook; Micro Hooks; reparatory, spoons, forks	A shaft-like surgical instrument, having a proximal end forming a handle which can be formed in various embodiments. The shaft at the distal end merges into one or more parallel hooks that apply tension to the skin. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.
Rectal Specula	An instrument, usually composed of two articulated branches, available in a variety of sizes, shapes and contours. It has a blunt distal end and is used for insertion and opening to stretch or stretch the rectal opening (canal). It is used to study or introduce other instruments. It is usually made of stainless steel. The portion of the branches in the proximal area forms the handle while the distal area opens when you squeeze the handle. It is a reusable, invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.
Plaster spreader	A hand instrument with specially designed branches for separating and spreading plaster casts on limbs or body parts. It is designed like large pliers with flat, wedge-shaped cutting blades serrated on the outer blade, which may be interchangeable. These blades are placed around the plaster cast to be cut and the plaster is separated by pressing the handles together. The instrument is usually made of stainless steel. It is a non-invasive, reusable product.
Retractors	A surgical hand instrument used to hold the edges of the wound during an opening of the abdominal cavity. It is commonly used during deep abdominal surgery, e.g. to make organs accessible or expose them. It consists of a blunt, flat, curved blade and a handle. The surgeon uses the instrument to keep the margins open manually. It is made of stainless steel and is available in different blade lengths. It is a reusable instrument intended for transient use.

Bone levers	A surgical hand instrument for lifting, positioning or levering of bone structures, other anatomical structures or surgical material during an orthopedic procedure. It is made of stainless steel and is available in a variety of sizes, styles and working ends. It is usually a sturdy instrument with a blunt working end that can be pointed, rounded, flat, curved or hooked. The handle at the proximal end is also available in different versions. Some models have special functions, e.g. for ankle surgery, or to treat a dislocation. It is a reusable instrument intended for transient use.
Elevators; Meniscus retractors	A hand-held surgical instrument used to lift, position, or lever open bone structures, other anatomical structures or surgical material during an orthopedic surgery. It is made of stainless steel and is available in a variety of sizes, styles and working ends. It is a reusable instrument intended for transient use.
Spinal retractors	A shaft-like surgical instrument whose proximal end forms a handle, which can be shaped in various designs. At the distal end, the shaft merges into one or more parallel hooks for the excision of tissue and tendons on the vertebral bone. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.
Trachea Dilators	A hand surgical instrument used to dilate tracheal structures e.g. to assist in the insertion of a tracheal cannula. It can have different designs and sizes, e.g. a slender hollow or solid body made of metal, plastic or other suitable materials in cylindrical form, or a scissors-like instrument that activates by closing the two ring handles two or three blunt branches in opposite directions. It is a reusable, invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.
Thyroid Gland Retractor; Trachea Retractors; Micro-laryngeal Instruments	A surgical instrument used to temporarily hold the trachea in a stable position during a tracheostomy and/or to temporarily separate/push away the edges of a tracheostomy to allow placement of a tracheostomy tube. The instrument can have different designs, e.g. single- or double-ended, with a wide retractor on one side and opposite a pointed hook or several parallel hooks, which can be angled differently. It is made of stainless steel. It is a reusable instrument intended for transient use.
Epiglottis Retractor; Foreign Body Levers	An instrument for gripping, holding or manipulating anatomical structures during an ENT procedure at e.g. the bronchi, esophagus, trachea, larynx, pharynx, nose, or ear. It is made of stainless steel. It is a reusable, invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.
Intestinal and Stomach Clamp; Intestinal Clamp; Gastrointestinal clamp; Pyloric spreader	A surgical instrument for atraumatically gripping, compressing, connecting or holding the intestine during gastrointestinal procedures. It is a reusable instrument intended for transient use.
Bladder Retractors; Prostatic Retractors	A frame-shaped device to be hooked into the retractor to hold the anatomical structures of the bladder. Retractors are either blunt or hook-shaped (sometimes with multiple hooks). The retractor blades are used to hold back the edges of e.g. incisions, an organ or tissue during a surgical procedure. It is a reusable instrument intended for transient use.

Lateral Vaginal Retractors; Vaginal Specula	A product used during the examination and treatment of the vagina to provide access to the inside. Usually this is achieved by stretching the opening. It is a reusable, invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.
Nasal specula (simple); ear specula	An instrument, usually composed of two articulated branches, available in a variety of sizes, shapes and contours. It has a blunt distal end and is used in insertion and opening to stretch or stretch tissue in the nostril. It is used for examination or during procedures. It is usually made of stainless steel with a self-opening mechanism or an adjustable mechanism. It is a reusable instrument intended for transient use. It is a reusable, invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.
Hooks for posterior nares	A shaft-like surgical instrument whose proximal end forms a handle that can vary in shape. The shaft tapers toward the distal end toward a hook or hooks used in plastic nose surgery. In the multi-hook variant, the hooks are arranged in parallel. It is a reusable instrument intended for transient use.
Lid Retractor; Eye Specula; Muscle and fixation hooks	A surgical eye instrument used to hold off the edges and associated tissue of an incision in the eye during a surgical procedure. It usually has a self-holding arm with curved hooks or blades at the distal end that grip and hold the tissue away. The proximal end is usually connected and has a lockable mechanism to hold the working end in an expanded position. It is made of stainless steel and is available in various sizes and designs. It is a reusable instrument intended for transient use.

3.2 Indication

Retractors are used in the following areas:

Rectal Specula	The instrument is used for anal and rectal examinations. The examination may only be carried out by suitably trained and qualified specialist personnel.
Epiglottis Retractor; Foreign Body Levers Trachea Dillators	The instrument is used for ENT examinations and treatments. It may only be used by trained and qualified professionals. The instrument is used in various procedures on the trachea. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Plaster spreader	The instrument is used in orthopedics. It may only be used by trained and qualified professionals.
Lateral Vaginal Retractors; Vaginal Specula Nasal specula (simple); ear specula Spinal retractors	The product is used in examinations and treatments in gynecology. It may only be used by trained and qualified personnel. The product is used in ENT examinations and treatments. It may only be used by trained and qualified professionals. The instrument is used for surgical procedures on the spine. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Intestinal and Stomach Clamp; Intestinal Clamp; Gastrointes-gastrointestinal clamp; Pyloric spreader	The instrument is used for surgical procedures in the area of the gastrointestinal tract. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Lid Retractor; Eye Specula; Muscle and fixation hooks	The instrument is used for surgical procedures in ophthalmology. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Bone levers	The instrument is used in various surgical procedures in orthopedics. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Bladder Retractors; Prostatic Retractors	The instrument is used for urological surgical procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Retractors	The instrument is used in various surgical procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Hooks; Exploration- and Coagulation Hooks; Root Retractor; Spreaders; Tissue hook; Skin Hook; Micro Hooks, rapsatory, spoons, forks	The instrument is used in various surgical procedures. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Chisels and osteotomes; rha-chiotome; flat chisel; Nasal osteotomes; Elevators; Meniscus retractors	The instrument is used in various surgical procedures in orthopedics. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Thyroid Gland Retractor; Trachea Retractors; Micro-laryngeal Instruments; Hooks for posterior nares	The instrument is used during surgical ENT procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.

3.3 Contraindications

No contraindications are known.

3.4 Patient target group

The products are suitable for all patient groups.

4 Warnings

	Medical products are delivered in a non-sterile condition and must be cleaned, disinfected and sterilized prior to their initial use.
	The use of faulty instruments is in principle forbidden and they have to go through the whole cleaning process before return.
	Please take into consideration that through higher power a bigger damage of the tissue can result: f.e. on forceps: the power at the end of the jaw is higher than at the tip of the jaw
	Remove all protective sleeves and films prior to first using or preparation for use.
	The safe combination of different products or of products with implants must be reviewed prior to clinical application by the user.
	Avoid improper throwing or dropping of instruments.
	Avoid mechanical overstressing of the instrument beyond the structural design, this can lead to breakage and deformation!
	A visual inspection of the instrument for damage and contamination must be carried out before each use!
	To prevent all contact corrosion, instruments with damaged surfaces must be separated immediately.
	If the products are used on patients with transmissible spongiform encephalopathy or HIV infection, we decline any responsibility for their reuse.
	After ophthalmical use, please pay attention to water quality during treatment (according to the specifications of AAMI TIR34 and the recommendations of the Rober-Koch-Institute on preparation of medical devices)!
	Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

5 Handling

The type of treatment must be determined in each individual case by the surgeon in cooperation with the internist and the anaesthetist.

For operational use in various surgical disciplines must be done by appropriately trained and qualified personnel.

6 Preparation

The person in charge of preparatory treatment is responsible for ensuring that the treatment is duly carried out using the relevant equipment, materials and personnel in the treatment facility and so achieves the desired result. This necessitates validation and routine monitoring of the process used. We urge you to take note of the national regulations dealing with instrument preparation.

The validated parameters refer to reusable surgical instruments. The validated parameters should also be observed for the other products described, unless a different procedure is explicitly described.

6.1 Information on instrument preparation

- Use cleaning and/or disinfection agents with a pH-value within 9-10.
- Please observe manufacturer instructions regarding dosage, exposure time and renewal of solutions.
- Use suitable brushes that do not remove material, i.e. **no** hard brushes (such as metal brushes and metal sponges) or coarse abrasive cleaners.
- Never leave instruments in cleaning or disinfection agents for longer than the specified time.
- Rinse and dry carefully through channels and pipes.
- Sensitive instruments must be cleaned in a storage or clamping fixture.
- Observe manufacturer instructions of cleaning – and sterilizing equipment.

6.1.1 Water quality

For cleaning, neutralization and rinsing, we recommend the use of demineralized water in accordance with the "Guideline DGKH, DGSV, AKI for the validation and routine monitoring of automated cleaning and thermal disinfection processes for medical devices and on the principles of device selection" (the guideline refers to DIN EN ISO 15883-1 point 6.4.2).

According to the guideline, the following values are recommended:

- Conductivity: $\leq 15 \mu\text{s}/\text{cm}$ (deviating from the table in DIN EN 285)
- pH value: 5 - 7
- Total hardness: $\leq 0.02 \text{ mmol CaO}/\text{l}$
- Salt content: $\leq 10 \text{ mg}/\text{l}$
- Phosphate (as P_2O_5): $\leq 0.5 \text{ mg}/\text{l}$
- Silicate (as SiO_2): $\leq 1 \text{ mg}/\text{l}$

- Chloride: ≤ 2 mg/l

6.2 Preparation at the place of use

Immediately after use, remove coarse dirt from the instruments with a suitable brush (see chapter 6.1) and rinse out the working cannulas. Do not use fixing agents or hot water ($> 40^{\circ}\text{C}$), as this results in residues becoming fixed and can affect the success of the subsequent cleaning operation

Dismantle and/or open instruments as far as possible. Within short time after use the instruments clean the instruments for reducing a drying of the residues.

This enables an easier cleaning. If instruments come into contact with corroding medicines or cleaning agents, wash these up with water immediately after use.

Longer drying times, e.g. for dry disposal are not validated and therefore not recommended.

V The drying time during validation was 1 hour.

6.3 Ultrasound bath (optional)

All instruments must be opened, dismantled and any cavities rinsed through.

Place instruments in the screen basket in such a way that overlaps and contact between instruments are avoided. Add cleaning agent to the water and adjust the temperature of the solution in line with the cleaning agent manufacturer's instructions.

The cleaning in the ultrasound bath should be at **35-40 kHz, 5 minutes** at least.

V To validate cleaning in an ultrasonic bath, the test items were ultrasonically treated in Neodisher Mediclean forte 0,5 % for 5 minutes.

Subsequently rinse instruments including all cavities before cleaning and disinfection.

The ultrasonic bath should not be used for medical devices where material could be removed in the ultrasonic bath, e.g. products made of soft plastic or coated instruments.

6.4 Manual cleaning


Since mechanical processes can be standardized, reproduced and therefore validated, mechanical cleaning/disinfection should be preferred to manual processes. Manual cleaning and disinfection process is not validated and therefore needs to be validated additionally by the end user.

6.5 Mechanical cleaning



Due to international standards (EN ISO 15883) and national guidelines, only validated mechanical cleaning and disinfection procedures should be used. For automated cleaning, we recommend a standard program for surgical instruments, e.g. instruments from Miele.

Observe the following by loading:

- Place the disassembled or opened instruments securely in the tray.
- Instruments with openings and recesses must be placed with the open side facing downwards so that they can be cleaned and no water from the cleaning process can collect in them.
- If available, use a coordinated rinsing device.

- Do not overload trays, avoid creating any overlaps

Preliminary rinsing is followed by chemical.

The chemical cleaning should take place at **40°C -60°C for at least 5 minutes**.

We recommend products with a **pH-value within 9-10**, e.g. Neodisher MediClean forte from Dr. Weigert. The cleaning agents used should be selected depending on the material and properties of the instruments and in accordance to national regulations: If there is a high chloride concentration in the water, pitting and tension crack corrosion can occur on the instruments. The occurrence of this type of corrosion is minimized by using alkaline cleaning agents and demineralized water. By adding an acid-based neutralization medium, the rinsing off of alkaline cleaning agent residues is facilitated during the first intermediate rinsing process.

In order to prevent the formation of deposits, it is advisable to use neutral cleaners where the water quality is unfavourable.

After the second intermediate rinsing process, thermal disinfection takes place.

The thermal disinfection should be carried out with demineralized water that corresponds to microbiological drinking water quality (< 100 CFU/ml drinking water) at **80 to 95°C** and an expo-

sure time in accordance with **EN ISO 15883**.

After the finishing of the programme take the good out of the machine because corrosion can arise if the instrument remains in the machine.

V Parameters used for the validation of preparation	
Pre-rinsing	1 minute with cold tap water
Cleaning	Temperature: 55 °C
	Soaking Time: 5 minutes (worst case)
	Neodischer Mediclean forte 0,4% (worst case)
Neutralization	Temperature: cold DI water
	Soaking Time: 2 minutes
	Neodisher Z 0,1%
Post-rinsing	2 minutes with cold DI water
Disinfection	Temperature: 90 °C (A ₀ 3000)
	Soaking Time: 5 minutes

6.6 Drying

Ensure adequate drying by the cleaning and disinfection device or using other suitable measures.

V Drying was omitted in the validation (worst case condition).

7 Maintenance, inspection

After cooling to room temperature, the instruments must be visually inspected for protein residues and other contamination. Slits, barriers, locks, tubes and other areas that are difficult to access must be thoroughly inspected. Instruments that are not residue-free must be repeatedly subjected to the entire reprocessing process.

To ensure that surgical instruments can be used for their intended purpose after reprocessing, it is necessary to carry out a functional test after cleaning, disinfection and drying of the visual inspection and care measures. Carry out the functional tests described in chapter 7.1.

Instruments that are stained, blunt, bent, no longer functional, damaged in any other way or the instrument labeling is no longer recognizable (e.g. labeling, UDI code no longer legible) must be discarded!

To help identify faulty instruments that need to be sorted out, we recommend the brochure "Ins-

trument Reprocessing" from the Working Group "Instrumenten Aufbereitung". This includes Chapter 8 "Checks and Care" and Chapter 12 "Surface Changes: Deposits, Discoloration, Corrosion, Aging, Swelling and Stress Cracks".

7.1 Functional check

A newly purchased product must be subjected to a thorough visual and function check after its delivery and before each use.

Products must be checked for irregularities. Paying attention to cracks, fractures and the occurrence of corrosion.

If there are joints, the instruments should be oiled with a care product before the functional test. We recommend a medical white oil based on paraffin oil.

Check instruments with joints for ease of movement.

The labeling must also be checked for legibility.

Carry out a function check in accordance with the intended application of the instrument.

Essential tests for instruments with hooks are among others:

- No bent tips and hooks
- If two-piece: Correct opening and closing (smooth movement)

Defective products must not be used and must have undergone the complete preparatory treatment process again before being returned.

8 Sterilization

Prior to sterilization, products must undergo cleaning and disinfection, be rinsed off without residue using demineralized water and subsequently dried. HEBUmedical recommends using a validated steam sterilization process (e.g. sterilizer in compliance with EN 285 and validated in accordance with DIN EN ISO 17665-1).

The validated parameters refer to reusable surgical instruments. The validated parameters should also be observed for the other products described, unless a different procedure is explicitly described.

On using the fractionated vacuum method, **sterilization** must be performed with at least **134°C (USA 132°C)** with a **minimum dwell period of 3 minutes**. Vacuum drying must then be carried out for at least 20 minutes.

V Parameters used for the validation of steam sterilization	
Prevacuum	3 times
Sterilization temperature	132 °C
Sterilization time	1,5 minutes (half cycle method)
Drying time	Drying time

The vapour must be free of ingredients, recommended limiting values of feed water and vapour condensate are determined through EN 285.

Other sterilization processes are compatible but not validated from HEBUmedical. When loading, observe the recommended total weight. After the sterilization, check the sterile product packaging for damage, and inspect the sterilization indicators.

8.1 Packaging

Compliant packaging of products for sterilization in line with ISO 11607. Packaging used must be suitable for the instruments and protect them from microbiological contamination during storage. The seal must not be under tension. HEBUmedical recommends container or hospital common sterilization paper/film packaging for sterilized packaging.

V During validation the instruments were packaged in hospital common sterilization packaging (paper/film packaging) and steam sterilized.

9 Lifetime

The steam sterilization procedure was validated by laboratory tests. The products were sterile validated at a pre-vacuum of at least 5min duration and a temperature of 134°C for a lifetime of 50 cycles.

You can continue to use the instruments at your own responsibility over this cycle value if the tests described in chapter 7 have been successfully completed.

10 Storage

Store products in a dry, clean and dust-free environment at moderate temperatures from 5°C to 40°C.

Protect from the effects of the sun's rays and artificial light.

11 Warranty / Repair

Our products are manufactured from high-grade materials and carefully checked prior to dispatch. However, even if used properly in accordance with their intended purpose they are subject to a greater or lesser degree of wear depending on their intensity of use.

This wear is technically induced and unavoidable.

Should faults occur independently of wear, please contact our customer services. Defective products should no longer be used.

They must undergo the complete preparatory treatment process before being returned.

12 Disposal

Before disposal, the product must undergo the complete reprocessing and sterilization process in accordance with the procedure described in these instructions for use.

Risk of infection from contaminated products!

- When disposing of the product, its components and their packaging, the national regulations must be observed.

Risk of injury from sharp-edged and/or pointed products!

- When disposing of the product, ensure that the packaging prevents injury from the product.

13 Service and manufacturer address

Should you require the instructions for use in paper form, please use the contact details below. The instructions for use in paper form will be made available to you within seven calendar days of receipt of the request.

Alternatively, you can print out the electronic instructions for use yourself.

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen / Germany
Tel. +49 7461 94 71 - 0
Fax +49 7461 94 71 - 22
eMail: service@HEBUmedical.de
Website: www.HEBUmedical.de

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