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# Contents

1	Symbol descriptions .....	10
2	Introduction .....	10
3	Scope .....	11
3.1	Intended Use .....	11
3.2	Contraindications.....	11
4	Warning.....	11
5	Handling .....	11
6	Preparation.....	12
6.1	Reutilization restrictions .....	12
6.2	Information on instrument preparation .....	12
7	Maintenance, inspection .....	12
8	Packaging .....	12
9	Storage.....	12
10	Warranty / Repair .....	12
11	Manufacturer and servicing address.....	13

# 1 Symbol descriptions

Symbol	Definition
	CE-labelling
	Attention
	Validated Parameters
	Manufacturer
	Lot-description
	Reference code
	Medical device / FDA Prescription device
	Medical device
	Non sterile
	Keep away from sunlight
	Dry storage required
	(Electronic) instruction for use

# 2 Introduction

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them. By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following. In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

### 3 Scope

The use of headlamps, head magnifiers and reflex mirrors takes place in almost every medical field.

#### 3.1 Intended Use


Head Mirrors and head band: A medical device (a lamp) that the doctor can wear on his head. The lamp is mounted on a band or helmet and placed on the forehead to bring light directly into view during a surgical procedure, examination or therapeutic treatment. The lamp usually consists of a magnifying glass, a reflector and a connector for the cold-light fiber optic cable, or a power supply from a battery. It is a reusable product.

Loupes binoculaires: An ophthalmic product worn on the head, with a round, concave mirror mounted on a headband, for projecting a beam of light to allow examination of the eye and its surrounding structures. It is a reusable product.

#### 3.2 Contraindications

No contraindications are known.

### 4 Warning

	Medical products are delivered in a non-sterile condition and must be cleaned and disinfected prior to their initial use.
	The use of faulty instruments is in principle forbidden and they must go through the whole cleaning process before return.
	Please observe the additional information enclosed with the products.
	Remove all protective sleeves and films prior to first using or preparation for use.
	The safe combination of different products must be reviewed prior to clinical application by the user.
	Avoid mechanical overstressing of the instrument beyond the structural design, this can lead to breakage and deformation!
	Before each use, the instrument must be checked for function, damage and contamination!
	Avoid improper throwing or dropping of products.
	If the products are used on patients with transmissible spongiform encephalopathy or HIV infection, we decline any responsibility for their reuse.

### 5 Handling

The product may only be used by appropriately trained and qualified personnel.

## 6 Preparation

The person in charge of preparatory treatment is responsible for ensuring that the treatment is duly carried out using the relevant equipment, materials and personnel in the treatment facility and so achieves the desired result. This necessitates validation and routine monitoring of the process used. We urge you to take note of the national regulations dealing with instrument preparation.

### 6.1 Reutilization restrictions

Frequent repeat preparatory treatment has minimal effects on the product. The end of the product life is normally determined by wear and damage due to use

### 6.2 Information on instrument preparation

- Use cleaning and/or disinfection agents with a pH-value within 9-10. Please observe manufacturer instructions regarding dosage, exposure time and renewal of solutions.
- Do not use hard brushes or coarse abrasive cleaners.
- Never leave instruments in cleaning or disinfection agents for longer than the specified time.
- Only used demineralized water for rinsing
- The preparation is carried out by a wipe disinfectant.

## 7 Maintenance, inspection

After cooling to room temperature, the instruments must be visually inspected for protein residues and other contamination, paying particular attention to cavities, blocks, inclusions, pipes, and other inaccessible areas. Instruments which are not free of residues must be returned for a complete retreatment process.

Carry out the functional check mentioned above. Instruments with stains, which are blunt, bent, no longer function or which are otherwise damaged must be segregated!

## 8 Packaging

Compliant packaging of products for sterilization in line with ISO 11607. Packaging used must be suitable for the instruments and protect them from microbiological contamination during storage. The seal must not be under tension. HEBUmedical recommends container for sterilized packaging.

## 9 Storage

Store products in a dry, clean and dust-free environment at moderate temperatures from 5°C to 40°C. Protect from the effects of the sun's rays and artificial light.

## 10 Warranty / Repair

Our products are manufactured from high-grade materials and carefully checked prior to dispatch. However, even if used properly in accordance with their intended purpose they are subject to a greater or lesser degree of wear depending on their intensity of use.

This wear is technically induced and unavoidable.

Should faults occur independently of wear, please contact our customer services. Defective products should no longer be used.

They must undergo the complete preparatory treatment process before being returned.

## 11 Manufacturer and servicing address

Should you require the instructions for use in paper form, please use the contact details below. The instructions for use in paper form will be made available to you within seven calendar days of receipt of the request.

Alternatively, you can print out the electronic instructions for use yourself.

HEBUmedical GmbH  
Badstraße 8  
78532 Tuttlingen / Germany  
Tel. +49 7461 94 71 - 0  
Fax +49 7461 94 71 - 22  
eMail: [service@HEBUmedical.de](mailto:service@HEBUmedical.de)  
Web: [www.HEBUmedical.de](http://www.HEBUmedical.de)

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**Our vast range of products covers over 10 000 various instruments.  
They can be found in Internet or request our catalog.**

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