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Suctions tubes and cannulas / Tuyaux d'aspiration et canettes / Tuberías y latas de succión / Tubi e lattine di aspirazione / Смукателни тръби и кутии / Sugrør og dåser / Iminaitorud ja purgid / Imurputket ja tõlkit / Σωλήνες και δοχεία αναρρόφησης / Usisne cijevi i limenke / Sūkšanas caurules un kārbas / Siurbimo vamzdžiai ir skardinės / Zuigpijpen en blikjes / Sugerør og kanyler / Rury ssące i puszki / Tubos de sucção e latas / Tuburi de aspirație și canule / Sugrör och burkar / Sacie potrubia a plechovky / Sesalne cevi in pločevinke / Sací potrubí a plechovky / Szívócsövek és kannák

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Seite 5

### **Operating Manual**

Page 19

### **Mode d'emploi**

Page 33

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Página 47

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Pagina 61

DEUTSCH

ENGLISH

FRANÇAIS

ESPAÑOL

ITALIANO

---

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# 1 Symbol descriptions

Symbols	Definition
	CE-labelling
	Attention
	Validated Parameters
	Manufacturer
	Lot-description
	Reference code
	Medical device / FDA Prescription device
	Medical device
	Non sterile
	Keep away from sunlight
	Dry storage required
  Hinweis auf eIFU	(Electronic) instruction for use

## 2 Introduction

By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following.

In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

Our products are intended exclusively for professional use in accordance with their intended purpose by appropriately trained and qualified personnel and may only be purchased by such personnel.

## 3 Intended use

The list of products for which these instructions for use apply can be found in the appendix.

### 3.1 Intended Purpose

Rectal probes; Suction Tubes, cannula	A product for periodic or continuous irrigation and/or suction of body cavities or wounds during a surgical procedure, without any significant change in energy, substances, or parameters. It is used to remove loose tissue, fluids, and other debris from the surgical field to keep it clear and improve visibility and access. It is not suitable for the delivery of parenteral fluids, administration of drugs, or other life-sustaining measures. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system. It is a reusable product and is intended for temporary use.
Biopsy Cannulas	A narrow, tapered metal tube used to take tissue samples for analysis in order to determine possible follow-up measures that could prevent suffering or prolong the patient's life. It is made in different diameters, the cannula tip is made in different versions (for example, straight or angled). It is a reusable instrument intended for transient use.
Metal Catheters	A cylindrical, rod-like surgical instrument for probing and/or dilating the urethra to prevent patient suffering. It has a semi-rigid, wedge-shaped shaft that ends at the distal end in a narrow, beaded tip. It can be hollow to pass other narrow instruments. It is a reusable instrument intended for transient use.

### 3.2 Indication

The use of suction tubes and cannulas take place in almost every area of surgery, including: Abdominal surgery, anesthesia, ENT, Hematology, Urology

### 3.3 Contraindications

No contraindications are known.

### 3.4 Patient target group

The products are suitable for all patient groups.

## 4 Warnings

	Medical products are delivered in a non-sterile condition and must be cleaned, disinfected and sterilized prior to their initial use.
	The use of faulty instruments is in principle forbidden and they have to go through the whole cleaning process before return.
	Please take into consideration that through higher power a bigger damage of the tissue can result: f.e. on forceps: the power at the end of the jaw is higher than at the tip of the jaw
	Remove all protective sleeves and films prior to first using or preparation for use.
	The safe combination of different products or of products with implants must be reviewed prior to clinical application by the user.
	Avoid improper throwing or dropping of instruments.
	Avoid mechanical overstressing of the instrument beyond the structural design, this can lead to breakage and deformation!
	A visual inspection of the instrument for damage and contamination must be carried out before each use!
	To prevent all contact corrosion, instruments with damaged surfaces must be separated immediately.
	If the products are used on patients with transmissible spongiform encephalopathy or HIV infection, we decline any responsibility for their reuse.
	After ophthalmical use, please pay attention to water quality during treatment (according to the specifications of AAMI TIR34 and the recommendations of the Rober-Koch-Institute on preparation of medical devices)!
	Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## 5 Usage

### 5.1 Handling

The type of treatment must be determined in each individual case by the surgeon in cooperation with the internist and the anaesthetist.

For operational use in various surgical disciplines must be done by appropriately trained and qualified personnel.

### 5.2 Applicaton / function

Rectal probes; suction tubes, cannulas:

- Connect the appropriate tube to the rectal probe/suction cannula.
- Insert the rectal probe/suction cannula into the area to be suctioned.
- Create negative pressure to enable suction
- For models with thumb plate and breaker hole:
  - Suction: Cover the suction interrupter on the handle plate with your thumb
  - Do not suction: Remove your thumb from the suction interrupter on the handle plate
- Remove the rectal probe/suction cannula from the area and return it to the reprocessing unit as quickly as possible

Biopsy cannulas:

- Insert the cannula into the corresponding area where the tissue sample is to be taken
- Remove the tissue sample
- Remove the biopsy cannula from the area and return it to the reprocessing unit as quickly as possible

Metal catheters:

- Insert the metal catheter into the urethra
- Create negative pressure to enable suction
- Probe the urethra, drain urine if necessary
- Remove the metal catheter from the area and reprocess it as quickly as possible

### 5.3 Variants / performance features

The products have the following performance features:

- Total length and / or working length
- Diameter
- With or without handle plate
- With or without suction interrupter
- Different hose connections (luer lock, luer cone, cone)
- Material: Stainless steel

The products are available in various versions. For detailed information, please refer to our current "HEBU general catalog" (<https://www.hebumedical.com>).

## 6 Preparation

The person in charge of preparatory treatment is responsible for ensuring that the treatment is duly carried out using the relevant equipment, materials and personnel in the treatment facility and so achieves the desired result. This necessitates validation and routine monitoring of the process used. We urge you to take note of the national regulations dealing with instrument preparation.

The validated parameters refer to reusable surgical instruments. The validated parameters should also be observed for the other products described, unless a different procedure is explicitly described.

### 6.1 Information on instrument preparation

- Use cleaning and/or disinfection agents with a pH-value within 9-10. Please observe manufacturer instructions regarding dosage, exposure time and renewal of solutions.
- Use suitable brushes that do not remove material, i.e. **no** hard brushes (such as metal brushes and metal sponges) or coarse abrasive cleaners.
- Never leave instruments in cleaning or disinfection agents for longer than the specified time.
- Rinse and dry carefully through channels and pipes.
- Sensitive instruments must be cleaned in a storage or clamping fixture.
- Observe manufacturer instructions of cleaning – and sterilizing equipment.

#### 6.1.1 Water quality

For cleaning, neutralization and rinsing, we recommend the use of demineralized water in accordance with the "Guideline DGKH, DGSV, AKI for the validation and routine monitoring of automated cleaning and thermal disinfection processes for medical devices and on the principles of device selection" (the guideline refers to DIN EN ISO 15883-1 point 6.4.2).

According to the guideline, the following values are recommended:

- Conductivity:  $\leq 15 \mu\text{s}/\text{cm}$  (deviating from the table in DIN EN 285)
- pH value: 5 - 7
- Total hardness:  $\leq 0.02 \text{ mmol CaO}/\text{l}$
- Salt content:  $\leq 10 \text{ mg}/\text{l}$
- Phosphate (as  $\text{P}_2\text{O}_5$ ):  $\leq 0.5 \text{ mg}/\text{l}$
- Silicate (as  $\text{SiO}_2$ ):  $\leq 1 \text{ mg}/\text{l}$
- Chloride:  $\leq 2 \text{ mg}/\text{l}$

### 6.2 Preparation at the place of use

Immediately after use, remove coarse dirt from the instruments with a suitable brush (see chapter 6.1) and rinse the working cannulas with a syringe or clean them using the wire. Do not use

fixing agents or hot water (> 40°C), as this results in residues becoming fixed and can affect the success of the subsequent cleaning operation

Dismantle and/or open instruments as far as possible. Within short time after use the instruments clean the instruments for reducing a drying of the residues.

This enables an easier cleaning. If instruments come into contact with corroding medicines or cleaning agents, wash these up with water immediately after use.

Longer drying times, e.g. for dry disposal are not validated and therefore not recommended.

**V** The drying time during validation was 1 hour.

### 6.3 Ultrasound bath (optional)

All instruments must be opened, dismantled and any cavities rinsed through.

Place instruments in the screen basket in such a way that overlaps and contact between instruments are avoided. Add cleaning agent to the water and adjust the temperature of the solution in line with the cleaning agent manufacturer's instructions.

The cleaning in the ultrasound bath should be at **35-40 kHz, 5 minutes** at least.

**V** To validate cleaning in an ultrasonic bath, the test items were ultrasonically treated in Neodisher Mediclean forte 0,5% for 5 minutes.

Subsequently rinse instruments including all cavities before cleaning and disinfection.

The ultrasonic bath should not be used for medical devices where material could be removed in the ultrasonic bath, e.g. products made of soft plastic or coated instruments.

### 6.4 Manual cleaning

Since mechanical processes can be standardized, reproduced and therefore validated, mechanical cleaning/disinfection should be preferred to manual processes. Manual cleaning and disinfection process is not validated and therefore needs to be validated additionally by the end user.


### 6.5 Mechanical cleaning

Due to international standards (EN ISO 15883) and national guidelines, only validated mechanical cleaning and disinfection procedures should be used. For automated cleaning, we recommend a standard program for surgical instruments, e.g. instruments from Miele.

Observe the following by loading:

- Place the disassembled or opened instruments securely in the tray.
- Instruments with openings and recesses must be placed with the open side facing downwards so that they can be cleaned and no water from the cleaning process can collect in them.
- Do not overload trays, avoid creating any overlaps

- If available, use a coordinated rinsing device.

Preliminary rinsing is followed by chemical. The chemical cleaning should take place at **40°C - 60°C** for **at least 5 minutes**.

We recommend products with a **pH-value within 9-10**, e.g. Neodisher MediClean forte from Dr. Weigert. The cleaning agents used should be selected depending on the material and properties of the instruments and in accordance to national regulations: If there is a high chloride concentration in the water, pitting and tension crack corrosion can occur on the instruments. The occurrence of this type of corrosion is minimized by using alkaline cleaning agents and demineralized water. By adding an acid-based neutralization medium, the rinsing off of alkaline cleaning agent residues is facilitated during the first intermediate rinsing process.

In order to prevent the formation of deposits, it is advisable to use neutral cleaners where the water quality is unfavourable.

After the second intermediate rinsing process, thermal disinfection takes place.

The thermal disinfection should be carried out with demineralized water that corresponds to microbiological drinking water quality (< 100 CFU/ml drinking water) at **80 to 95°C** and an exposure time in accordance with **EN ISO 15883**.

After the finishing of the programme take the good out of the machine because corrosion can arise if the instrument remains in the machine.

<b>V Parameters used for the validation of preparation</b>	
<b>Pre-rinsing</b>	1 minute with cold tap water
<b>Cleaning</b>	Temperature: 55 °C
	Soaking Time: 5 minutes (worst case)
	Neodischer Mediclean forte 0,4% (worst case)
<b>Neutralization</b>	Temperature: cold DI water
	Soaking Time: 2 minutes
	Neodisher Z 0,1%
<b>Post-rinsing</b>	2 minutes with cold DI water
<b>Disinfection</b>	Temperature: 93°C ( $A_0 > 3000$ )
	Soaking Time: 5 minutes

## 6.6 Drying

Ensure adequate drying by the cleaning and disinfection device or using other suitable measures.

<b>V</b>	During validation, the drying temperature was 100°C and the drying time was 20 minutes.
----------	---

## 7 Maintenance, inspection

After cooling to room temperature, the instruments must be visually inspected for protein residues and other contamination. Slits, barriers, locks, tubes and other areas that are difficult to access must be thoroughly inspected. Instruments that are not residue-free must be repeatedly subjected to the entire reprocessing process.

To ensure that surgical instruments can be used for their intended purpose after reprocessing, it is necessary to carry out a functional test after cleaning, disinfection and drying of the visual inspection and care measures. Carry out the functional tests described in chapter 7.1.

Instruments that are stained, blunt, bent, no longer functional, damaged in any other way or the instrument labeling is no longer recognizable (e.g. labeling, UDI code no longer legible) must be discarded!

To help identify faulty instruments that need to be sorted out, we recommend the brochure "Instrument Reprocessing" from the Working Group "Instrumenten Aufbereitung". This includes Chapter 8 "Checks and Care" and Chapter 12 "Surface Changes: Deposits, Discoloration, Corrosion, Aging, Swelling and Stress Cracks".

### 7.1 Functional check

A newly purchased product must be subjected to a thorough visual and function check after its delivery and before each use.

Products must be checked for irregularities. Paying attention to cracks, fractures and the occurrence of corrosion.

If there are joints, the instruments should be oiled with a care product before the functional test. We recommend a medical white oil based on paraffin oil.

Check instruments with joints for ease of movement.

The labeling must also be checked for legibility.

Carry out a function check in accordance with the intended application of the instrument.

Defective products must not be used and must have undergone the complete preparatory treatment process again before being returned.

## 8 Sterilization

Prior to sterilization, products must undergo cleaning and disinfection, be rinsed off without residue using demineralized water and subsequently dried. HEBUmedical recommends using a validated steam sterilization process (e.g. sterilizer in compliance with EN 285 and validated in accordance with DIN EN ISO 17665-1).

The validated parameters refer to reusable surgical instruments. The validated parameters should also be observed for the other products described, unless a different procedure is explicitly described.

On using the fractionated vacuum method, **sterilization** must be performed with at least **134°C (USA 132°C)** with a **minimum dwell period of 3 minutes**. Vacuum drying must then be carried out for at least 20 minutes.

<b>V</b> Parameters used for the validation of steam sterilization	
<b>Prevacuum</b>	3 times
<b>Sterilization temperature</b>	132 °C
<b>Sterilization time</b>	1,5 minutes (half cycle method)
<b>Drying time</b>	Drying time

The vapour must be free of ingredients, recommended limiting values of feed water and vapour condensate are determined through EN 285.

Other sterilization processes are compatible but not validated from HEBUmedical.

When loading, observe the recommended total weight. After the sterilization, check the sterile product packaging for damage, and inspect the sterilization indicators.

### 8.1 Packaging

Compliant packaging of products for sterilization in line with ISO 11607. Packaging used must be suitable for the instruments and protect them from microbiological contamination during storage. The seal must not be under tension. HEBUmedical recommends container or hospital common sterilization paper/film packagings for sterilized packaging.

<b>V</b> During validation the instruments were packaged in hospital common sterilization packagings (paper/film packagings) and steam sterilized.
--

## 9 Lifetime

The steam sterilization procedure was validated by laboratory tests. The products were sterile validated at a pre-vacuum of at least 5min duration and a temperature of 134°C for a lifetime of 50 cycles.

You can continue to use the instruments at your own responsibility over this cycle value if the tests described in chapter 7 have been successfully completed.

## 10 Storage

Store products in a dry, clean and dust-free environment at moderate temperatures from 5°C to 40°C.

Protect from the effects of the sun's rays and artificial light.

## 11 Warranty / Repair

Our products are manufactured from high-grade materials and carefully checked prior to dispatch. However, even if used properly in accordance with their intended purpose they are subject to a greater or lesser degree of wear depending on their intensity of use.

This wear is technically induced and unavoidable.

Should faults occur independently of wear, please contact our customer services. Defective products should no longer be used.

They must undergo the complete preparatory treatment process before being returned.

## 12 Disposal

Before disposal, the product must undergo the complete reprocessing and sterilization process in accordance with the procedure described in these instructions for use.

**Risk of infection from contaminated products!**

- When disposing of the product, its components and their packaging, the national regulations must be observed.

**Risk of injury from sharp-edged and/or pointed products!**

- When disposing of the product, ensure that the packaging prevents injury from the product.

## 13 Service and manufacturer address

Should you require the instructions for use in paper form, please use the contact details below. The instructions for use in paper form will be made available to you within seven calendar days of receipt of the request.

Alternatively, you can print out the electronic instructions for use yourself.

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## **Article list**

### Rectal probes; Suction Tubes, cannula:

HB1245-02, HB2590-53, HB2591-43, HB2591-44, HB2591-53, HB2591-65, HB2592-33, HB2593-22, HB2593-33, HB2594-12, HB2594-33, HB2594-54, HB2594-63, HB2594-65, HB2594-66, HB2596-32, HB2596-33, HB2596-34, HB2596-53, HB2596-54, HB2596-55, HB2596-63, HB2596-64, HB2596-65, HB2597-44, HB2597-54, HB2597-64, HB2597-65, HB2665-03, HB2665-04, HB2665-05, HB2665-06, HB2665-07, HB2665-08, HB2665-09, HB2665-10, HB2665-11, HB2665-12, HB2665-13, HB2665-14, HB2665-15, HB2665-16, HB2665-17, HB2665-18, HB2665-19, HB2665-20, HB2665-22, HB2665-23, HB2665-24, HB2665-25, HB2665-26, HB2665-27, HB2665-28, HB2665-29, HB2665-30, HB2665-32, HB2665-33, HB2665-34, HB2665-35, HB2665-36, HB2665-37, HB2665-38, HB2665-39, HB2665-40, HB2665-42, HB2667-00, HB2667-05, HB2667-07, HB2667-10, HB2667-15, HB2667-17, HB2667-20, HB2667-22, HB2667-23, HB2667-25, HB2667-27, HB2667-30, HB2667-32, HB2667-33, HB2667-35, HB2667-37, HB2667-40, HB2667-42, HB2667-45, HB2667-47, HB2667-50, HB2667-55, HB2667-57, HB2667-60, HB2667-65, HB2667-67, HB2668-00, HB2668-04, HB2668-08, HB2668-12, HB2668-13, HB2668-15, HB2668-38, HB2668-40, HB2668-42, HB2669., HB2669-00, HB2669-05, HB2669-07, HB2669-10, HB2669-12, HB2669-20., HB2669-21, HB2669-27, HB2669-29, HB2669-85, HB2669-86, HB2669-87, HB2669-88, HB2669-89, HB2669-90, HB2669-91, HB2669-92, HB2669-95, HB2670-03, HB2670-04, HB2670-05, HB2670-06, HB2670-07, HB2670-07L, HB2670-08, HB2670-09, HB2670-10, HB2670-10L, HB2670-11, HB2670-12, HB2670-12L, HB2670-14, HB2670-15, HB2670-25, HB2670-26, HB2670-27, HB2670-28, HB2670-29, HB2670-30, HB2670-31, HB2670-32, HB2670-35, HB2670-45, HB2670-46, HB2670-47, HB2670-48, HB2670-49, HB2670-50, HB2670-51, HB2670-52, HB2670-55, HB2670-65, HB2670-66, HB2670-67, HB2670-68, HB2670-69, HB2670-70, HB2670-71, HB2670-72, HB2670-74, HB2670-75, HB2670-86, HB2670-87, HB2670-88, HB2670-89, HB2670-90, HB2670-91, HB2670-92, HB2670-95, HB2672-06, HB2672-07, HB2672-08, HB2672-09, HB2672-10, HB2672-11, HB2672-12, HB2672-14, HB2672-15, HB2672-17, HB2672-20, HB2672-25, HB2672-27, HB2673-05, HB2673-06, HB2673-07, HB2673-08, HB2673-09, HB2673-10, HB2673-18, HB2673-20, HB2673-22, HB2673-23, HB2673-24, HB2673-25, HB2673-26, HB2673-27, HB2673-33, HB2673-35, HB2673-37, HB2673-39, HB2674, HB2674-01, HB2674-02, HB2674-30, HB2674-40, HB2674-50, HB2674-60, HB2675-06, HB2678-01, HB2678-02, HB2678-03, HB2678-04, HB2678-05, HB2679-00, HB2680, HB2680-21, HB2680-33, HB2681, HB2682, HB2683, HB2683-24, HB2684-01, HB2684-02, HB2685, HB2685-01, HB2686-01, HB2686-02, HB2686-03, HB2687-01, HB2687-09, HB2687-10, HB2687-21, HB2687-22, HB2687-31, HB2687-32, HB2687-41, HB2687-42, HB2688-12, HB2688-21, HB2688-22, HB2688-23, HB2688-41, HB2688-42, HB2688-51, HB2688-52, HB2688-55, HB2688-61, HB2688-71, HB2688-75, HB2688-81, HB2700-20, HB2700-21, HB2700-22, HB2700-23, HB2700-24, HB2700-25, HB4272-00, HB4277-20, HB4277-25, HB4277-30, HB4277-40, HB4278-00, HB4278-23, HB5648-05

### Biopsy Cannulas

HB2676-01, HB2676-02, HB2676-03, HB2676-04, HB2676-05, HB2676-06, HB2676-07, HB2676-08, HB5060-01, HB5060-02, HB5060-03, HB5060-04, HB5060-05, HB5067

### Metal Catheters

HB4599-01, HB4599-02, HB4599-03, HB4600-08, HB4600-10, HB4600-12, HB4600-14, HB4600-16, HB4600-18, HB4600-20, HB4600-22, HB4600-24, HB4602-06, HB4602-08,

HB4602-10, HB4602-12, HB4602-14, HB4602-16, HB4602-18, HB4602-20, HB4602-22,  
HB4602-24

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HEBUmedical GmbH