Endoskope /

Endoscopes

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Seite 3

Operating Manual

Side 23

DEUTSCH

ENGLISH

Inhalt

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Table of contents

1	Symbol descriptions	25
2	Introduction	26
3	Intended use	26
3.1	Intended purpose	26
3.2	Indication	26
3.3	Contraindications	26
3.4	Patient target group	26
4	Warnings	27
4.1	Risk of infection/complications	28
5	Handling	29
6	Assembly and disassembly	29
6.1	Rigid endoscopes	29
6.2	Fibre-optic connection	30
6.3	Stopcocks	31
6.4	Camera head	32
6.5	Instrument coupling	32
7	Preparation	33
7.1	Information on instrument preparation	33
7.2	Basic principles	34
7.3	Cleaning and disinfection	34
7.4	Pre-treatment	35
7.5	Manual cleaning	35
7.6	Mechanical cleaning and disinfection	36
7.7	Checking	37
8	Maintenance, inspection	38
8.1	Functional check	38
9	Sterilisation	39
9.1	Preparation and packaging for sterilisation	39
9.2	Steam sterilisation	39
9.1	Storage	40

10 Life time	40
11 Storage	40
12 Warranty / Repair	41
13 Disposal	41
14 Service and manufacturer address	41

1 Symbol descriptions

Symbols	Definition
	Attention
	Validated Parameters
	Manufacturer
	Lot-description
	Reference code
	Medical device / FDA Prescription device
	Medical device
	Non sterile
	Keep away from sunlight
	Dry storage required
  Hinweis auf eIFU	(Electronic) instruction for use

2 Introduction

By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following.

In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

Our products are intended exclusively for professional use in accordance with their intended purpose by appropriately trained and qualified personnel and may only be purchased by such personnel.

3 Intended use

3.1 Intended purpose

The endoscopes are intended for temporary or short-term illumination and visualization during diagnostic and therapeutic endoscopic, minimally invasive, surgical procedures. It is a reusable instrument and is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.

3.2 Indication

The instrument is used in various surgical procedures during arthroscopy, nephroscopy, spino-scopy, laparoscopy or thoracoscopy. The procedure must be performed by trained and qualified specialists.

3.3 Contraindications

No contraindications are known. There are no known contraindications directly related to using an endoscope. In principle, the use of rigid endoscopes for a single endoscopic procedure are contra-indicated if endoscopic procedures in general are contraindicated.

Only materials suitable for use in the medical field are used for the rigid endoscopes. In rare cases, hypersensitised patients may experience pseudo-allergic reactions during prolonged contact. Having this in mind, you must ensure the patient does not have an allergy to silicone, nickel and possibly brass before the endoscopic procedure.

3.4 Patient target group

The products are suitable for all patient groups.

4 Warnings

	Medical products are delivered in a non-sterile condition and must be cleaned, disinfected and sterilized prior to their initial use.
	The use of faulty instruments is in principle forbidden and they have to go through the whole cleaning process before return.
	Please take into consideration that through higher power a bigger damage of the tissue can result: f.e. on forceps: the power at the end of the jaw is higher than at the tip of the jaw
	Remove all protective sleeves and films prior to first using or preparation for use.
	The safe combination of different products or of products with implants must be reviewed prior to clinical application by the user.
	Avoid improper throwing or dropping of instruments.
	Avoid mechanical overstressing of the instrument beyond the structural design, this can lead to breakage and deformation!
	A visual inspection of the instrument for damage and contamination must be carried out before each use!
	To prevent all contact corrosion, instruments with damaged surfaces must be separated immediately.
	When using the endoscope in a trocar, avoid any bending stress during insertion and removal.
	Some parts of the endoscope may heat up significantly during operation, there is a risk of burns. Appropriate protective equipment must be worn to avoid injury. The light source should only be used to the extent necessary.
	Using the endoscope in a clinical setting increases the risk of infection. Therefore, all precautions must be taken to prevent infections.
	As light energy and heat energy are released, the distal end and the light guide terminal can heat up significantly during use. Direct contact with fabric and highly flammable materials must be avoided. If possible, use only the brightness level that is really required and not the maximum lighting setting.
	When using RF electrodes, it is important to ensure that the active electrode is always in the user's line of sight and does not come into contact with the endoscope or other metal parts of the instrument
	During laser surgery, the use of reflective objects in the work area should be avoided, and the laser beam should never be directed at the endoscope.

	Endoscopes must be prepared in accordance with the instructions for use before the first use and before each subsequent use.
	To avoid infections, all employees must wear the following personal protective equipment: a protective hood to completely cover the hair, as well as eye, mouth and nose protection, gloves, a protective gown and suitable, moisture-proof shoes. Further hygiene guidelines of the facility must be followed.
	Damaged endoscopes should never be used. (see section 2.8 – “Visual and functional inspection”)
	When using chemical cleaning agents, the occupational health and safety instructions on the respective safety data sheet must be observed.
	The products must not be cleaned in an ultrasonic bath.
	Never use hot air sterilisation, flash sterilisation, or radiation sterilisation.
	If the products are used on patients with transmissible spongiform encephalopathy or HIV infection, we decline any responsibility for their reuse.
	After ophthalmical use, please pay attention to water quality during treatment (according to the specifications of ANSI/AAMI ST108 and the recommendations of the Rober-Koch-Institute on preparation of medical devices)!
	Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

4.1 Risk of infection/complications

Depending on the area of application, complications that may arise due to the result of the procedure must be considered.

The risk of infection during endoscopic examinations is particularly important.

Risk factors causing an increased risk of infection may be divided into two categories:

Procedure-related risks:

- the nature and extent of tissue damage during therapeutic interventions.
- circumstances surrounding the endoscopic intervention (emergency or planned procedure).
- incompetence and low experience of the examiner/user.
- incorrect cleaning and disinfection of the endoscopes and accessories.

Patient-related risks:

- patients with impaired immune status or immunosuppression (HIV, leukaemia, lymphoma, immunosuppressive therapy, advanced liver or kidney disease, old age).
- presence of specific sources of infection or anatomical conditions.
- conditions facilitating the adhesion of bacteria in the organism (heart valve defect, heart valve replacement, endoprotheses, permanent intravenous catheters).

Endoscopic examinations may cause endogenous displacement of the body's own microorganisms, followed by bacteraemia. In this context, the national and international recommendations regarding prophylactic administration of antibiotics before certain interventions must be observed (ESGE Guidelines).

5 Handling

The type of treatment must be determined in each individual case by the surgeon in cooperation with the internist and the anaesthetist.

For operational use in various surgical disciplines must be done by appropriately trained and qualified personnel.

6 Assembly and disassembly

6.1 Rigid endoscopes

Models with working channels particularly, the (inner lumen) require very thorough cleaning to prevent the build-up of deposits in the delicate channels.

During cleaning and disinfection, all detachable endoscopic attachments must be removed to uncover any concealed surfaces. This is the only way to ensure adequate reprocessing.

6.2 Fibre-optic connection

Disassembly:

Unscrew adapters (1) and (2) and remove them from the endoscope

In case of working channels – if applicable –:


- Remove sealing cap (5)
- Unscrew valve body (4)
- Remove the valve (3)

(see picture above)

Ventile valves and valve caps are disposable items for single use. Multiple use is not intended and is not validated.


Assembly:

- Screw on adapters (1) and (2) adapter.
- In case of working channels (if applicable)
 - Insert a new valve.
 - Screw on the valve body (4)
 - Attach the sealing cap (5) (see picture above)

(1) Storz® - / Aesculap® - / Olympus® - adapter,
(2) Wolf® - adapter,
(3) ACMI® - connection.

6.3 Stopcocks

- (6) Stopcock,
(7) Conical valve,
(8) Spring cap.


Disassembly:

- Unscrew the spring cap (8) and remove the conical valve (7) from the stopcock (6).

Assembly:

- Prior to every sterilisation, the conical valve (7) must be treated with a lubricant that is approved for the purpose in question and the respective reprocessing technique. This is to provide corrosion protection and to retain the functional capacity.
- Insert the conical valve (7) into the stopcock (6) and screw on the spring cap (8). Check the stopcocks to see whether they function correctly.
- When inserting the conical valve, make sure that the guide pin runs smoothly in the guide and that the lever points towards the opening when it is open.

6.4 Camera head


- Unlock the camera head (1) as shown in the picture and put it on the eyepiece funnel (2) and lock it (3).

Due to the different compatible systems, the instructions for use from the manufacturer must be followed. The eyepiece funnel for coupling the camera head complies with the ISO/TS 18339 specifications.

6.5 Instrument coupling

- Unlocking the camera head (1), placing it on the eyepiece funnel (2) and locking the camera head (3). as shown in the picture.

7 Preparation

The person in charge of preparatory treatment is responsible for ensuring that the treatment is duly carried out using the relevant equipment, materials and personnel in the treatment facility and so achieves the desired result. This necessitates validation and routine monitoring of the process used. We urge you to take note of the national regulations dealing with instrument preparation.

7.1 Information on instrument preparation

- Detergents and disinfectants can cause considerable damage to endoscopes. These agents must not contain the following components:
 - organic acids, mineral acids, and oxidising acids (minimum pH 5)
 - strong alkalis (maximum pH 10)
 - phenols or halogens (such as chlorine, iodine, bromine)
 - aromatic/ halogenated hydrocarbons
- Agents used in combination must be compatible with each other.
- Use cleaning and/or disinfection agents with a pH-value within 9-10
- Please observe manufacturer instructions regarding dosage, exposure time and renewal of solutions.
- Use suitable brushes that do not remove material, i.e. **no** hard brushes (such as metal brushes and metal sponges) or coarse abrasive cleaners.
- Never leave instruments in cleaning or disinfection agents for longer than the specified time.
- Rinse and dry carefully through channels and pipes.
- Sensitive instruments must be cleaned in a storage or clamping fixture.
- Observe manufacturer instructions of cleaning – and sterilizing equipment.
- Never forcefully accelerate the endoscope cooling process (for instance with cold water); sudden changes in temperature may destroy the optical components.
- The endoscopes may only be exposed to temperatures below 137°C (279°F).
- Never clean endoscopes in an ultrasonic bath (damage to the optical system).

7.1.1 Water quality

For cleaning, neutralization and rinsing, we recommend the use of demineralized water in accordance with the "Guideline DGKH, DGSV, AKI for the validation and routine monitoring of automated cleaning and thermal disinfection processes for medical devices and on the principles of device selection" (the guideline refers to DIN EN ISO 15883-1 point 6.4.2).

According to the guideline, the following values are recommended:

- Conductivity: $\leq 15 \mu\text{S}/\text{cm}$ (deviating from the table in DIN EN 285)
- pH value: 5 - 7

- Total hardness: ≤ 0.02 mmol CaO/l
- Salt content: ≤ 10 mg/l
- Phosphate (as P_2O_5): ≤ 0.5 mg/l
- Silicate (as SiO_2): ≤ 1 mg/l
- Chloride: ≤ 2 mg/l

7.2 Basic principles

The endoscopes must be cleaned, disinfected, and sterilised before and after every use, especially brand-new devices, since all endoscopes are delivered in a non-sterile condition (cleaning and disinfection after removing the transport packaging; (sterilisation in suitable sterilisation packaging) is very important.

The following conditions must be met to ensure effective reprocessing:

- Cleaning and reprocessing immediately after use (max. 1h after the last use) to avoid encrustation and corrosion during dry removal. If they are expected to be left for a while, the optical systems are to be kept in a detergent/disinfection solution to stop them drying up.
- Wet removal should be done only in non-corrosive solutions at the specified concentration.
- Determining the configurations for the handling of the used devices and compliance with the manufacturer's instructions.
- Regular maintenance and checking of the used devices.
- Validated procedures for all reprocessing steps.
- Adherence to the standardised parameters for each reprocessing cycle.
- Verifying the disinfection and sterilisation efficiency based on corresponding indicators.

In addition, the applicable national hygiene regulations and the local guidelines of the doctor's practice or hospital must be observed, especially the various requirements regarding effective prior inactivation.

7.3 Cleaning and disinfection

The mechanical method described here for the washer-disinfector should be given priority. Effective cleaning and disinfection are an essential prerequisite for effective sterilisation. The pre-treatment must always be carried out separate from this (see chapter 7.4).

Notes:

- Only use cleaning products that have been tested and approved in accordance with the respective national hygiene regulations and local guidelines.
- In endoscopes with a channel system (irrigation and/or working channels), the inner lumina must be carefully cleaned and disinfected to prevent fixation and preservation of organic residues by aldehyde.

- If cleaning is not carried out properly, there is a risk of infection. Contamination by germs must be prevented.
- To avoid infections, all employees must wear personal protective equipment. Further hygiene guidelines of the facility must be followed.

7.4 Pre-treatment

Working steps immediately after application:

- Remove all fibre-optic adapters before reprocessing! If applicable, disassemble all stopcocks (see chapter 6.3).
- All cleaning steps must be performed below the surface of the liquid to prevent squirting of contaminated liquids.
- Rinse thoroughly with cold running purified water (max. 20°C) to remove coarse impurities from the endoscopes.
- Remove adhering contaminations with a mild detergent solution approved for medical endoscopes (see chapter 7.1).
- All empty channels must be rinsed at least five times using a disposable syringe (with a volume of at least 50 ml).
- Never use abrasive cleaners or metal brushes and avoid excessive force when manually removing contaminants.
- Finally rinse the endoscopes with purified water (demineralised water in accordance with DIN EN ISO 15883-1) for one minute to prevent staining and pitting.
- Dry fully with compressed air (especially cavities) or a lint-free cloth.

Evidence of the basic suitability of the endoscopes for effective cleaning was provided by an independent accredited test laboratory using Neodisher Mediclean detergent (Dr. Weigert GmbH & Co. KG, 20539 Hamburg) with an immersion time of 5 minutes in a 0.5% solution as part of the validation of the mechanical reprocessing technique. In doing so, the procedure described above was used.

7.5 Manual cleaning

Since mechanical processes can be standardized, reproduced and therefore validated, mechanical cleaning/disinfection should be preferred to manual processes. Manual cleaning and disinfection process is not validated and therefore needs to be validated additionally by the end user.

7.6 Mechanical cleaning and disinfection

Requirements for suitable cleaning and disinfecting equipment:

- Program selection for optimised endoscope cleaning with enough rinsing cycles. The manufacturer's specifications regarding loading of the receiving wash basins must not be exceeded.
- The device must be equipped with suitable racks and connections that allow safe cleaning and disinfection with the enabled program. This must be verified by the operator of the washer-disinfector.
- Controlled thermal disinfection program (A0 value > 3000 or at least 5 min. at 90°C) with proven efficacy.
- Regular maintenance and proven efficacy. The devices must meet the EN ISO 15883-1:2014 requirements. Country-specific requirements must be observed.
- Final rinsing with cold demineralised water (in accordance with DIN EN ISO 15883-1) for at least 120 sec.
- Controlled drying phase. The program must be verified by the operator of the washer-disinfector.

With chemo-thermal disinfection there is a risk of disinfectant residues remaining on the endoscopes.

A quarterly microbiological check for quality assurance of the reprocessing method is recommended.

Requirements for suitable detergents/disinfectants:

- Approval for the cleaning of endoscopic instruments with proven efficacy
- Compatibility of the used detergents/disinfectants (especially in the case of chemo-thermal disinfection).
- Listed chemicals (see chapter 7.1) should be avoided.
- If powders are used, it must be ensured that they have been safely and completely dissolved before cleaning. Potential residues may, for instance, block the inner lumens.
- Use enzyme-based agents with a neutral pH value.

Increased chloride concentrations in the feed water cycle may damage the material (pitting). The rinsing water must be carefully reprocessed to prevent recontamination.

The manufacturer's instructions for the cleaning agent and, if applicable, disinfectant regarding concentration, temperature and exposure time must be observed.

Procedure:

- Safely attach the endoscopes to the inserts of the disinfection unit. Avoid rinse residue and do not allow the endoscopes to be in contact with other instruments.
- Open the stopcocks – if applicable – and connect all lumens of the endoscopes to the special inserts with a rinsing system to guarantee complete and thorough rinsing of all cavities.
- Do not overload the disinfection unit.
- Start the program.
- After completing the program, check whether the program was run in accordance with the specifications and whether all control parameters are met.
- As soon as the automatic cleaning cycle is completed, immediately remove the endoscopes from the disinfection unit to avoid corrosion. Make sure your hands are disinfected or you are wearing fresh disposable gloves. Avoid accelerated cooling (for instance in water).
- Dry the tubes and channels with compressed air and, if necessary, wipe the devices dry with a lint-free cloth.
- Checking and maintenance (see chapter 7.7 and 8).
- Packaging the endoscopes (see chapter 9.1).

V

Evidence of the basic suitability of the endoscopes for effective automatic cleaning and disinfection was provided by an independent accredited test laboratory using the Neodisher Mediclean cleaner (0.5% solution, 5min immersion) for pre-cleaning and then the Lautenschläger ZentraCert disinfectant (thermal disinfection), and the Neodisher Mediclean detergent (Dr. Weigert GmbH & Co. KG, Hamburg) in the Miele PG 8535 washer-disinfectant (cleaning, Miele Professional, Princeton), as part of automatic reprocessing validation. In doing so, the procedure described above was used.

7.7 Checking

After cleaning and disinfection, the endoscope can then be assembled (see chapter 6). The following checks must be carried out:

- Visual check of the three optical surfaces (see chapter 8.1), and, if necessary, cleaning with an alcohol-soaked swab (70%). Deposits on the light guide can cause considerable loss of light and may also affect the optical system.
- Checking the surfaces for corrosion, wear, general damage, sharp edges, or chipping in the distal area.
- If there are still any residues or contaminations left, repeat the disinfection after manual pre-treatment.
- Damaged endoscopes must be discarded.

8 Maintenance, inspection

After cooling to room temperature, the instruments must be visually inspected for protein residues and other contamination. Slits, barriers, locks, tubes and other areas that are difficult to access must be thoroughly inspected. Instruments that are not residue-free must be repeatedly subjected to the entire reprocessing process.

After each cleaning and disinfection cycle, the stopcocks must be lubricated before sterilisation (see chapter 6.3). Only lubricants with verified biocompatibility should be used. The lubricant must be suitable for this application and approved for steam sterilisation.

Cleaning the optical surfaces with 70% alcohol (ethanol, isopropanol) prevents deposits from setting/burning in.

To ensure that surgical instruments can be used for their intended purpose after reprocessing, it is necessary to carry out a functional test after cleaning, disinfection and drying of the visual inspection and care measures. Carry out the functional tests described in chapter 8.1.

Instruments that are stained, blunt, bent, no longer functional, damaged in any other way or the instrument labeling is no longer recognizable (e.g. labeling, UDI code no longer legible) must be discarded!

To help identify faulty instruments that need to be sorted out, we recommend the brochure "Instrument Reprocessing" from the Working Group "Instrumenten Aufbereitung". This includes Chapter 8 "Checks and Care" and Chapter 12 "Surface Changes: Deposits, Discoloration, Corrosion, Aging, Swelling and Stress Cracks".

8.1 Functional check

A newly purchased product must be subjected to a thorough visual and function check after its delivery and before each use.

Products must be checked for irregularities. Paying attention to cracks, fractures and the occurrence of corrosion.

If there are joints, the instruments should be oiled with a care product before the functional test. We recommend a medical white oil based on paraffin oil.

Check instruments with joints for ease of movement.

The labeling must also be checked for legibility.

Carry out a function check in accordance with the intended application of the instrument.

Examination for:

- external damage (shaft deformation, dents, or sharp edges).
- residues of detergents or disinfectants. Condition of the three optical surfaces – 1. Lens window, 2. Ocular window, 3. Fibre-optic connection - using reflecting light or a magnifier (smooth, clean, and intact).
- optimum image quality (a bright and clear image with high definition).
- loss-free light transmission from the fibre-optic connection to the light output (Possible comparison with a new device).
- free passage of the endoscope's working channels.
- material changes on the metal and plastic surfaces.
- legibility of the device labelling.

Defective products must not be used and must have undergone the complete preparatory treatment process again before being returned.

9 Sterilisation

Proceed as follows before sterilising the components.

9.1 Preparation and packaging for sterilisation

- Open all stopcocks if applicable.
- Only use disposable sterilisation packaging and/or sterilisation containers that are suitable for steam sterilisation (sufficient temperature resistance, air, and steam permeability according to DIN EN ISO 11607).
- The packaging must guarantee optimum protection of the sterile endoscopes during transport and storage. The storage location must be dust-free, low on germs, dry, dark, and free from temperature fluctuations.
- Reusable sterilisation containers must be maintained in accordance with the manufacturer's specifications; they must offer secure fixation of the endoscopes and protect them from damage.

Since the suitability of the packaging has a significant impact on the sterilisation result, this must be checked when defining the sterilisation parameters.

Make sure that only completely cleaned, well maintained, dry and disinfected devices are being sterilised in the user's process cycle.

The following sterilisation procedure was validated for its germicidal effect:

9.2 Steam sterilisation

- Fractionated vacuum process (with 3x pre-vacuum) for endoscopes with and without an empty channel.

Steam sterilisation conditions:

- Sterilisation temperature: from min. 132°C (270°F) to max. 134°C (273°F); in accordance with DIN EN ISO 17664
- Sterilisation time at sterilisation temperature: 3-5 min at 132°C (270°F)
- Steam steriliser, approved according to DIN EN 13060 and DIN EN 285 and tested according to DIN EN ISO 17665
- Observe the cooling time. Accelerated cooling (e.g., with cold water) may destroy the endoscopes.
- According to the KRINKO/BfArM-RKI recommendations, saturated steam sterilization at 134

The steriliser manufacturer's instructions regarding preparation, packaging and execution of the sterilisation process must be carefully observed.

Other sterilisation procedures are prohibited!

V

vidence of the basic suitability of the endoscopes for effective steam sterilisation was provided by an independent accredited test laboratory using Lautenschläger Zentra-Cert (F. & M. Lautenschläger GmbH & Co. KG, Köln) in a half-cycle process (1.5 min/132°C) within the scope of the validation of the sterilisation technique. Typical conditions in clinics and medical practices were assumed and the endoscopes were tested under laboratory conditions in a half-cycle process. This also pertains to other parameters with a longer holding time and/or higher temperature.

9.1 Storage

The packaging must guarantee optimum protection of the sterile endoscopes during transport and storage. The requirements of the relevant standard DIN EN ISO 11607 must always be complied with.

10 Life time

No damage is to be expected if used as intended.

If the instructions for use are observed, the number of reprocessing cycles has only a minor effect on the service life of the products. The reprocessing methods as well as the products used for this purpose have a significant impact on the service life. When used as intended, no damage to the product is expected. A visual inspection must therefore be carried out before each use (siehe chapter 8.1).

11 Storage

Store products in a dry, clean and dust-free environment at moderate temperatures from 5°C to 40°C.

Protect from the effects of the sun's rays and artificial light.

12 Warranty / Repair

Our products are manufactured from high-grade materials and carefully checked prior to dispatch. However, even if used properly in accordance with their intended purpose they are subject to a greater or lesser degree of wear depending on their intensity of use.

This wear is technically induced and unavoidable.

Should faults occur independently of wear, please contact our customer services. Defective products should no longer be used.

They must undergo the complete preparatory treatment process before being returned.

13 Disposal

Before disposal, the product must undergo the complete reprocessing and sterilization process in accordance with the procedure described in these instructions for use.

Risk of infection from contaminated products!

- When disposing of the product, its components and their packaging, the national regulations must be observed.

Risk of injury from sharp-edged and/or pointed products!

- When disposing of the product, ensure that the packaging prevents injury from the product.

14 Service and manufacturer address

Should you require the instructions for use in paper form, please use the contact details below. The instructions for use in paper form will be made available to you within seven calendar days of receipt of the request.

Alternatively, you can print out the electronic instructions for use yourself.

HEBUmedical GmbH
Badstraße 8

78532 Tuttlingen / Germany
Tel. +49 7461 94 71 - 0
eMail: service@HEBUmedical.de
Website: www.HEBUmedical.de

HEBUmedical GmbH

Badstraße 8 • 78532 Tuttlingen / Germany
Tel. +49 (0) 7461 94 71 - 0
info@HEBUmedical.de • www.HEBUmedical.de